




Patient-friendly Compression



Catalogue & Price List | Issue 10






Code Finder App

Quickly and easily find information for prescribing.
Works on smartphone, tablet or PC.

- No need to use a form - get a **variant identifier** for your configured garment
- EMIS & SystmOne compatible data
- Generate and download a referral letter for your patient to take to their GP



Scan here...





Head-to-toe compression

Which area of the body do you need to treat?

Chest

Page

Comfiwave Breast Band	118 ►
Eto Bras & Vests	130 ►
Prairie Bras	142 ►
Kinesio Tape	166 ►

Abdomen

Abdominal Wrap	168 ►
Comfiwave Abdominal Band	106 ►

Genitals

Comfiwave Shorts	104 ►
Eto Shorts	136 ►
Eto Scrotal Supports	135 ►
Veni Shorts	26 ►

Legs

Circular Knit Socks	24 ►
Flat Knit Socks	18 ►
Comfiwave Below Knee	98 ►
easywrap Leg	74 ►
Clinistretch Bandages	162 ►
Lymphflow IPC	152 ►

Head & Neck

Page

Comfiwave Neck	120 ►
Eto Garments	133 ►
Kinesio Tape	166 ►

Arms

Comfiwave Sleeves	108 ►
Armsleeves	48 ►
easywrap Arm	88 ►
Lymphflow IPC	154 ►

Hands


Microfine Glove	64 ►
Flat Knit Gloves	52 ►
Comfiwave Glove	112 ►
easywrap Hand	86 ►

Thigh

Circular Knit Stockings	24 ►
Flat Knit Stockings	18 ►
Comfiwave Thigh	100 ►
easywrap Thigh + Knee	78 ►
Eto Grace	138 ►

Feet

Microfine Toecap	34 ►
Custom Made Toecaps	35 ►
Comfiwave Toecap	102 ►
easywrap Foot	76 ►



Patient-friendly Compression

Haddenham Healthcare is a family-owned, UK-based medical device company with a global outlook. We specialize in compression therapy for lymphoedema and related conditions, designing products that integrate seamlessly into patients' daily lives.



**Supporting
Therapists**

Page 4 ►



**Supporting
Patients**

Page 5 ►



Our Ethos

Page 6 ►



**Conditions
We Treat**

Page 8 ►

About us

How we work

We listen to patients and healthcare professionals, using real-life feedback for research and development.

We innovate with purpose, creating lasting solutions. We partner long-term, ensuring everyone is supported.



Who we are

We foster a friendly culture where ideas can be shared freely.

We care about our communities, donating generously from our retained earnings to charities annually.

We grow responsibly, committed to our employees, partners, and the environment.



Why it matters

Our ethos drives us to deliver life-changing products. We measure success by whether we improve people's lives.

“We strive to be innovative, trustworthy and sincere.”

Supporting therapists...

Sales

E-mail

Send your order/query to sales@hadhealth.com

Call us

We'd be happy to take your call 9am - 5pm, Mon-Fri.
[01844 208842](tel:01844 208842).

Portal

Order via our user friendly e-commerce portal. Contact us for more information.
portal.haddenham.com

Support

Clinical support service

Book an appointment with our clinical team. They have a wealth of knowledge to support compression therapy solutions.

www.hadhealth.com/contact

Learning

Online modules, webinars and face to face training on our products and more. www.hadhealth.com/learning

Resources

Order forms

Download measurement forms and charts.

www.hadhealth.com/forms

Codes

Lookup codes for prescribing
www.hadhealth.com/code-finder

Social

Connect to our social media channels and mailing list to stay up to date with the latest news and products.



hadhealth.com/subscribe



facebook.com/hadhealth



instagram.com/hadhealth



youtube.com/@haddenham

...and their patients

www.Lymphshop.com



Compression


Online store for patient ordering of **Extra Garments** - encouraging supervised self-management. For example spare pairs of stockings for taking on holiday, or extra colours to suit their wardrobe.

Accessories

Garment glue, kinesiology tape and more.

Resources

Patient-friendly books on Lymphoedema, "Compression Matters" blog and mailing list with offers and promotions.



A patient guide for the recently diagnosed
Copies available to request for your waiting room or reception.

lymphshop.com/what-is-lymphoedema

Our Ethos

Approachability & sincerity

We put people at the heart of everything we do - from having a real person answer the phone, to offering free training sessions led by experienced clinicians.

Our open-door policy ensures every employee or customer feels acknowledged and valued.



Innovation for real-world needs

We develop products that are thoughtfully engineered to address patient and clinical needs. Feedback from healthcare professionals and patients keep us responsive and forward-thinking.


Long-term relationships

We regard our customers and suppliers as partners, collaborating to refine solutions that enhance patient care and outcomes.

Trust is built on reliability, honesty, and genuine collaboration—the values that define us as a family business.

A place people want to work

We foster a friendly, encouraging culture, where continuous improvement is valued. Employees benefit from work-life balance and opportunities for growth, making Haddenham a pleasant place to work.



Social & environmental responsibility

We contribute a percentage of our retained earnings to charitable causes, both locally and internationally to help build a better future for all.

We adhere to recognised best practices in environmental management, ensuring we minimise our impact throughout our operations.

Why it matters

Our goal is to be a trusted innovator in compression therapy, delivering products that patients prefer to wear and which healthcare professionals can depend on.

“By combining innovation, trust and sincerity, we aim to leave a lasting, positive mark on the world benefiting patients, communities, and our shared future.”

Conditions we treat

Lymphoedema / Chronic Oedema

What is the lymphatic system

The lymphatic system is a network of vessels and glands throughout the body that helps fight infection and remove excess fluid. It is responsible for collecting all fluid in the tissues and returning it to the circulatory system, failure of the lymphatics leads to chronic oedema.

The terms chronic oedema and lymphoedema are sometimes used interchangeably, both describe failure of the lymphatics but there are some key differences.




What can go wrong?


Chronic oedema is an umbrella term used to describe swelling of all causes and occurs when fluid accumulates in the tissues and is present for at least 3 months. This can be caused by an underlying problem with the venous or lymphatic systems preventing fluid balance in the tissues.

Lymphoedema is a long-term (chronic) condition that causes swelling in the body's tissues. It can affect any part of the body, but usually develops in the arms or legs.

Unhealthy lymphatics



Healthy/compressed lymphatics



How can compression help?

When the lymphatic system is compromised, pressure within the system causes the vessels to dilate and results in backflow of lymph into the tissues.

The external pressure created when compression is applied, moves fluid from the tissues into the initial lymphatics, increasing lymph transport. Where there is lymphatic failure compression will help to push fluid through the tissues towards areas where drainage is functional.

Vascular Disorders


What is the Venous/Circulatory System?

It comprises of the heart, arteries supplying oxygenated blood to our tissues and organs, the veins return deoxygenated blood back to the lungs and heart. Capillaries are the intersection between arteries and veins where high pressure arterial blood transfers to the venous system causing fluid to leak out of the blood vessels into the interstitium containing nutrients, other substances and debris.

What can go wrong?

When the veins become damaged the valves within are unable to close and return the blood back to the lungs & heart. This results in venous reflux, and venous hypertension.

Prolonged venous hypertension leads to leakage of proteins and fluid into the skin, creating inflammation and subsequent tissue breakdown, a venous leg ulcer.



How can compression help?

Compression works by applying external pressure on the skin, thus compressing the veins. This force helps to improve the function of the valves within the veins, which reduces the amount of fluid that leaks out of the veins and into the tissues.

Compression also supports the muscles, which helps to create a force to push blood through the venous system.

CEAP Classification of venous disorders

Clinical Manifestations, Etiology, Anatomic Distribution, Pathophysiology

C0	No visible or palpable signs of venous disease
C1	Telangiectasias or reticular veins
C2	Varicose Veins
C2r	Recurrent varicose veins
C3	Oedema
C4a	Pigmentation or eczema
C4b	Lipodermatosclerosis or atrophie blanche
C4c	Corona phlebectatica
C5	Healed venous ulcer
C6	Active venous ulcer
C6r	Recurrent active venous ulcer

Lipoedema

What is Lipoedema

Lipoedema is a condition that causes a symmetrical build-up of adipose (fat) often in the lower limbs but can also affect the arms, usually affecting women. It typically spares the hands and feet and is often mistaken for lymphoedema or obesity. The affected areas may feel tender, bruise easily, and appear disproportional. While some people experience pain or discomfort, others may have few symptoms.

How can compression help?

Compression therapy is a key part of conservative management for lipoedema. While it does not reduce the underlying fat tissue, compression can help alleviate symptoms such as pain, heaviness, and tenderness. It supports the tissues, promotes comfort, and can improve mobility. When secondary lymphatic involvement occurs, compression also assists in managing oedema by helping to reduce fluid accumulation and support lymphatic return.

General guidance

You should always aim to treat your patient in line with best practice guidelines. The recommendations in this document are based on best practice guidelines from national and international expert groups, and are within published documents.

Remember that this guide does not replace your clinical judgement. If you are unsure then do seek advice from your peers, or directly from the clinical team at Haddenham Healthcare (clinical@hadhealth.com).


There may be instances where you are unable to confer. In such cases the compression prescribed should be one that is tolerated well by the patient and applies a suitable level of compression for the symptoms being treated.

Hosiery



Medical Compression Hosiery

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard, graduated compression to the lower limb.



Which fabric do I need?

Use our pathway to guide your clinical judgement [Page 16](#) ►



What makes it so patient-friendly?

Turn the page to discover all features [Page 12](#) ►

Sized Flat Knits



Pertex Light

[Page 18](#) ►

Sized Circular Knits



Veni

[Page 24](#) ►



Stride

[Page 32](#) ►

Custom Made



Hosiery

[Page 36](#) ►



Pertex

[Page 20](#) ►



Star Cotton

[Page 28](#) ►



Doktus

[Page 30](#) ►



Toecaps

[Page 35](#) ►



Goldpunkt

[Page 22](#) ►



Microfine Toecap

[Page 34](#) ►

Hosiery Appendix



Sizing information, measuring guide, and full range of options & add-ons.

[Page 39](#) ►



Freedom to choose

Haddenham medical compression hosiery is made from high-quality, long-lasting yarns that deliver RAL standard graduated compression.

Graduated compression

Works with the body's natural circulation - concentrating the highest compression at the ankle and gradually decreasing up the limb.

The perfect fit

The widest range of sized garments available, with only a few simple measurements required.

- 3 leg lengths
- 3 foot lengths
- Up to 12 sizes
- Extra wide option for every size



Sized

Our sized Hosiery is **made-to-order**, with the care and quality that we believe each patient deserves.

We also retain **stock** of frequently ordered variants for an even faster turnaround.

or

Custom Made

We aim to be flexible wherever possible by providing a made-to-measure service that offers unparalleled customisation and fit.

7 versatile fabrics

Flat and Circular knit garments suitable for treatment of all stages of Lymphoedema, Venous disorders and Lipoedema.

Fabric selection guide

See page 16 ►




Quality

We only use high-quality, durable, long-lasting yarns which are guaranteed for 6 months of use.

Compression levels are standardised using the RAL classification system.

Reliable

Stockings that can keep up with the most active of lifestyles, thanks to conforming fabrics and a variety of grip top options.



Natural colours

A wide variety of colours to choose, from a versatile palette of natural tones.



Subtle patterns

Our most versatile fabric Veni also offers stylish pattern options for added texture and additional aesthetic combinations.



Veni

Configurable

Comprehensive options to tailor the garment to the patient's needs. Style like active-wear with **Footless** finish, add a **Fly** to tights, or enhance donnability with **Velcro closures**.

For full range of options see [Hosiery appendix Page 39](#) ►



Doktus

Which fabric?

Please use the below maintenance phase pathway as a guide to which of our fabrics would be most suitable for your patient based on their presentation and symptoms.

How would you describe the level of your patient's oedema?

Early accumulation of fluid which subsides with elevation, with or without pitting.

ISL Stage 0-1

Mild

Tissue swelling rarely reduces with elevation and pitting is present.

ISL Stage 2

Moderate

Pitting is absent with skin thickening fibrosis and excess subcutaneous fat.

ISL Late Stage 2 - Stage 3

Severe

MILD

What is the level of shape distortion?

EXTENSIVE

MODERATE

Is the oedema soft and watery?

YES

NO

How would you describe the texture of the patient's tissue?

WOODY

PUTTY LIKE

What is the condition of the patient's skin?



What is the level of pitting?



Are warty overgrowths with fat lobules present?

NO

YES

HEALTHY SKIN

SENSITIVE

FRAGILE SLIGHT

MODERATE

DEEP

veni

STAR
Cotton

Perlex
Light

STRIDE
Below Knee

Doktus

Perlex

GOLDPUNKT

Page 24 ►

Page 28 ►

Page 18 ►

Page 32 ►

Page 30 ►

Page 20 ►

Page 22 ►

Signs and symptoms

▲ Skin conditions

Fragile Thin skin prone to skin tears

Sensitive the presence of skin conditions where the skin is prone to reactions such as redness and itching.

◆ Textures

Watery soft pliable tissues.

Putty putty consistency represents early fibrotic changes, will pit with firm pressure.

Woody firm woody texture is consistent with advanced fibrotic tissue changes and does not pit.

● Pitting Levels


Slight 1+ 2mm barely detectable immediate rebound

Moderate 2+ 4mm deep <10 second rebound

Deep 3+ 6 mm deep <30 second rebound

Very deep 4+ up to 8mm deep >30 second rebound

● Shape distortion



Mild

Moderate

Extensive

How much compression?

Graduated compression* is highest at the ankle and gradually decreases as the yarns travel up the limb.

When deciding how much compression to apply also consider these 3 factors:

1. The nature of the fabric yarns and the way in which they are knitted together.
2. The fabric's ability to retain its structure under pressure. Measured as the increase in compression, per centimetre increase in the circumference of the limb (mmHg/cm).
3. The resultant force exerted on the limb measured in mmHg.

* Read more about how pressure is distributed on [Page 45](#) ►



Stiffness mmHg/cm

✓ Low = 0-2 | ✓ Medium = 2-4 | ✓ High = 4-6

	Circular Knit				Flat Knit		
	veni	STAR Cotton	Doktus	STRIDE	Pertex Light	Pertex	GOLDPUNKT
Pressure							
Class 1 18 - 21mmHg	✓ 1.1	✓ 1.4		✓ 1.2	✓ 1.5		
Class 2 23-32mmHg	✓ 1.3	✓ 1.5	✓ 1.6	✓ 1.7		✓ 2.1	✓ 2.8
Class 3 34-46mmHg	✓ 1.8		✓ 1.8			✓ 4.3	✓ 3.3
Class 4 49+ mmHg							✓ 4.7

[Page 24](#) ►

[Page 28](#) ►

[Page 30](#) ►

[Page 32](#) ►

[Next page](#) ►

[Page 20](#) ►

[Page 22](#) ►

Pertex Light

Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane

A unique stretchy flat knit that stretches with ease and hugs gently—ideal for fragile or sensitive skin.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.

Effortless

Easy application
Especially for elderly or carers.

Conforming
Great at bridging limb shape discrepancies

Skin kind
Gentle on fragile or sensitive skin



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
Class 1	18-21	Low	Soft watery, slightly pitting oedema with minimal shape distortion. Dependency oedema.	No visible or palpable signs of venous disease	Smooth skin with an increase of enlarged subcutaneous fat tissue

Codes Available on **Drug Tariff** as Custom Made. See www.hadhealth.com/code-finder


	Anklet*	Below Knee	Thigh High	Thigh High + Waistband	Tights	Tights	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
Class 1	090	901 900	902 903	905 904 906	908 907 909		
Unit of Sale	Pair	Pair Pair	Pair Pair	Single Single Single	Single Single Single	Single Single Single	Single Single Single

Sizing

Extra wide sizes are shown in green

Circumferences cm	4	5	6	7	8	9	10	11	12	13	14	15
t Waist	65-75	70-80	70-80	75-85	75-85	80-90	80-90	85-95	85-95	90-105	90-105	95-110
j Buttocks	85-105	90-110	90-110	95-115	95-115	100-120	100-120	105-125	105-125	110-130	110-130	115-135
g Thigh Top	46-65 44-53	48-68 46-56	50-71 48-59	52-74 50-62	54-77 52-65	56-81 54-69	58-85 56-73	60-89 58-77	62-93 60-81	64-98 62-86	66-103 64-91	68-108 66-96
d Below Knee	29-40 28-32	31-42 30-34	32-44 31-36	34-46 33-38	35-48 34-40	37-50 36-42	38-52 37-44	40-54 39-46	41-56 40-48	43-59 41-51	44-62 43-54	46-65 44-57
c Widest Calf	31-39 29-33	33-41 31-35	34-43 33-37	36-45 35-39	37-47 36-41	39-49 38-43	40-51 39-45	42-53 41-47	43-55 42-49	45-58 44-52	46-61 45-55	48-64 47-58
b Ankle	19-21	20-22	21-23	22-24	23-25	24-26	25-27	26-28	27-29	28-31	29-33	30-36

Measurements falling outside of the above chart?
See our Custom Made solutions on Page 36 ►



Length

Straight Distance cm

A-g Floor to Thigh Top

SHORT LEG

REGULAR LEG

LONG LEG

57-67

68-80

81-90

A-d Floor to Below Knee

33-37

38-43

44-48

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

Heel to End of Longest Toe

SHORT FOOT

REGULAR FOOT

LONG FOOT

21-23

23-26

26-29

Standard Options

Standard options must be specified when placing an order.

FOOT			
	Closed toe	Open toe	Footless
GRIP TOP		3cm Strong Plain Below Knee Default	
			5cm Strong Plain Thigh High Default
COLOUR			
	Beige	Black	Light Beige

Special Options



All our special options are fully documented in the Hosiery Appendix Page 39 ►

Increase **panty compression**, switch to **male gusset**, add **fastenings** to aid donning, change the **grip top style**, or adjust the **heel and knee stitching** based on your patients activity levels.

Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane

Firm knit with soft touch making it a go-to option for reliable medium compression needs.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.

Strong

Robust yet soft

Suitable for moderate to severe oedema, but with a soft, medium-density structure.

Post-surgical solution

Ideal post-decongestive or post-liposuction garment



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
 Class 2	23-32	Medium	Moderate shape distortion, deep to very deep pitting.	Telangiectasia, thread, reticular or broken veins. Varicose veins/ recurrent varicose veins.	Uneven skin with indentations in the fat tissue and lipomas.

 Class 3	34-46	High	Obesity-related lymphoedema or lymphoedema where there is excess subcutaneous fibrosis and fat, post liposuction for lymphoedema.	Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica. Healed Venous Ulcer.	Large deformations of fat, lipolymphoedema, large fat nodules and lipomas.
--	-------	------	---	--	--

Codes

Available on Drug Tariff as Custom Made. See www.hadhealth.com/code-finder


											
	Anklet*	Below Knee		Thigh High		Thigh High + Waistband			Tights		
		+ GRIP TOP		+ GRIP TOP		LEFT LEG	BOTH	RIGHT LEG	LEFT LEG	BOTH	RIGHT LEG
 Class 2	092	921	920	922	923	925	924	926	928	927	929
 Class 3	093	931	930	932	933	935	934	936	-	-	-
Unit of Sale	Pair	Pair	Pair	Pair	Pair	Single	Single	Single	Single	Single	Single
UNISEX	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN	WOMEN

Sizing

Extra wide sizes are shown in green

Circumferences cm		4	5	6	7	8	9	10	11	12	13	14	15
t	Waist	65-75	70-80	70-80	75-85	75-85	80-90	80-90	85-95	85-95	90-105	90-105	95-110
j	Buttocks	85-105	90-110	90-110	95-115	95-115	100-120	100-120	105-125	105-125	110-130	110-130	115-135
g	Thigh Top	46-48 44-46	48-50 46-48	50-52 48-50	52-54 50-52	54-56 52-54	56-58 54-56	58-60 56-58	60-62 58-60	62-64 60-62	64-66 62-64	66-68 64-66	68-70 66-68
d	Below Knee	29-31 28-30	31-33 30-31	32-34 31-33	34-36 33-34	35-37 34-36	37-39 36-37	38-40 37-39	40-42 39-40	41-43 40-41	43-45 41-43	44-47 43-45	46-48 44-47
c	Widest Calf	31-33 29-31	33-35 31-33	34-36 33-35	36-38 35-36	37-39 36-38	39-41 38-39	40-42 39-41	42-44 41-42	43-45 42-44	45-47 44-45	46-49 45-47	48-50 47-49
b	Ankle	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29	29-30	30-31

Measurements falling outside of the above chart?
See our Custom Made solutions on Page 36 ►



Length

Straight Distance cm

	SHORT LEG	REGULAR LEG	LONG LEG
A-g Floor to Thigh Top	60-67	68-75	76-82
A-d Floor to Below Knee	33-37	38-42	43-46

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

	SHORT FOOT	REGULAR FOOT	LONG FOOT
Heel to End of Longest Toe	21-23	23-26	26-29

Standard Options

Standard options must be specified when placing an order.

FOOT			
	Closed toe	Open toe	Footless
GRIP TOP		3cm Strong Plain Below Knee Default	
			5cm Strong Plain Thigh High Default
COLOUR			
	Beige	Black	Light Beige Class 2 Only

Special Options



All our special options are fully documented in the Hosiery Appendix Page 39 ►

Increase panty compression, switch to male gusset, add fastenings to aid donning, change the grip top style, or adjust the heel and knee stitching based on your patients activity levels.

Goldpunkt

Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane

Made for stubborn and rebound oedema support, with silky low friction on the inside to make it easy to apply.

Tough

High density body yarns

For high stiffness and oedema control

Smooth inner texture

For easier donning at high compression



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
Class 2	23-32	Medium	Obesity-related lymphoedema or lymphoedema with moderate shape distortion and deep pitting.	Telangiectasia, thread, reticular or broken veins. Varicose veins/ recurrent varicose veins.	Uneven skin with indentations in the fat tissue and lipomas.
Class 3	34-46	Medium	Obesity-related lymphoedema or lymphoedema with moderate to extensive shape distortion,woody and/or very deep pitting tissues.	Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica. Healed venous ulcer.	Large deformations of fat.
Class 4	49+	High	Obesity-related lymphoedema or lymphoedema with extensive shape distortion, fat depositions, fibrosis and warty overgrowths, indicating elephantiasis.	Suitable for use in venous disorders where oedema and shape distortion is extensive.	Large deformations of fat, lipolymphoedema, large fat nodules and lipomas.

Codes

Available on **Drug Tariff** as Custom Made. See www.hadhealth.com/code-finder


	Anklet*	Below Knee	Thigh High	Thigh High + Waistband	Tights							
		+ GRIP TOP	+ GRIP TOP			LEFT LEG	PAIR	RIGHT LEG	LEFT LEG	SINGLE	RIGHT LEG	
Class 2	094	941	940	942	943	945	944	946				
Class 3	096	961	960	962	963	965	964	966				Available as Custom Made
Class 4	098	981	980	982	983	985	984	986				
Unit of Sale	Pair	Pair	Pair	Pair	Pair	Single	Single	Single	Single	Single	Single	

Sizing

Extra wide sizes are shown in green

Circumferences cm	4	5	6	7	8	9	10	11	12	13	14	15
t Waist	65-75	70-80	70-80	75-85	75-85	80-90	80-90	85-95	85-95	90-105	90-105	95-110
j Buttocks	85-105	90-110	90-110	95-115	95-115	100-120	100-120	105-125	105-125	110-130	110-130	115-135
g Thigh Top	46-48 44-46	48-50 46-48	50-52 48-50	52-54 50-52	54-56 52-54	56-58 54-56	58-60 56-58	60-62 58-60	62-64 60-62	64-66 62-64	66-68 64-66	68-70 66-68
d Below Knee	29-31 28-30	31-33 30-31	32-34 31-33	34-36 33-34	35-37 34-36	37-39 36-37	38-40 37-39	40-42 39-40	41-43 40-41	43-45 41-43	44-47 43-45	46-48 44-47
c Widest Calf	31-33 29-31	33-35 31-33	34-36 33-35	36-38 35-36	37-39 36-38	39-41 38-39	40-42 39-41	42-44 41-42	43-45 42-44	45-47 44-45	46-49 45-47	48-50 47-49
b Ankle	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29	29-30	30-31

Measurements falling outside of the above chart?
See our **Custom Made** solutions on Page 36 ►



Length

Straight Distance cm

	SHORT LEG	REGULAR LEG	LONG LEG
A-g Floor to Thigh Top	60-67	68-75	76-82
A-d Floor to Below Knee	33-37	38-42	43-46

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

	SHORT FOOT	REGULAR FOOT	LONG FOOT
Heel to End of Longest Toe	21-23	23-26	26-29

Standard Options

Standard options must be specified when placing an order.

FOOT			
	Closed toe	Open toe	Footless
GRIP TOP		3cm Strong Plain Below Knee Default	
			5cm Strong Plain Thigh High Default
COLOUR			
	Beige	Black	

Special Options

All our special options are fully documented in the **Hosiery Appendix** Page 39 ►

Increase **panty compression**, switch to **male gusset**, add **fastenings** to aid donning, change the **grip top style**, or adjust the **heel and knee stitching** based on your patients activity levels.

Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane

Our most versatile fabric—soft, stretchy, and easy to don and doff, with wide colour and pattern options.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.

Versatile

Conforming

Conforms to limb shape for comfort

Stylish

Available in 10 colours and 5 patterns



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
 Class 1	18-21	Low	Subclinical Oedema. Mild watery oedema. Post lymphatic surgery.	Tired aching legs	Smooth skin with an increase of enlarged subcutaneous fat tissue.
 Class 2	23-32	Low	Swelling that subsides overnight. Slight tissue pitting.	Telangiectasia, thread, broken or reticular veins, varicose veins, recurrent varicose veins.	Uneven skin with indentations in the fat tissue and lipomas.
 Class 3	34-46	Low	Swelling rarely reduces with elevation. Mild shape distortion with healthy skin.	Healed venous ulcer.	N/A

Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder



Anklet*

Below Knee

Thigh High

Thigh High + Waistband

Tights

+ GRIP TOP

+ GRIP TOP

LEFT LEG

BOTH

RIGHT LEG

LEFT LEG

BOTH

RIGHT LEG

 Class 1	050	501	500	505	509	517	513	521	527	525	529
 Class 2	053	533	532	537	541	549	545	553	559	557	561
 Class 3	056	565	564	569	573	581	577	585	591	589	593
Unit of Sale	Pair	Pair	Pair	Pair	Pair	Single	Single	Single	Single	Single	Single

Sizing

Extra wide sizes are shown in green

Stock Options


All other variants made-to-order in 8-10 days



Class	Grip Top	Sizes	Width	Leg & Foot	Toes	Colour
1-2	None	1-8	Standard	Regular	Closed	Beige Black

Circumferences cm	1	2	3	4	5	6	7	8
t Waist	65-75	70-80	75-85	80-90	85-95	90-105	95-110	100-120
j Buttocks	85-105	90-110	95-115	100-120	105-125	110-130	115-135	120-140
g Thigh Top	46-55 43-48	48-59 45-51	52-63 49-55	56-66 53-59	60-69 57-63	64-72 61-66	67-75 64-69	70-78 67-72
c Widest Calf	31-35 29-33	33-37 31-35	35-39 33-37	37-42 35-39	39-45 37-42	41-47 39-45	43-52 41-48	45-56 43-51
b Ankle	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-34

Measurements falling outside of the above chart?
See our Custom Made solutions on Page 36 ►



Length

Straight Distance cm

	SHORT LEG	REGULAR LEG	LONG LEG
A-g Floor to Thigh Top	57-67	68-80	81-90
A-d Floor to Below Knee	33-37	38-43	44-48

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

	SHORT FOOT	REGULAR FOOT	LONG FOOT
Heel to End of Longest Toe	21-23	23-26	26-29

Standard Options

Standard options must be specified when placing an order.

FOOT			
	Closed toe	Open toe	Footless

GRIP TOP		3cm Strong Plain Below Knee Default		5cm Strong Plain Thigh High Default
-----------------	--	-------------------------------------	--	-------------------------------------

COLOUR					
	Beige	Black	Burgundy	Caramel	Mocha
	Navy	Pearl	Slate	Smoke	White

PATTERNS Class 1 & 2 only						
	None	Diamond	Honeycomb	Narrow Rib	Pinstripe	Zig-zag

Special Options

All our special options are fully documented in the **Hosiery Appendix Page 39** ►

Increase **panty compression**, change **waistband type**, change the **grip top style**, and browse reference **photos of colours and patterns**.





Compression Shorts

Material Composition: Polyamide, Elastane

Veni compression shorts are made from specially sized Veni RAL Class 2 fabric, to deliver a stiffer and higher compression than standard compression tights.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.


Versatile

Compressive Panty

For effective treatment of genital oedema

Graduated Leg Compression

Keeps legs comfortable



Compression

Class	mmHg	Stiffness	Recommended for
RAL Class 2	23-32	Low	Truncal and genital oedema, lipoedema affecting the hips and buttocks

Codes



Cycling Shorts




Capri Shorts

RAL Class 2	558	560
Unit of Sale	Single	Single

Sizing

Circumferences cm		XS	S	M	L	XL	XXL	XXXL
t	Waist	65-75	75-85	85-95	95-105	105-115	115-125	125-135
j	Buttocks	85-95	95-105	105-115	115-125	125-135	135-145	145-155
g	Thigh Top	43-51	47-55	49-56	55-63	59-67	63-71	67-74
f	Mid-Thigh	43-48	47-52	49-54	55-61	59-65	63-69	67-72
d	Below Knee	25-28	27-30	29-33	31-35	35-39	37-41	40-44
c	Widest Calf	27-30	29-32	33-37	36-40	39-43	42-46	45-51



Options

Colour

Beige

Black

Star Cotton

Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane, Cotton

100% cotton fabric in contact with skin — perfect for sensitive skin with superb comfort and support.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.

Comfort

High quality cotton
Soft and breathable for all day comfort

Opaque body yarns
Similar in appearance to active-wear or socks



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
 Class 1	18-21	Low	Subclinical Oedema. Mild watery oedema with sensitive skin conditions	Subclinical venous disorders / at risk of venous disorders with sensitive skin conditions.	Smooth skin with an increase of enlarged subcutaneous fat tissue.
 Class 2	23-32	Low	Swelling rarely reduces with elevation. Mild shape distortion with sensitive skin conditions.	Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Uneven skin with indentations in the fat tissue and lipomas.

Codes Available on **Drug Tariff**. See www.hadhealth.com/code-finder


											
	Anklet*	Below Knee	Thigh High	Thigh High + Waistband			Tights				
		+ GRIP TOP		+ GRIP TOP		LEFT LEG	BOTH	RIGHT LEG	LEFT LEG	BOTH	RIGHT LEG
 Class 1	034	341	340	345	349	357	353	361	367	365	369
 Class 2	074	741	740	742	743	747	746	748	750	749	751
Unit of Sale	Pair	Pair	Pair	Pair	Pair	Single	Single	Single	Single	Single	Single

Sizing

Extra wide sizes are shown in green

Circumferences cm	1	2	3	4	5	6	7	8
t Waist	65-75	70-80	75-85	80-90	85-95	90-105	95-110	100-120
j Buttocks	85-105	90-110	95-115	100-120	105-125	110-130	115-135	120-140
g Thigh Top	46-55 43-48	48-59 45-51	52-63 49-55	56-66 53-59	60-69 57-63	64-72 61-66	67-75 64-69	70-78 67-72
c Widest Calf	31-35 29-33	33-37 31-35	35-39 33-37	37-42 35-39	39-45 37-42	41-47 39-45	43-52 41-48	45-56 43-51
b Ankle	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-34

Measurements falling outside of the above chart?
See our **Custom Made** solutions on Page 36 ►



Length

Straight Distance cm

A-g Floor to Thigh Top

SHORT LEG

REGULAR LEG

LONG LEG

60-67

68-75

76-82

A-d Floor to Below Knee

33-37

38-42

43-46

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

Heel to End of Longest Toe

SHORT FOOT

REGULAR FOOT

LONG FOOT

21-23

23-26

26-29

Standard Options

Standard options must be specified when placing an order.


FOOT



GRIP TOP



COLOUR



Special Options

All our special options are fully documented in the **Hosiery Appendix Page 39** ►

Increase **panty compression**, change **waistband type**, change the **grip top style**.



Below Knee, Thigh High & Tights

Material Composition: Polyamide, Elastane

Robust circular knit for moderate lymphoedema that's soft and comfortable to wear, with the appearance of active-wear.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.

Durable

Thicker, opaque body yarns

Similar in appearance to active-wear or leggings

Stiffer circular knit

Suitable for patients needing stronger compression



Compression

Class	mmHg	Stiffness	Well suited for		
			LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
Class 2	23-32	Low	Moderate oedema which rarely reduces overnight with moderate pitting tissues.	Telangiectasia, thread, broken or reticular veins. Varicose veins, recurrent varicose veins.	Uneven skin with indentations in the fat tissue and lipomas.
Class 3	34-46	Low	Moderate oedema which rarely reduces overnight and is prone to rebound.	Healed venous ulcer.	N/A

Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder


											
	Anklet*	Below Knee		Thigh High		Thigh High + Waistband			Tights		
		+ GRIP TOP		+ GRIP TOP		LEFT LEG	BOTH	RIGHT LEG	LEFT LEG	BOTH	RIGHT LEG
 Class 2	070	701	700	702	703	707	706	708	717	709	718
 Class 3	072	721	720	722	723	727	726	728	737	729	738
Unit of Sale	Pair	Pair	Pair	Pair	Pair	Single	Single	Single	Single	Single	Single

Sizing

Extra wide sizes are shown in green

Circumferences cm	1	2	3	4	5	6	7	8
t Waist	65-75	70-80	75-85	80-90	85-95	90-105	95-110	100-120
j Buttocks	85-105	90-110	95-115	100-120	105-125	110-130	115-135	120-140
g Thigh Top	46-55 43-48	48-59 45-51	52-63 49-55	56-66 53-59	60-69 57-63	64-72 61-66	67-75 64-69	70-78 67-72
c Widest Calf	31-35 29-33	33-37 31-35	35-39 33-37	37-42 35-39	39-45 37-42	41-47 39-45	43-52 41-48	45-56 43-51
b Ankle	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-34

Measurements falling outside of the above chart?
See our **Custom Made** solutions on Page 36 ►



Length

Straight Distance cm

	SHORT LEG	REGULAR LEG	LONG LEG
A-g Floor to Thigh Top	60-67	68-75	76-82
A-d Floor to Below Knee	33-37	38-42	43-46

* Leg length options are unavailable for anklets.

Foot

Straight Distance cm

	SHORT FOOT	REGULAR FOOT	LONG FOOT
Heel to End of Longest Toe	21-23	23-26	26-29

Standard Options

Standard options must be specified when placing an order.

Foot			
	Closed toe	Open toe	Footless
Grip Top		3cm Strong Plain Below Knee Default	
			5cm Strong Plain Thigh High Default
Colour			
	Beige	Black	Navy
			
			Light Beige

Special Options

All our special options are fully documented in the **Hosiery Appendix Page 39** ►

Increase **panty compression**, switch to **male gusset**, change the **grip top style**, and browse reference **photos of colours and patterns**.



Socks

Material Composition: Polyamide, Elastane

Stride is a circular knit below knee garment knitted in the likeness of a sports sock, combining a pile knit sole with medical grade compression yarns up the leg.

Haddenham medical compression hosiery utilises high-quality, long-lasting yarns that deliver RAL standard graduated compression to the lower limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	VENOUS DISEASE
 Class 1	18-21	Low	Subclinical or Mild swelling with slight pitting and healthy skin.	Subclinical venous disorders /at risk of venous disorders with healthy skin.
 Class 2	23-32	Low	Moderate swelling with slight pitting and healthy skin.	Telangiectasia, thread, broken or reticular veins. Varicose veins, recurrent varicose veins.

Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder



Below Knee

 Class 1	820
 Class 2	821
Unit of Sale	Pair

Stock Options

All other variants
made-to-order in
8-10 days



Class

1-2

Sizes

XXS - XXL

Leg

Regular

Foot

Short
Regular
Long
Extra Long

Colour

Black

Active

Pile knit sole

Comfortable cushioning for people on their feet, or those keeping active.

Maximise movement

Helps maximise the benefits of movement to the lymphatic system.


Comfort

Four foot lengths for the perfect fit.



Sizing

Circumferences cm	XXS	XS	S	M	L	XL	XXL
c Widest Calf	29-33	31-35	33-37	35-39	37-42	39-45	41-48
b Ankle	18-20	20-22	22-24	24-26	26-28	28-30	30-32



Length

Straight Distance cm

A-d Floor to Below Knee

SHORT LEG

33-37

REGULAR LEG

38-43

Foot

Straight Distance cm

Heel to End of
Longest Toe

SHORT FOOT

21-23

REGULAR FOOT

23-26

LONG FOOT

26-29

EXTRA LONG

29-31

Options

Colour



Black



Grey



Navy



White

Microfine

Toecap

Material Composition: Polyamide, Elastane

Haddenham Microfine is an ultra-thin and smooth fabric that trims to fit without fraying.

Microfine Toecaps are an easy to measure, time-saving solution for oedema that extends to the toes.



Versatile

Time saving

No need to bandage toes

Trim to fit

Easily modified to desired toe length

Compression

mmHg

Well suited for

20-36


LYMPHOEDEMA

Positive Stemmer sign, fibrosis in the toes.
Mild to moderate shape distortion.

VENOUS DISEASE

In patients who present with oedema and are unable to wear closed toe hosiery. Prevention of toe oedema.

Sizing



Circumferences cm

XXS

XS

S

M

L

XL

a Ball of Foot

19-24

21-23

23-25

25-27

27-29

29-31

Codes

Available on **Drug Tariff**. See below or www.hadhealth.com/code-finder

	XXS	XS	S	M	L	XL
 Sand	218S-0 -	218S-1 347-9649	218S-2 347-9656	218S-3 347-9664	218S-4 347-9672	218S-5 386-4949
 Black	218B-0	218B-1 376-6284	218B-2 376-6292	218B-3 376-6300	218B-4 376-6318	218B-5 386-4956
Unit of Sale	Single	Single	Single	Single	Single	Single

Custom Made

Toecaps & Toemits

We aim to be flexible wherever possible by providing a made-to-measure service that offers unparalleled customisation and fit.



Bespoke

Made to Measure

Custom made toecaps are knitted to a patient's exact measurements, for comfort and conformity.

Download the [Custom Made Toecap](http://www.hadhealth.com/F0006) form
www.hadhealth.com/F0006 ►

Special Options



All our special options are fully documented in the [Hosiery Appendix Page 39](#) ►

Add a **heal strap**, or **customise toes**.

Codes

Available on [Drug Tariff](#). See www.hadhealth.com/code-finder



Toecaps



Toemits

	PERTEX LIGHT	PERTEX	GOLDPUNKT	MICROFINE	PERTEX LIGHT	PERTEX	GOLDPUNKT	MICROFINE
RAL Class 1	PX-TOE-1 PX-11	-	-	-	PX-TOEMIT-1 PX-11	-	-	-
RAL Class 2	-	PX-TOE-2 PX-12	GP-TOE-2 GP-12	-	-	PX-TOEMIT-2 PX-12	GP-TOEMIT-2 GP-12	-
RAL Class 3	-	PX-TOE-3 PX-13	GP-TOE-3 GP-13	-	-	PX-TOEMIT-3 PX-13	GP-TOEMIT-3 GP-13	-
20-36mmHg	-	-	-	MF-TOE	-	-	-	MF-TOEMIT
Unit of Sale	Single	Single	Single	Single	Single	Single	Single	Single

Custom Made

Medical Compression Hosiery

We aim to be flexible wherever possible by providing a made-to-measure service that offers unparalleled customisation and fit.

Bespoke

The widest choice

Select from our widest range of styles, fabric and options



Special Options



All our special options are fully documented in the **Hosiery Appendix Page 39** ►

Increase **panty compression**, change **waistband type**, add **fastenings** to aid donning, change the **grip top style**, or adjust the **heel and knee stitching** based on your patients activity levels.




Download the Custom Made Hosiery form
www.hadhealth.com/F0001 ►


Custom Made Flat Knit Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder

Pertex Light

					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 1	091CM	901CM PX-21 900CM PX-21 + CM90	902CM PX-31 903CM PX-31 + CM90	905CM PX-41 904CM PX-41 906CM PX-41	908CM PX-51 907CM PX-51 909CM PX-51
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Pertex


					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 2	092CM	921CM PX-22 920CM PX-22 + CM90	922CM PX-32 923CM PX-32 + CM90	925CM PX-42 924CM PX-42 926CM PX-42	928CM PX-52 927CM PX-52 929CM PX-52
RAL Class 3	093CM	931CM PX-23 930CM PX-23 + CM90	932CM PX-33 933CM PX-33 + CM90	935CM PX-43 934CM PX-43 936CM PX-43	938CM PX-53 937CM PX-53 939CM PX-53
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Goldpunkt


					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 2	094CM	941CM GP-22 940CM GP-22 + CM90	942CM GP-32 943CM GP-32 + CM90	945CM GP-42 944CM GP-42 946CM GP-42	948CM GP-52 947CM GP-52 949CM GP-52
RAL Class 3	096CM	961CM GP-23 960CM GP-23 + CM90	962CM GP-33 963CM GP-33 + CM90	965CM GP-43 964CM GP-43 966CM GP-43	968CM GP-53 967CM GP-53 969CM GP-53
RAL Class 4	098CM	981CM GP-24 980CM GP-24 + CM90	982CM GP-34 983CM GP-34 + CM90	985CM GP-44 984CM GP-44 986CM GP-44	988CM GP-54 987CM GP-54 989CM GP-54
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Custom Made Circular Knit Codes

Veni

					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 1	050CM	501CM 500CM	505CM 509CM	517CM 513CM 521CM	527CM 525CM 529CM
RAL Class 2	053CM	533CM 532CM	537CM 541CM	549CM 545CM 553CM	559CM 557CM 561CM
RAL Class 3	056CM	565CM 564CM	569CM 573CM	581CM 577CM 585CM	591CM 589CM 593CM
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Star Cotton

					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 1	034CM	341CM 340CM	345CM 349CM	357CM 353CM 361CM	367CM 365CM 369CM
RAL Class 2	074CM	741CM 740CM	742CM 743CM	747CM 746CM 748CM	750CM 749CM 751CM
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Doktus

					
	Anklet	Below Knee	Thigh High	Thigh High + Waistband	Tights
		+ GRIP TOP	+ GRIP TOP	LEFT LEG BOTH RIGHT LEG	LEFT LEG BOTH RIGHT LEG
RAL Class 2	070CM	701CM 700CM	702CM 703CM	707CM 706CM 708CM	717CM 709CM 718CM
RAL Class 3	072CM	721CM 720CM	722CM 723CM	727CM 726CM 728CM	737CM 729CM 738CM
Unit of Sale	Single	Single Single	Single Single	Single Single Single	Single Single Single

Appendix

Contents

Full information regarding sizing and the numerous special options available.

Unless stated otherwise all options available on **Drug Tariff**.

Measuring



Sizing


Page 40 ►

Special Options



Foot Options

Page 41 ►



Waistbands

Page 42 ►

Veni Exclusives



Patterns

Page 44 ►



Closures

Page 41 ►



Panty Options

Page 43 ►




Colours

Page 44 ►



Grip Tops

Page 42 ►



Other Options

Page 43 ►

Special Option Codes

Page 45 ►

Sizing



Foot Length

		SHORT FOOT	REGULAR FOOT	LONG FOOT
HEEL TO TOP OF LONGEST TOE	Centimetres	21-23	23-26	26-29
	UK Shoe Size	2-5	5-9	9-11
	EU Shoe Size	35-38	39-42	43-46

Reference sizing for Open Toe

HEEL TO BALL OF FOOT AT BIG TOE	Centimetres	13-15cm	15-18cm	18-21cm
---------------------------------	-------------	---------	---------	---------

Foot Options



Open Toe

Open toe finish on the foot.



Closed Toe

Closed toe finish on the foot.



Slant Toe

Open toe slant option. Please include inside and outside foot measurements.

Flat Knit Only
Custom Made Only



Footless

You can choose to finish your footless tights at any length. Keep in mind that without anchorage, footless garments might ride up so choose your finishing point with consideration of the leg shaping.



y-Heel

A y-heel is the **default fit**, offering a soft, contoured shape that accommodates natural movement around the ankle and heel.




t-Heel

A t-heel design redirects fabric tension across the front of the ankle to reduce pressure and irritation during upward foot movement while ensuring a secure, contoured fit around the heel.

Flat Knit Only

Closures



Zips

Zips can be added in any location but remember that this will alter the fabric properties at that location. For example, it won't stretch longitudinally where the zip is sewn in.

Flat Knit Only



Velcro Fastenings

The Velcro runs vertically down the length of the garment to oppose the two sides and close the garment.

Flat Knit Only



Velcro Straps

Velcro straps can be added in any location along the length of the leg. This however does not allow the garment to perform the same function as a wrapping device. The two edges need to be carefully aligned to prevent puckering.

Flat Knit Only

Grip Tops



Strong Plain

3cm / 5cm



Strong Lace

5cm




Fine Lace

5cm



Non-Silicone

4cm




Strong Plain (Veni)

3cm / 5cm

Microdot design, only available on Veni garments.

Waistbands



Drawstring Waist

The adjustable drawstring waist is a good option if you need to accommodate wider hips and waist.



Elastic Waistband


Flat knit default, non-silicon elastic, 4cm wide waistband is made to fit the waist circumference. This option may not be suitable for wider hips and waist.



Knitted Waistband

Circular knit default, knitted ribbed waist designed to stretch.
Circular Knit Only

Panty Finishing



No Gusset

Open pubis.



Fly/Male Gusset

Closed pubis with horizontal or vertical fly at front.

Larger gusset to increase comfort for male genitals.



Loose Fit

Fully panty but made more generously with two seams to accommodate larger abdomens or to allow for extra comfort if someone is intolerant of compression over the abdomen.

Panty Compression



No Compression

Unless another option is specified panty section will be made without compression.



Half Compression

The panty section will be knitted with half compression. That is, every second row will have a compression inlay threads.



Full Compression

The garments will continue to be made with graduated compression all the way to the waist with the same thread and tension as what you have asked for in the legs.

Other Options



Pocket/Padding/Lining

To provide additional protection to delicate tissues.

Custom Made Only

Not Available on Drug Tariff



Malleolus Pad

Padding at either side of the ankle for additional support and comfort.

Custom Made Only




Knee Shaping Flat Knit Only

No knee shaping is default.

90° knee reduced bunching when the leg remains in a fixed position, such as those who are wheelchair bound.

60° knee offering reduced bunching for those who are active but have day to day roles where sitting a lot.



Slant Top

Flat Knit Only

Slant finish to accommodate more garment coverage on the lateral side of the limb. Please include inside and outside measurements.

Custom Made Only

Veni Patterns

Available in Class 1 and Class 2 only.



Honeycomb

Zig Zag

Pin Stripe

Narrow Rib

Diamond

Veni Colours



Beige

Black

Pearl

Mocha*

*Class 1 and 2 only

Burgundy*

*Plain Grip Top Only
*Class 1 and 2 only



Slate

White

Smoke

Navy


Special Option Codes

Available on **Drug Tariff** *. See www.hadhealth.com/code-finder

Category	Description	Sized		Custom Made	
		Code		Code	
Foot	Footless	NF *		NF-CM *	
	T-Heel	TH		TH-CM *	
Padding	Partial Lining	SPECIAL		LIN-CM	
	Pocket	SPECIAL		POC-CM	
	Malleolus Pads	SPECIAL		MAL-CM CM-91	
	Pad	SPECIAL		PAD-CM	
	Non-Standard/Additional	SPECIAL		GT-CM CM-90	
Grip Top	Zip	SPECIAL		ZIP-CM CM-92	
	Velcro Fastening	SPECIAL		VELF-CM CM-93	
Closures	Velcro Strap	SPECIAL		VELS-CM CM-94	
	Slant Foot	SF		SF-CM *	
Finishing	Slant Top	ST		ST-CM *	
	Heel Strap	HEEL		HEEL-CM	
Tights	No Gusset	NG *		NG-CM *	
	Fly Opening	FY *		FY-CM *	
	Open Front	OF *		OF-CM *	
	Loose Fit	LO *		LO-CM *	
	Male Gusset	MG *		MG-CM *	
	Elastic Waistband	EWB		-	
	Half Compression Panty	HC *		HC-CM *	
	Full Compression Panty	FC *		FC-CM *	

Compression

Residual pressure refers to the gradually decreasing pressure above the ankle, expressed as a percentage of the ankle pressure, and varies by compression class.



Compression Class	mmHg	Residual Pressure Ratio		
		b¹ Start of calf muscle	c Mid-calf	f / g Mid-thigh / Thigh-top
Class 1	18-21	70-100%	50-80%	20-60%
Class 2	23-32	70-100%	50-80%	20-50%
Class 3	34-46	70-100%	50-70%	20-40%
Class 4	49+	70-100%	50-70%	20-40%



Armsleeves



Medical Compression Armsleeves

Haddenham medical compression armsleeves are made from high-quality, long-lasting yarns that deliver RAL standard graduated compression.



What makes it so patient-friendly?

Turn the page to discover all features [Page 48](#) ►

Sized Flat Knits



Pertex Light

Page 52 ►

Sized Circular Knits



Venex

Page 58 ►



Star Cotton

Page 60 ►

Custom Made



Armsleeves

Page 65 ►



Pertex

Page 54 ►



Doktus

Page 62 ►



Microfine Glove

Page 64 ►



Goldpunkt

Page 56 ►

Haddenham medical compression armsleeves are made from high-quality, long-lasting yarns that deliver RAL standard graduated compression.

Freedom to choose



Pertex

6 versatile fabrics


Flat and Circular knit garments suitable for treatment of all stages of Lymphoedema, Venous disorders and Lipoedema.

[Fabric selection guide](#)
See page 50 ►

Reliable

Sleeves that can keep up with the most active of lifestyles, thanks to conforming fabrics and a variety of grip top options.





Venex

Graduated compression

Works with the body's natural circulation - concentrating the highest compression at the wrist and gradually decreasing up the limb.

Comfort

Elbow stretch zone and gentle thumb gusset for freedom in motion and maintaining dexterity.

Wide range of sizing to find the perfect fit.

Sized

Our sized sleeves are **made-to-order**, with the care and quality that we believe each patient deserves.

We also retain **stock** of frequently ordered variants for an even faster turnaround.

_____ or _____

Custom Made

We aim to be flexible wherever possible by providing a made-to-measure service that offers unparalleled customisation and fit.



Microfine
Glove

Which fabric?

Please use the below maintenance phase pathway as a guide to which of our fabrics would most suitable for your patient based on their presentation and symptoms.

How would you describe the level of your patient's oedema?

Early accumulation of fluid which subsides with elevation, with or without pitting.
ISL Stage 0-1

Mild

Tissue swelling rarely reduces with elevation and pitting is present.
ISL Stage 2

Moderate

Pitting is absent with skin thickening fibrosis and excess subcutaneous fat.
ISL Late Stage 2 - Stage 3

Severe

What is the level of shape distortion? ●

MILD

EXTENSIVE

MODERATE

Is the oedema soft and watery?

YES

NO

How would you describe the texture of the patient's tissue? ◆

WOODY

PUTTY LIKE

What is the condition of the patient's skin? ▲

▲

What is the level of pitting? ●

●

Are warty overgrowths with fat lobules present?

NO

YES

HEALTHY SKIN

SENSITIVE

FRAGILE

SLIGHT

MODERATE

DEEP

Venex

STAR
Cotton

Pertex
Light

Dokus

Pertex

GOLDPUNKT

References

▲ Skin conditions

Fragile Thin skin prone to skin tears

Sensitive the presence of skin conditions where the skin is prone to reactions such as redness and itching.

◆ Textures

Watery soft pliable tissues.

Putty putty consistency represents early fibrotic changes, will pit with firm pressure.

Woody firm woody texture is consistent with advanced fibrotic tissue changes and does not pit.

● Pitting Levels


Slight 1+ 2mm barely detectable immediate rebound

Moderate 2+ 4mm deep <10 second rebound

Deep 3+ 6 mm deep <30 second rebound

Very deep 4+ up to 8mm deep >30 second rebound

● Shape distortion



Mild

Moderate

Extensive

How much compression?

Graduated compression is highest at the wrist and gradually decreases as the yarns travel up the limb.

When deciding how much compression to apply also consider these 3 factors:

1. The nature of the fabric yarns and the way in which they are knitted together.
2. The fabric's ability to retain its structure under pressure. Measured as the increase in compression per centimetre increase in the circumference of the limb (mmHg/cm).
3. The resultant force exerted on the limb measured in mmHg.



Stiffness mmHg /cm

✓ Low = 0-2 | ✓ Medium = 2-4 | ✓ High = 4-6

	Circular Knit			Flat Knit		
	Venex	STAR Cotton	Dokus	Pertex Light	Pertex	GOLDPUNKT
Pressure						
Class 1 18 - 21mmHg	✓ 1.2	✓ 1.4		✓ 1.5		
Class 2 23-32mmHg	✓ 1.3	✓ 1.5	✓ 1.6		✓ 2.1	✓ 2.8
Class 3 34-46mmHg			✓ 1.8		✓ 4.3	✓ 3.3

Page 58 ►

Page 60 ►

Page 62 ►

Page 52 ►

Page 54 ►

Page 56 ►

Pertex Light

Armsleeve, Glove and Mitten

Material Composition: Polyamide, Elastane

A unique stretchy flat knit that stretches with ease and hugs gently—ideal for fragile or sensitive skin.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for
			LYMPHOEDEMA
 Class 1	18-21	Low	Mild soft watery oedema with slightly pitting tissues and fragile skin. Mild shape distortion. Palliative care.
			LIPOEDEMA
			Smooth skin with an increase of enlarged subcutaneous fat tissue.

Codes

Available on **Drug Tariff** as Custom Made. See www.hadhealth.com/code-finder

					
Mitten		Glove		Sleeve	
		LEFT	RIGHT	+ GRIP TOP	
 Class 1	106	105L	105R	143	145
Unit of Sale	Single	Single	Single	Single	Single

Effortless

Uniquely stretchy

A stretchy knit makes for effortless donning and doffing, which means better patient compliance and self management.




Sizing

EXTRA WIDE shown in green

Circumferences cm	XS	S	M	L
g Axilla	26-32 20-26	28-35 21-28	34-48 24-34	40-50 30-40
e Elbow Fold	18-22	19-24	22-28	26-32
c Wrist	14-16	16-19	20-24	22-26
a Base of Fingers	17-18	18-20	20-22	22-24

Length cm	REGULAR	LONG
c-g Wrist to Axilla	40-44	44-50



Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black



Light Beige*

*Light Beige
not available in
Class 3

Shoulder Cap

Please state g-h measurement
and left/right arm



With Velcro® Bra Attachment
£20.00

Please state g-h measurement,
g-t circumference and left/right arm



With Adjustable Strap
£25.00

Pertex

Armsleeve, Glove and Mitten

Material Composition: Polyamide, Elastane

Firm knit with soft touch making it a go-to option for reliable medium compression needs.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	LIPOEDEMA
 Class 2	23-32	Medium	Moderate oedema with moderate shape distortion and deep pitting tissues. Post liposuction for lymphoedema.	Uneven skin with indentations in the fat tissue and lipomas.
 Class 3	34-46	High	Obesity-related lymphoedema or lymphoedema with extensive shape distortion, fat depositions, fibrosis and warty overgrowths, indicating elephantiasis.	N/A

Codes

Available on **Drug Tariff** as Custom Made. See www.hadhealth.com/code-finder

Mitten		Glove		Sleeve		Sleeve with Mitten	
		LEFT	RIGHT		+ GRIP TOP		+ GRIP TOP
Class 2	206	205L	205R	333	203	234	204
Class 3	306	305L	305R	331	311	332	312
Unit of Sale	Single	Single	Single	Single	Single	Single	Single

Conforming

Firm knitting


Firm consistent hold for moderate oedema



Sizing

EXTRA WIDE shown in green

Circumferences cm	XS	S	M	L
g Axilla	33-36 28-32	36-40 32-36	40-44 36-40	44-48 40-44
e Elbow Fold	23-25	25-27	28-30	30-32
d Mid Forearm	23-24	24-26	26-29	29-32
c Wrist	17-18	18-20	20-22	22-24
a Base of Fingers	18-19	19-20	21-23	23-25



Length cm	REGULAR	LONG
c-g Wrist to Axilla	40-44	44-50

Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black



Light Beige*

*Light Beige
not available in
Class 3

Shoulder Cap

Please state g-h measurement
and left/right arm



With Velcro® Bra Attachment
£20.00

Please state g-h measurement,
g-t circumference and left/right arm



With Adjustable Strap
£25.00

Armsleeve, Glove and Mitten

Material Composition: Polyamide, Elastane

Made for stubborn and rebound oedema support, with silky low friction on the inside to make it easy to apply.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	LIPOEDEMA
Class 2	23-32	Medium	Obesity-related lymphoedema or lymphoedema with moderate shape distortion and very deep pitting. Post liposuction for lymphoedema.	Uneven skin with indentations in the fat tissue and lipomas.
Class 3	34-46	High	Obesity-related lymphoedema or lymphoedema where there is excess subcutaneous fibrosis and fat.	N/A

Codes

Available on **Drug Tariff** as Custom Made. See www.hadhealth.com/code-finder



Mitten



Glove



Sleeve



Sleeve with Mitten

		LEFT	RIGHT		+ GRIP TOP		+ GRIP TOP
 Class 2	209	207L	207R	315	213	316	214
 Class 3	309	307L	307R	317	313	318	314
Unit of Sale	Single	Single	Single	Single	Single	Single	Single

Tough

High Compression


Effective treatment for the most stubborn oedema



Sizing

EXTRA WIDE shown in green

Circumferences cm	XS	S	M	L
g Axilla	33-36 28-32	36-40 32-36	40-44 36-40	44-48 40-44
e Elbow Fold	23-25	25-27	28-30	30-32
d Mid Forearm	23-24	24-26	26-29	29-32
c Wrist	17-18	18-20	20-22	22-24
a Base of Fingers	18-19	19-20	21-23	23-25



Length cm

REGULAR

LONG

c-g Wrist to Axilla	40-44	44-50
----------------------------	-------	-------

Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black

Shoulder Cap

Please state g-h measurement
and left/right arm



With Velcro® Bra Attachment

Please state g-h measurement,
g-t circumference and left/right arm



With Adjustable Strap

Material Composition: Polyamide, Elastane

A semi-sheer, lightweight and breathable fabric with microfibre softness.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	LIPOEDEMA
 Class 1	18-21	Low	Mild soft watery oedema with healthy skin.	Smooth skin with an increase of enlarged subcutaneous fat tissue.
 Class 2	23-32	Low	Mild soft watery oedema with slightly pitting tissues, mild shape distortion and healthy skin.	Uneven skin with indentations in the fat tissue and lipomas.

Codes

Available on Drug Tariff. See www.hadhealth.com/code-finder



Sleeve



Sleeve with Mitten

		+ GRIP TOP		+ GRIP TOP
 Class 1	103	101	104	102
 Class 2	223	221	224	222
Unit of Sale	Single	Single	Single	Single

Stock Options

All other variants
made-to-order in
8-10 days



Class

1-2

Grip Top

3cm Plain

Sizes

1-8

Widths

Standard

Length

Regular

Colour

Beige

Black

Gentle

Soft and stretchy
for easy application

Conforms
to limb shape for comfort




Sizing

EXTRA WIDE shown in green

Circumferences cm	XXS	XS	S	M	L
g Axilla	29-33 26-28	33-36 28-32	36-40 32-36	40-44 36-40	44-48 40-44
e Elbow Fold	21-23	23-25	25-27	28-30	30-32
c¹ 4cm Above Wrist	14-15	15-17	17-19.5	19.5-22.5	22.5-26
a Base of Fingers	13-16	14-18	16-21	18-24	21-27

Length cm	REGULAR	LONG
c-g Wrist to Axilla	38-44	44-50



Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black



Caramel



Pearl



Slate



Smoke

Shoulder Cap

Please state g-h measurement
and left/right arm



With Velcro® Bra Attachment

Please state g-h measurement,
g-t circumference and left/right arm



With Adjustable Strap

Star Cotton

Armsleeve

Material Composition: Cotton, Polyamide, Elastane

Star Cotton is a circular knit fabric comprised of cotton covered yarns, which makes it an ideal choice for those with an intolerance to non-natural fibres.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	LIPOEDEMA
 Class 1	18-21	Low	Mild soft watery oedema with sensitive skin.	Smooth skin with an increase of enlarged subcutaneous fat tissue.
 Class 2	23-32	Low	Mild soft watery oedema with slightly pitting tissues, mild shape distortion and sensitive skin.	Uneven skin with indentations in the fat tissue and lipomas.

Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder



Sleeve

+ GRIP TOP



Sleeve with Mitten

+ GRIP TOP

 Class 1	161	163	162	164
 Class 2	261	263	262	264
Unit of Sale	Single	Single	Single	Single

Comfort

100% Cotton against skin

Minimises skin irritation for those with an intolerance to non-natural fibres.

Soft wide cuff at wrist

To increase comfort and reduce irritation




Sizing

EXTRA WIDE shown in green

Circumferences cm	XS	S	M	L
g Axilla	33-36 28-32	36-40 32-36	40-44 36-40	44-48 40-44
e Elbow Fold	23-25	25-27	28-30	30-32
c¹ 4cm Above Wrist	15-17	17-19.5	19.5-22.5	22.5-26
a Base of Fingers	14-18	16-21	18-24	21-27

Length cm	REGULAR	LONG
c-g Wrist to Axilla	40-44	44-50



Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black

Shoulder Cap

Please state g-h measurement and left/right arm



With Velcro® Bra Attachment

Please state g-h measurement, g-t circumference and left/right arm



With Adjustable Strap

Doktus

Armsleeve

Material Composition: Polyamide, Elastane

Doktus is a circular knit fabric with thick, opaque body yarns.

Haddenham medical compression armsleeves utilise high-quality, long-lasting yarns that deliver RAL standard graduated compression to the upper limb.



Compression

Class	mmHg	Stiffness	Well suited for	
			LYMPHOEDEMA	LIPOEDEMA
 Class 2	23-32	Medium	Moderate oedema which rarely reduces overnight with moderate pitting tissues and shape distortions.	Uneven skin with indentations in the fat tissue and lipomas.
 Class 3	34-46	High	Moderate oedema which rarely reduces overnight and is prone to rebound.	N/A

Codes

Available on Drug Tariff. See www.hadhealth.com/code-finder



Sleeve



Sleeve with Mitten

		+ GRIP TOP		+ GRIP TOP
 Class 2	215	201	217	202
 Class 3	216	301	219	302
Unit of Sale	Single	Single	Single	Single

Strength

Stiff Compression

Tackles stubborn rebound oedema

Sturdy body yarns


Distribute pressure and contain the oedema

Sizing

EXTRA WIDE shown in green

Circumferences cm	XS	S	M	L
g Axilla	33-36 28-32	36-40 32-36	40-44 36-40	44-48 40-44
e Elbow Fold	23-25	25-27	28-30	30-32
c¹ 4cm Above Wrist	15-17	17-19.5	19.5-22.5	22.5-26
a Base of Fingers	14-18	16-21	18-24	21-27

Length cm	REGULAR	LONG
c-g Wrist to Axilla	40-44	44-50



Options

Grip Top



3cm
Strong Plain



5cm
Strong Plain



5cm
Fine Lace



5cm
Strong Lace

Colour



Beige



Black



Navy



Light Beige

Shoulder Cap

Please state g-h measurement
and left/right arm



With Velcro® Bra Attachment

Please state g-h measurement,
g-t circumference and left/right arm



With Adjustable Strap

Microfine

Glove

Material Composition: Polyamide, Elastane

Microfine Gloves are an easy to measure, time-saving solution for oedema that extends to the fingers.

Haddenham Microfine is an ultra-thin and smooth fabric that trims to fit without fraying.

Finger Saving

No need to bandage fingers
Time saving and finger saving



Compression

Class	mmHg	Well suited to
-	20-36	Positive Stemmer sign. Mild to moderate shape distortion and slight to moderate tissue pitting

Sizing



Circumferences cm	XXS	XS	S	M	L
a Base of Fingers	14-16	16-18	18-20	20-22	22-24
b Palm at Fold of Thumb	15-17	17-19	19-21	21-23	24-26

Codes Available on Drug Tariff

	XXS	XS	S	M	L
Sand	208S-0	208S-1 347-9631	208S-2 347-9607	208S-3 347-9615	208S-4 347-9623
Black	208B-0	208B-1 376-6219	208B-2 376-6227	208B-3 376-6243	208B-4 376-6268
Unit of Sale	Single	Single	Single	Single	Single

Custom Made for Upper Limb

We aim to be flexible wherever possible by providing a made-to-measure service that offers unparalleled customisation and fit.

Bespoke

The widest choice

Select from our widest range of styles, fabric and options

Download the [Custom Made Armsleeves](http://www.hadhealth.com/F0002) form
www.hadhealth.com/F0002 ►



Special Options

Available on [Drug Tariff](#) *. See www.hadhealth.com/code-finder

Padding



Partial Lining

UL-LIN-CM



Pocket

UL-POC-CM

Closures



Zips

UL-ZIP-CM
CM-81



Velcro Fastenings

UL-VELF-CM



Velcro Straps

UL-VELS-CM

Finishing



Slant Top

UL-ST-CM

Shoulder Cap



Shoulder Cap & Velcro® Bra Attachment

UL-SCBA-CM
CM-83



Shoulder Cap & Adjustable Strap

UL-SCAS-CM
CM-83

Custom Made Codes

Available on **Drug Tariff**. See www.hadhealth.com/code-finder

Pertex Light

											
	Mitten	Mitten	Glove	Sleeve	Sleeve with Mitten	Sleeve with Mitten	Sleeve with Glove				
	No Thumb					No Thumb					
				+ GRIP TOP	+ GRIP TOP	+ GRIP TOP	+ GRIP TOP				
RAL Class 1	150CM	106CM	105CM	143CM	145CM	180CM	181CM	182CM	183CM	184CM	185CM
	PX-61	PX-61	PX-64	PX-71	+CM-80	PX-74	+CM-80	PX-74	+CM-80	PX-77	+CM-80
Unit of Sale	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single

Pertex

											
	Mitten No Thumb	Mitten	Glove	Sleeve	Sleeve with Mitten	Sleeve with Mitten No Thumb	Sleeve with Glove				
				+ GRIP TOP	+ GRIP TOP	+ GRIP TOP	+ GRIP TOP				
RAL Class 2	250CM PX-62	206CM PX-62	205CM PX-65	333CM PX-72	203CM +CM-80	234CM PX-75	204CM +CM-80	282CM PX-75	283CM +CM-80	284CM PX-78	285CM +CM-80
RAL Class 3	350CM PX-63	306CM PX-63	305CM PX-66	331CM PX-73	311CM +CM-80	332CM PX-76	312CM +CM-80	382CM PX-76	383CM +CM-80	384CM PX-79	385CM +CM-80
Unit of Sale	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single

Goldpunkt

											
	Mitten No Thumb	Mitten	Glove	Sleeve		Sleeve with Mitten		Sleeve with Mitten No Thumb		Sleeve with Glove	
				+ GRIP TOP		+ GRIP TOP		+ GRIP TOP		+ GRIP TOP	
 Class 2	251CM GP-62	209CM GP-62	207CM GP-65	315CM GP-72	213CM +CM-80	316CM GP-75	214CM +CM-80	286CM GP-75	287CM +CM-80	288CM GP-78	289CM +CM-80
 Class 3	351CM GP-63	309CM GP-63	307CM GP-66	317CM GP-73	313CM +CM-80	318CM GP-76	314CM +CM-80	386CM GP-76	387CM +CM-80	388CM GP-79	389CM +CM-80
Unit of Sale	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single	Single

Custom Made Codes

Venex

	 Mitten No Thumb	 Sleeve	 Sleeve with Mitten No Thumb
		+ GRIP TOP	+ GRIP TOP
RAL Class 1	110CM	103CM	101CM
RAL Class 2	210CM	223CM	221CM
Unit of Sale	Single	Single	Single

Star Cotton

	 Mitten No Thumb	 Sleeve	 Sleeve with Mitten No Thumb
		+ GRIP TOP	+ GRIP TOP
RAL Class 1	170CM	161CM	163CM
RAL Class 2	270CM	261CM	263CM
Unit of Sale	Single	Single	Single

Doktus

	 Mitten No Thumb	 Sleeve	 Sleeve with Mitten No Thumb
		+ GRIP TOP	+ GRIP TOP
RAL Class 2	230CM	215CM	201CM
RAL Class 3	330CM	216CM	301CM
Unit of Sale	Single	Single	Single



easywrap





An adjustable compression wrap system comprised of single-layer knitted fabric bands which are fastened in an overlapping pattern.



What makes it so patient-friendly?

Turn the page to discover all features [Page 70 ▶](#)



Leg
[Page 74 ▶](#)



Foot
[Page 76 ▶](#)



Thigh
[Page 78 ▶](#)



Knee
[Page 80 ▶](#)



Fusion
[Page 82 ▶](#)



Liners
[Page 84 ▶](#)



Extenders
[Page 85 ▶](#)




Hand
[Page 86 ▶](#)



Arm
[Page 88 ▶](#)



Additional Bands
[Page 90 ▶](#)



Custom Made
[Page 91 ▶](#)

Haddenham easywrap is a single layer, compression wrap system comprised of overlapping bands which mimic the therapeutic effect of short-stretch bandages.

A cost-effective alternative to bandaging during the decongestive phase of treatment, when tackling stubborn or rebound oedema, as well as for the treatment of venous leg ulcers.

Compression made easy





Reduce

High working pressures mean that the more active you are, the better your lymphatic system will work to reduce oedema.

Maintain

Low resting pressures and high static-stiffness maintain limb reduction when you are less active.

Two classes of compression



*Compression values delivered within $\pm 15\%$ + 1mmHg

Low profile

Can be worn comfortably under most everyday clothing and footwear.





Easy to apply

A distinctive “lock-out” or end-stretch makes them easy to apply at the correct tension.



Position



Stretch



Fasten

Liners

easywrap for lower limb comes with a pair of protective liners - one to wear and one to wash.

For more details on our liner policy see Liners Page 84 ►

Effective

“Patients who have previously been non concordant with treatment, whether that be due to the psychosocial aspects of bandages or the inability to apply compression garments have all experienced and increase in quality of life when using easywrap.”

- British Journal of Community Nursing



Figure 4. Lymphorrhoea and venous in chronic bilateral lower-limb oedema

Case with Specialist Interest in Chronic and Complex Wounds, Hywel Dda University Health Board

William, presented with two wounds on his lower limbs, but he reported no lymph fluid leakage or

(ABPI) as 1.2, with biphasic waveforms bilaterally. He had no ischaemic rest pain, intermittent claudication, a capillary refill of less than 3 seconds, and a limb that was warm to touch. These findings indicated that William had venous insufficiency, and this was the likely cause of his oedema, lymphorrhoea and ulceration.

William reported a great deal of discomfort in his legs as a result of the high volume of lymphorrhoea causing excoriation and burning. When asked to score his pain using a numeric rating scale from 0 (no pain) to 3 (worst pain ever), he reported that it was at least 1 all the time, on a scale of 0 to 3 if he had

applying compression therapy was first discussed with the heart-failure nurse. The heart-failure nurse agreed that, as William's symptoms were stable, he could undertake compression therapy with close monitoring. Therefore, the benefits of compression therapy were communicated to William, who agreed to trying it.

However, William was reluctant to use compression bandages that would need to be regularly changed in clinic, when his work commitments meant he could not attend clinic appointments more than once a week. He was very keen for an alternative compression system that



Figure 2. Week 6: left anterior (a) and posterior (b) lower leg

periwound skin, and an antimicrobial cream with a non-adherent dressing to protect and prevent trauma to the skin affected. He also received gloves and an information leaflet on leg care during dressing changes. William and his wife were shown how to apply and remove easywrap and the dressings and given advice on hand and leg-washing hygiene, such as cleaning the shower after use and washing towels after they had been used to dry the leg. William was delighted that he would be able to shower and wash his leg independently, reducing any exudate present on the skin,




Figure 3. Week 12: left anterior (a) and posterior (b) lower leg

had reduced by 1 cm. His pain score had reduced, occasionally reaching 1, such as when he had been walking up stairs or driving. Otherwise, William felt that his quality of life had improved significantly, because he was able to mobilise better

with the pain changed every 2 or 3 days. William's left leg was less painful, more often than not scoring 0 and only occasionally scoring 1. The anterior wound had reduced in size to 3 cm long and 1 cm wide and comprised 100% granulation tissue (Figure 3). The posterior wound had healed. The periwound skin was even softer to touch. His calf and ankle circumferences had reduced by 3 cm and 2 cm respectively in total over the 12 weeks.

Discussion

William felt much more optimistic since using easywrap Strong. He found that it was comfortable to wear, did not increase his pain levels and allowed him to see positive results, such as exudate reduction and pain reduction, fairly quickly. Despite initially worrying that it would be complex, his wife had found easywrap to be easy to apply and remove. They both felt that it was a better alternative to bandaging for them,



Mobility

Featuring a patented knee section which maintains freedom of movement and enhances everyday comfort.

“This is very easy to put on, comfortable and really helpful for my lymphoedema.”

Patient review, lymphshop.com

Time-saving

A simple and cost-effective alternative to bandaging that promotes self-care and independence, and reduces waste.

Guaranteed

6 months of use.

Material composition

Polyester, Polyurethane

Availability

All variants held in stock.



easywrap

Leg

Material Composition: polyamide, elastane

The easywrap compression wrap system for lower limb comes in three parts, which work in isolation or combination.

easywrap Leg is comprised of short-stretch bands. Each band overlaps by 50% to ensure therapeutic, graduated compression and can be adjusted as swelling fluctuates.

“The most comfy wrap I have ever worn.”


Patient review, lymphshop.com



Compression

Class	mmHg	Well suited for		
		LYMPHATIC	VENOUS	LIPOEDEMA
Light	20-30	Mild oedema, Fragile skin, palliative care	Mixed aetiology ulcers. Telangiectasia, thread,spider or broken veins. Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.
Strong	30-40	Moderate oedema with moderate shape distortion and deep pitting tissues. Post liposuction for lymphoedema.	Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica. Healed Venous Ulcer.	Large deformations of fat, lipolymphoedema, large fat nodules and lipomas.

Sizing



Circumferences cm		XS	S	M	L	XL
c	Widest point of calf	38-45	44-52	50-60	55-65	60-70
b	Ankle	23-27	28-32	33-38	39-44	45-52
Contoured Lengths cm		REGULAR		TALL		
b-d	Ankle to below knee	30-34		34-38		

Easy to wear

Low profile

Can be worn under everyday clothing, even jeans

Velcro® Fasteners

For a strong hold and an effortless release

Adaptive

Highest static stiffness in class*

Stops oedema rebounding after reduction

High working pressure

for improved venous return



Codes

Available on **Drug Tariff**

*Phillips, N. and Wright, T (2024)
'Adjustable Compression Wraps, Background,
evidence and case studies on use in lymphatic
and venous disease' *Journal of Woundcare*, 15 July.

Light 20-30mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	REGULAR	EL-LR1 400-4917	EL-LR2 400-4925	EL-LR3 400-4933	EL-LR4 400-4941	EL-LR5 400-4958
	TALL	EL-LT1 400-4966	EL-LT2 400-4974	EL-LT3 400-4982	EL-LT4 400-4990	EL-LT5 400-5005
 BLACK	REGULAR	ELB-LR1 408-9397	ELB-LR2 408-9405	ELB-LR3 408-9413	ELB-LR4 408-9421	ELB-LR5 408-9439
	TALL	ELB-LT1 408-9355	ELB-LT2 408-9389	ELB-LT3 408-9363	ELB-LT4 408-9348	ELB-LT5 408-9314

Strong 30-40mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	REGULAR	ES-LR1 400-5518	ES-LR2 400-5526	ES-LR3 400-5534	ES-LR4 400-5542	ES-LR5 400-5559
	TALL	ES-LT1 400-5567	ES-LT2 400-5575	ES-LT3 400-5583	ES-LT4 400-5591	ES-LT5 400-5609
 BLACK	REGULAR	ESB-LR1 408-9009	ESB-LR2 408-8985	ESB-LR3 408-8993	ESB-LR4 408-9025	ESB-LR5 408-9041
	TALL	ESB-LT1 408-8951	ESB-LT2 408-8969	ESB-LT3 408-8944	ESB-LT4 408-8977	ESB-LT5 408-8928

easywrap

Foot

Material Composition: polyamide, elastane

The easywrap compression wrap system for lower limb comes in three parts, which work in isolation or combination, and are designed to be worn over a liner.

The easywrap Foot consists of a set of short-stretch bands extending from a central spine at the back of the garment.

“This is very easy to put on, comfortable and really helpful for my lymphoedema.”

Patient review, lymphshop.com




Compression

Class	mmHg	Well suited for		
		LYMPHATIC	VENOUS	LIPOEDEMA
Light	20-30	Mild oedema of the foot, Fragile skin, palliative care.	Mixed aetiology ulcers. Telangiectasia, thread,spider or broken veins. Recurrent/ Varicose veins. Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Feet are usually spared.
Strong	30-40	Moderate oedema of the foot with moderate shape distortion and deep pitting tissues. Post liposuction for lymphoedema.	Recurrent/Active venous Leg Ulcer. Healed venous ulcer.	Lipolymphoedema, large fat nodules and lipomas.

Sizing

Circumferences cm		XS	S	M	L	XL
a ¹	Mid-foot	22-24	25-27	28-30	31-34	35-40
Straight Distance cm						
a-heel Base of toes to heel	REGULAR	14.5-16.5	16-18	17.5-19.5	19-21	21-23
	LONG	16.5-18.5	18-20	19.5-21.5	21-23	23-25



Easy to wear

Low profile

Fits into most shoes

Single posterior layer

For maximum comfort
and smooth contouring

Versatile

The perfect partner

Pair with an easywrap Leg and/
or a Microfine Toecap for graduated
compression from toes to knee.



Codes

Available on **Drug Tariff**

Light 20-30mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	REGULAR	EL-FR1 400-4867	EL-FR2 400-4875	EL-FR3 400-4883	EL-FR4 400-4891	EL-FR5 400-4909
	LONG	EL-FL1 400-4818	EL-FL2 400-4826	EL-FL3 400-4834	EL-FL4 400-4842	EL-FL5 400-4859
 BLACK	REGULAR	ELB-FR1 408-9553	ELB-FR2 408-9546	ELB-FR3 408-9538	ELB-FR4 408-9520	ELB-FR5 408-9512
	LONG	ELB-FL1 408-9447	ELB-FL2 408-9454	ELB-FL3 408-9488	ELB-FL4 408-9496	ELB-FL5 408-9504

Strong 30-40mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	REGULAR	ES-FR1 400-5468	ES-FR2 400-5476	ES-FR3 400-5484	ES-FR4 400-5492	ES-FR5 400-5500
	LONG	ES-FL1 400-5419	ES-FL2 400-5427	ES-FL3 400-5435	ES-FL4 400-5443	ES-FL5 400-5450
 BLACK	REGULAR	ESB-FR1 409-0064	ESB-FR2 409-0049	ESB-FR3 409-0072	ESB-FR4 408-9090	ESB-FR5 408-9116
	LONG	ESB-FL1 409-0023	ESB-FL2 409-0056	ESB-FL3 409-0031	ESB-FL4 408-9074	ESB-FL5 408-9082

easywrap

Thigh with Knee

Material Composition: polyamide, elastane

The easywrap compression wrap system for lower limb comes in three parts, which work in isolation or combination, and are designed to be worn over a liner.


The easywrap Thigh with Knee is a two part solution for the upper leg.



Compression

Class	mmHg	Well suited for		
		LYMPHATIC	VENOUS	LIPOEDEMA
Light	20-30	Mild oedema extending into the thigh, Fragile skin, palliative care.	Telangiectasia, thread, spider or broken veins.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.
Strong	30-40	Moderate oedema extending into the thigh, moderate shape distortion, deep pitting tissues. Post liposuction for lymphoedema.	Recurrent/ Varicose veins.	Large deformations of fat especially on the thighs and knees. Lipo-lymphoedema.

Sizing



Circumferences cm	XS	S	M	L	XL
g Just below the groin	60-70	65-75	70-80	75-90	85-100
e¹ Bottom of the thigh	47-52	53-59	60-66	67-73	74-80
Contoured Lengths cm	SHORT	REGULAR	TALL		
e¹-g Bottom of thigh to groin	16-20	20-24	25-30		

Comfortable

Patented zone behind knee

Padding and elliptical cut prevent pinching or cutting in, enhancing comfort and mobility

Wide bands

For easy donning



Codes

Available on **Drug Tariff**

Light 20-30mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	SHORT	EL-TS1 400-5062	EL-TS2 400-5070	EL-TS3 400-5088	EL-TS4 400-5096	EL-TS5 400-5104
	REGULAR	EL-TR1 400-5013	EL-TR2 400-5021	EL-TR3 400-5039	EL-TR4 400-5047	EL-TR5 400-5054
	TALL	EL-TT1 400-5112	EL-TT2 400-5120	EL-TT3 400-5138	EL-TT4 400-5146	EL-TT5 400-9816
 BLACK	SHORT	ELB-TS1 408-9223	ELB-TS2 408-9215	ELB-TS3 408-9207	ELB-TS4 408-9181	ELB-TS5 408-9199
	REGULAR	ELB-TR1 408-9231	ELB-TR2 408-9256	ELB-TR3 408-9272	ELB-TR4 408-9280	ELB-TR5 408-9298
	TALL	ELB-TT1 408-9124	ELB-TT2 408-9132	ELB-TT3 408-9140	ELB-TT4 408-9165	ELB-TT5 408-9173

Strong 30-40mmHg

COLOUR	LENGTH	XS	S	M	L	XL
 SAND	SHORT	ES-TS1 400-5666	ES-TS2 400-5674	ES-TS3 400-5682	ES-TS4 400-5690	ES-TS5 400-5708
	REGULAR	ES-TR1 400-5617	ES-TR2 400-5625	ES-TR3 400-5633	ES-TR4 400-5641	ES-TR5 400-5658
	TALL	ES-TT1 400-5716	ES-TT2 400-5724	ES-TT3 400-5732	ES-TT4 400-5740	ES-TT5 400-5757
 BLACK	SHORT	ESB-TS1 408-8886	ESB-TS2 408-8910	ESB-TS3 408-8902	ESB-TS4 408-8845	ESB-TS5 408-8860
	REGULAR	ESB-TR1 408-8779	ESB-TR2 408-8795	ESB-TR3 408-8811	ESB-TR4 408-8787	ESB-TR5 408-8837
	TALL	ESB-TT1 408-8712	ESB-TT2 408-8746	ESB-TT3 408-8753	ESB-TT4 408-8738	ESB-TT5 408-8704

easywrap

Knee

Material Composition: polyamide, elastane

The easywrap compression wrap system for lower limb comes in three parts, which work in isolation or combination, and are designed to be worn over a liner.


The easywrap Knee is a one-piece solution for the Knee.



Compression

Class	mmHg	Well suited for		
		LYMPHATIC	VENOUS	LIPOEDEMA
Light	20-30	Mild oedema extending into the thigh, Fragile skin, palliative care.	Telangiectasia, thread, spider or broken veins.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.
Strong	30-40	Moderate oedema extending into the thigh, moderate shape distortion, deep pitting tissues. Post liposuction for lymphoedema.	Recurrent/ Varicose veins.	Large deformations of fat especially on the thighs and knees. Lipo-lymphoedema.

Sizing



Circumferences cm	XS	S	M	L	XL
e¹ Bottom of the thigh	47-52	53-59	60-66	67-73	74-80
C Widest calf	38-45	44-52	50-60	55-65	60-70

Easy to wear

Wide bands

For easy donning

Comfortable

Patented zone behind knee

Padding and elliptical cut prevent pinching or cutting in, enhancing comfort and mobility



Codes

Light 20-30mmHg

COLOUR	XS	S	M	L	XL
 SAND	EL-KP1	EL-KP2	EL-KP3	EL-KP4	EL-KP5
 BLACK	ELB-KP1	ELB-KP2	ELB-KP3	ELB-KP4	ELB-KP5

Strong 30-40mmHg

COLOUR	XS	S	M	L	XL
 SAND	ES-KP1	ES-KP2	ES-KP3	ES-KP4	ES-KP5
 BLACK	ESB-KP1	ESB-KP2	ESB-KP3	ESB-KP4	ESB-KP5

easywrap

Fusion Liners

Material Composition: polyamide, elastane

easywrap fusion is a liner and foot piece in one, which makes it the perfect partner for an easywrap leg piece.

The liner provides compression to the foot, and half-compression at the ankle for overlap with an easywrap leg piece, offering a low-profile, smart alternative to regular foot wraps.

“Discovered Fusion liners when I found it difficult to use the foot wraps. Very comfortable to wear, the liner is soft and comfortable against my sensitive skin.”


Patient review, lymphshop.com



Compression

Class	mmHg	Well suited for		
		LYMPHATIC	VENOUS	LIPOEDEMA
Light	20-30	Mild to moderate oedema with slight to moderate pitting tissues, mild shape distortion.	Mixed aetiology ulcers. Telangiectasia, thread,spider or broken veins. Recurrent/ Varicose veins. Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.
Strong	30-40	Moderate to severe oedema with moderate shape distortion and deep pitting tissues. Where stemmer sign is negative.	Recurrent/Active venous Leg Ulcer. Healed venous ulcer.	Large deformations of fat or lipo-lymphoedema.

Sizing



Circumferences cm		XS	S	M	L	XL
c	Widest point of calf	≤ 45	≤ 50	≤ 60	≤ 65	≤ 70
b	Ankle	19-24	25-29	30-34	35-40	41-46
a¹	Mid-foot	19-24	25-29	30-34	35-40	41-46

Easy to wear

Compressive foot zone

Just like wearing socks

Half-compression zone at the ankle

To maintain graduated compression and 50% overlap when worn with an easywrap Leg.



Codes

Available on **Drug Tariff**

Open Toe

COLOUR	COMPRESSION	XS	S	M	L	XL
 SAND	Light 20-30mmHg	EL-FD1 415-3110	EL-FD2 415-3128	EL-FD3 415-3136	EL-FD4 415-3144	EL-FD5 415-3151
	Strong 30-40mmHg	ES-FD1 415-3243	ES-FD2 415-3235	ES-FD3 415-3276	ES-FD4 415-3268	ES-FD5 415-3250
 BLACK	Light 20-30mmHg	ELB-FD1 415-3219	ELB-FD2 415-3185	ELB-FD3 415-3193	ELB-FD4 415-3201	ELB-FD5 415-3227
	Strong 30-40mmHg	ESB-FD1 415-3284	ESB-FD2 415-3292	ESB-FD3 415-3318	ESB-FD4 415-3326	ESB-FD5 415-3300

Closed Toe

COLOUR	COMPRESSION	XS	S	M	L	XL
 SAND	Light 20-30mmHg	EL-FD1-CT	EL-FD2-CT	EL-FD3-CT	EL-FD4-CT	EL-FD5-CT
	Strong 30-40mmHg	ES-FD1-CT	ES-FD2-CT	ES-FD3-CT	ES-FD4-CT	ES-FD5-CT
 BLACK	Light 20-30mmHg	ELB-FD1-CT	ELB-FD2-CT	ELB-FD3-CT	ELB-FD4-CT	ELB-FD5-CT
	Strong 30-40mmHg	ESB-FD1-CT	ESB-FD2-CT	ESB-FD3-CT	ESB-FD4-CT	ESB-FD5-CT

easywrap

Standard Liners

Material Composition: polyamide, elastane

Forms a barrier between the limb and easywrap which absorbs sweat and traps odours, reducing the need to wash the easywrap.

easywrap Standard liners are made of a non-compressive, stretchy, soft knit for maximum comfort and skin protection.

“Soft and comfortable against my sensitive skin. They wash well and dry quickly, good quality liners.”

Patient review, lymphshop.com



Codes

Available on Drug Tariff

COLOUR	LEG SIZE		
 SAND	XS/S/M	E-L-AD 425-2649	E-L-AG 425-2631
	L/XL	E-L-ADW 425-2623	E-L-AGW 425-3951
 BLACK	XS/S/M	EB-L-AD 425-2664	EB-L-AG 425-2656
	L/XL	EB-L-ADW 425-3944	EB-L-AGW 425-2672

*One pair of Below Knee Liners will be included free of charge when purchasing an easywrap Foot and Leg together.

easywrap

Extenders

Extends the circumference of a single band of easywrap by 10cm. When oedema in this area is reduced, the extender can then be easily removed.

Useful where patients measure outside the standard circumference size of a single measurement point, reducing the need for ordering an entire custom made garment.



Codes

COLOUR	LENGTH	WIDTH		
		8cm	12cm	16cm
 SAND	10cm	EE-0810	EE-1210	EE-1610
 BLACK	10cm	EEB-0810	EEB-1210	EEB-1610

easywrap

Hand

Material Composition: polyamide, elastane

The easywrap compression wrap system for upper limb comes in two parts, which work in isolation or in combination.

easywrap Hand features two self-adjustable Velcro® tabs, which fasten around a silky smooth mitten. Can also be successfully paired with a *Microfine Glove* to add compression to fingers.



Compression

Class	mmHg	Well suited for
		LYMPHATIC
Light	20-30	Mild to moderate shape distortion and slight to moderate tissue pitting. Post liposuction.

Sizing



Circumferences cm	XS	S	M	L
b Palm at fold of thumb	<19	19-21	21-23	>24
a Palm at base of fingers	<18	18-20	20-22	>22

Easy to wear

Velcro® fasteners

Lengthened wrist section to ensure overlap with Easywrap Arm

4-way stretch soft fabric

for maximum comfort whilst maintaining range of motion.



Codes

Available on Drug Tariff

Left 20-30mmHg

COLOUR	XS	S	M	L
 SAND	EL-H1L 400-5377	EL-H2L 400-5385	EL-H3L 400-5393	EL-H4L 400-5401
 BLACK	ELB-H1L 408-9603	ELB-H2L 408-9637	ELB-H3L 408-9629	ELB-H4L 408-9611

Right 20-30mmHg

COLOUR	XS	S	M	L
 SAND	EL-H1R 400-5369	EL-H2R 400-5336	EL-H3R 400-5344	EL-H4R 400-5351
 BLACK	ELB-H1R 408-9595	ELB-H2R 408-9587	ELB-H3R 408-9579	ELB-H4R 408-9561

Arm

Material Composition: polyamide, elastane

The easywrap compression wrap system for upper limb comes in two parts, which work in isolation or combination.


easywrap Arm features self-adjustable vecro tabs, which fasten around a silky smooth sleeve.



Compression

Class	mmHg	Well suited for
		LYMPHATIC
Light	20-30	Mild to moderate shape distortion and slight to moderate tissue pitting. Post liposuction.
		LIPOEDEMA
		Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.

Sizing



Circumferences cm	S	M	L	XL
g Axilla	22-31	29-39	32-45	36-50
e Elbow Crease	20-27	25-34	30-40	32-43
c Wrist Crease	14-18	16-21	19-25	19-25
Contoured Lengths cm	SHORT	REGULAR	LONG	
c-g Axilla to wrist	40-44	44-48	48-52	

Easy to wear

Velcro® fasteners

Can be applied and adjusted with one hand

4-way stretch soft fabric

for maximum comfort whilst maintaining range of motion.



Codes

Available on **Drug Tariff**

Left 20-30mmHg

COLOUR	LENGTH	S	M	L	XL
 SAND	SHORT	EL-AS2L 400-5187	EL-AS3L 400-5195	EL-AS4L 400-5203	EL-AS5L 405-8236
	REGULAR	EL-AR2L 400-5245	EL-AR3L 400-5252	EL-AR4L 400-5260	EL-AR5L 405-8194
	LONG	EL-AL2L 400-5302	EL-AL3L 400-5310	EL-AL4L 400-5328	EL-AL5L 405-8186
 BLACK	SHORT	ELB-AS2L 408-8373	ELB-AS3L 408-8381	ELB-AS4L 408-8399	ELB-AS5L 408-8407
	REGULAR	ELB-AR2L 408-8290	ELB-AR3L 408-8308	ELB-AR4L 408-8316	ELB-AR5L 408-8324
	LONG	ELB-AL2L 408-8241	ELB-AL3L 408-8225	ELB-AL4L 408-8233	ELB-AL5L 408-8217

Right 20-30mmHg

COLOUR	LENGTH	S	M	L	XL
 SAND	SHORT	EL-AS2R 400-5153	EL-AS3R 400-5161	EL-AS4R 400-5179	EL-AS5R 405-8202
	REGULAR	EL-AR2R 400-5211	EL-AR3R 400-5229	EL-AR4R 400-5237	EL-AR5R 405-8210
	LONG	EL-AL2R 400-5278	EL-AL3R 400-5286	EL-AL4R 400-5294	EL-AL5R 405-8228
 BLACK	SHORT	ELB-AS2R 408-8415	ELB-AS3R 408-8423	ELB-AS4R 408-8431	ELB-AS5R 408-8449
	REGULAR	ELB-AR2R 408-8332	ELB-AR3R 408-8340	ELB-AR4R 408-8357	ELB-AR5R 408-8365
	LONG	ELB-AL2R 408-8258	ELB-AL3R 408-8266	ELB-AL4R 408-8274	ELB-AL5R 408-8282

easywrap

Additional Bands

Material Composition: polyamide, elastane

A single band of easywrap lower limb material with Velcro® fasteners attached at either end.

They can be used to provide additional support and compression to either the upper or lower limb.



Compression

Class	mmHg	Well suited for
Light	20-30	Mild to moderate shape distortion, moderate pitting.
Strong	30-40	Moderate shape distortion, massive localised lymphoedema, obesity related lymphoedema, moderate skin changes.

Codes

Light 20-30mmHg			Strong 30-40mmHg		
	8 x 40cm	8 x 60cm		8 x 40cm	8 x 60cm
 SAND	EL-0840	EL-0860	 SAND	ES-0840	ES-0860
 BLACK	ELB-0840	ELB-0860	 BLACK	ESB-0840	ESB-0860
Unit of Sale	Single	Single	Unit of Sale	Single	Single

	16 x 45cm	16 x 150cm		16 x 45cm	16 x 150cm
 SAND	EL-1645	EL-16150	 SAND	ES-1645	ES-16150
 BLACK	ELB-1645	ELB-16150	 BLACK	ESB-1645	ESB-16150
Unit of Sale	Single	Single	Unit of Sale	Single	Single

easywrap

Custom Made

Material Composition: polyamide, elastane

For patients who fall outside of the sizing chart ranges.

Easy-to-use, step-by-step measurement form that guides you through the process of ordering a custom made easywrap.

We aim to provide additional flexibility wherever possible. Please call or email us to discuss.

Codes

		Foot	Leg	Thigh	Knee	Arm	Hand
Light	20-30mmHg	EL-FC	EL-LC	EL-TC	EL-KC	EL-AC	EL-HC
Strong	30-40mmHg	ES-FC	ES-LC	ES-TC	ES-KC	-	-
Unit of Sale		Single	Single	Single	Single	Single	Single





Comfiwave



Comfiwave

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves grooves on a patient's skin that help soften the skin while promoting lymphatic drainage.



What makes it so patient-friendly?

Turn the page to discover all features [Page 94](#) ▶



Below Knee

[Page 98](#) ▶



Thigh

[Page 100](#) ▶



Toecap

[Page 102](#) ▶



Shorts

[Page 104](#) ▶



Abdominal Band

[Page 106](#) ▶



Sleeve

[Page 108](#) ▶



Sleeve + Mitten

[Page 110](#) ▶



Glove

[Page 112](#) ▶



Glove to Elbow

[Page 114](#) ▶



Glove to Axilla

[Page 116](#) ▶



Breast Band

[Page 118](#) ▶



Neck

[Page 120](#) ▶



Thigh + Hip

[Page 122](#) ▶



Kids

[Page 124](#) ▶

Comfiwave is a soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient's skin that help soften the skin while promoting lymphatic drainage.

Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and even overnight.

Compression to drift-off in



Innovative

Comfiwave garments use proprietary 'low-resistance' yarns which allow them to be effortlessly donned and doffed.

Conforming

3D knitting technology creates a seamless knit which delivers a gentle yet targeted, therapeutic pressure of 15-21_{mmHg}.



Breathable

High air permeability for coolness all year round, day or night.

Organic cotton

Soft, organic cotton against the skin ensures maximum comfort.





For night time

Continue compression therapy overnight in supreme comfort with Comfiwave's fabric providing gentle compression while you sleep.

Making waves

Comfiwave's wave patterning mimics the pressure differential massage effect that is ordinarily found in cut-foam, chip-bag style garments.

Effective

“Can be used for stand-alone treatment or as an adjunct to other treatment modalities, and for 24 hours a day if needed.”

- British Journal of Community Nursing

“Patients reported that Comfiwave was comfortable and easy to use.”

D	F	G	Limb volume A	Unaffected	% difference
25	29	31.5	2329	1954	+16%
24	27.5	31	2139	1953	+8%

through her teenage years with bilateral lower leg oedema and recurrent episodes of cellulitis.

On initial assessment, there was evidence of bilateral lymphoedema below the knee, with both feet showing positive Stemmer's sign, mild hyperkeratosis and post-cellulitis skin discolouration on the left gaiter area. The left leg was 30% bigger than the right, with an ankle circumference of 51cm and calf of 65cm. There were distinct skin folds distorting the shape around her left ankle, as well as foot oedema, which was restricting footwear (Figure 4). The limb size also restricted clothing, and the patient was only able to wear very loose trousers or long skirts.

Initial advice of skin care, exercise, elevation and weight management were discussed, and 2 weeks of daily treatment were booked. Miss H was signed off work by her GP to enable her to attend the sessions. Daily treatment commenced in October 2018 and involved alternating days of MLD or intermittent compression using Haddenham

Figure 5. Miss H wearing Comfiwave on the left leg

Figure 6. Miss H's legs after using Comfiwave for 2 weeks

PRODU



ments
to the
y were
improve
ment of
that this
ment or as
24 hours
efficacy in
decreasing
edema and
ient, family

ent: a prospective
3-541. <https://doi.org/10.1016/j.cnn.2017.09.001>

staker J. Examination
in tissue oxygenation
429-435. <https://doi.org/10.1016/j.cnn.2017.09.001>



Simple sizing

Easy measuring with colour coded sizing system that eliminates the need for uncomfortable sewn-in labels. 3 red dots = Medium Regular.

“ I love this garment, I have been wearing it overnight for the last week and it's making a massive difference!

- Patient review, lymphshop.com



Kid-friendly

A one-of-a-kind addition to treatment options for children with lymphoedema.

Guaranteed

6 months of use

Material composition

Cotton, Polyester, Polyurethane

Availability

All variants held in stock

Comfiwave

Below Knee

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for		
	LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA
15-21	Mild to moderate oedema with or without pitting and/or fibrosis. Palliative care.	Telangiectasia, thread,spider or broken veins. Recurrent/ Varicose veins. Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.

Sizing



Circumferences cm		S	M	L	XL
c	Widest Calf	25 - 32	30 - 38	36 - 46	44 - 56
b	Ankle	17 - 22	20 - 26	22 - 30	27 - 35
Lengths cm		SHORT	REGULAR	LONG	
A-d	Floor to Below Knee	32 - 37	37 - 42	42 - 48	

Comfort

Stretchy lateral toe zone

For enhanced fit at ball of foot

Shaped ankle zone

Increases comfort and prevents creasing

Support

Enhanced instep compression

Added support around the top of the foot

Contoured heel section

Helps keep the garment in place

“This has been a real game changer for me. My leg has never been as good.”

Patient review, lymphshop.com



Codes

Available on **Drug Tariff**

	S		M		L		XL	
SHORT	CW1-AD2S 422-1412	■	CW1-AD3S 422-1362	■	CW1-AD4S 422-1370	■	CW1-AD5S 422-1388	■
REGULAR	CW1-AD2R 422-1446	■	CW1-AD3R 422-1420	■	CW1-AD4R 422-1396	■	CW1-AD5R 422-1404	■
LONG	CW1-AD2L 422-1438	■	CW1-AD3L 422-1479	■	CW1-AD4L 422-1453	■	CW1-AD5L 422-1461	■

Unit of Sale

Single

Comfiwave

Thigh High

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient's skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for			
	LYMPHOEDEMA	VENOUS DISEASE	LIPOEDEMA	
15-21	Mild to moderate oedema with or without pitting and/or fibrosis. Palliative care.	Telangiectasia, thread,spider or broken veins. Recurrent/ Varicose veins. Pigmentation or eczema, Lipodermatosclerosis or atrophie blanche. Corona phlebectatica.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.	

Sizing



Circumferences cm	S	M	L	XL	XXL
g Thigh Top	40 - 48	46 - 56	54 - 65	62 - 75	72 - 92
c Widest Calf	25 - 32	30 - 38	36 - 46	44 - 56	52 - 64
b Ankle	17 - 22	20 - 26	22 - 30	27 - 35	27 - 35
Lengths cm	SHORT		REGULAR		LONG
A-g Floor to Thigh Top	64 - 70		70 - 77		77 - 84

Comfort

Stretchy lateral toe zone

For enhanced fit at ball of foot

Shaped ankle zone

Increases comfort and prevents creasing

Support

Enhanced instep compression

Added support around the top of the foot

Contoured heel section

Helps keep the garment in place

“ Really comfortable to wear at night. Unlike other products, doesn’t make my leg hot when temperature is high.”

Patient review, lymphshop.com

Codes

Available on Drug Tariff

	S		M		L		XL		XXL	
SHORT	CW1-AG2S 422-1487	■	CW1-AG3S 422-1560	■	CW1-AG4S 422-1552	■	CW1-AG5S 422-1503	■	CW1-AG6S 429-5507	■
REGULAR	CW1-AG2R 422-1529	■	CW1-AG3R 422-1537	■	CW1-AG4R 422-1511	■	CW1-AG5R 422-1578	■	CW1-AG6R 429-5515	■
LONG	CW1-AG2L 422-1586	■	CW1-AG3L 422-1495	■	CW1-AG4L 422-1545	■	CW1-AG5L 422-1594	■	CW1-AG6L 429-5499	■

Unit of Sale

Single

Comfiwave

Toecap

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.

Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Toe swelling, base of toe swelling, positive Stemmer sign

Sizing



Circumferences cm	M	L
a Ball of Foot	20 - 27	26 - 34

Comfort

3D knitting of the forefoot

Creates anatomical shaping for a better fit

Wider circumference at small toe

Creates an expanded ‘pinky toe’ opening to increase comfort

Support

Sturdy knit

Delivers a firm and supportive fit



“ This is an amazing garment. My toes go from fat bloated sausages to toe-like toes after I have been wearing it, and it fits really well.”

Patient review, lymphshop.com

Codes

	M		L	
REGULAR	CW1-ZA3		CW1-ZA4	
Unit of Sale	Single			

Comfiwave

Shorts

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Female genital oedema, male genital oedema, suprapubic swelling, hip/buttock swelling, obesity related lymphoedema

Sizing



Circumferences cm		S	M	L	XL	XXL
t	WAIST	50 - 68	66 - 86	84 - 104	102 - 124	120 - 146
g	THIGH TOP	40 - 48	46 - 56	54 - 65	62 - 75	72 - 92

Comfort

3D Knitting

Conforming, seamless knit to target sensitive areas

Breathable cotton

Promotes airflow and reduces moisture and infection risk

Versatility

Pad and pull-up friendly

Compatible with incontinence products

Unisex design

Conforms for men or women



“ They felt amazing as soon as I put them on...so comfortable, soft and have the perfect amount of squeeze!”

Codes

	S		M		L		XL		XXL
REGULAR	CW1-GT2		CW1-GT3		CW1-GT4		CW1-GT5		CW1-GT6

Unit of Sale

Single

Comfiwave

Abdominal Band

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for	
	LYMPHOEDEMA	LIPOEDEMA
15-21	Male or female suprapubic oedema. Obesity related lymphoedema.	Lipoedema of the pelvis, hips and buttocks.

Sizing



Circumferences cm		M	L	XL
t	Waist	60 - 75	70 - 85	90 - 120
j	Hips	80 - 100	90 - 105	110 - 140

Comfort

3D mapped compression

Conforms to body shape where traditional compression garments can't.

Versatility

Use in combination

Can be used with the Comfiwave Thigh or Thigh+Hip garment.

Use post-operatively

Can be used over dressings for covering wounds, including fungating lesions.

“The first garment that offers help in this difficult to treat area. Very comfortable.”

Male lymphoedema patient, with secondary lymphoedema in the pubic area and upper thighs on both legs.

Codes

	M	L	XL
REGULAR	CW1-HT3 436-1556	CW1-HT4 436-1564	CW1-HT5 436-1549
Unit of Sale		Single	




Comfiwave

Sleeve

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for	
	LYMPHOEDEMA	LIPOEDEMA
15-21	Oedema which does not manifest in the hand, pitting and fibrosis to the forearm, palliative care.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.

Sizing



Circumferences cm		S	M	L
g	Axilla	25 - 30	30 - 40	40 - 50
d	Mid-point of forearm	23 - 26	26 - 30	30 - 38
c	Wrist	13 - 18	18 - 22	22 - 28
Lengths cm		SHORT	REGULAR	LONG
c-g	Wrist to Axilla	36 - 42	40 - 46	44 - 50

Comfort

Elbow stretch zone

Pressure is evenly distributed to prevent tightness when the elbow bends.

Conforms to forearm and elbow

Prevents oedema from pooling in these difficult areas

Versatility

Can be used with Comfiwave Glove

Don and doff the glove as required without having to remove compression from your arm



Codes

Also available on **Drug Tariff**

	S		M		L
SHORT	CW1-CG2S 422-1198	■	CW1-CG3S 422-1180	■	CW1-CG4S 422-1164
REGULAR	CW1-CG2R 422-1172	■	CW1-CG3R 422-1149	■	CW1-CG4R 422-1156
LONG	CW1-CG2L	■	CW1-CG3L	■	CW1-CG4L

Unit of Sale

Single

Comfiwave

Sleeve with Mitten

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Oedema extending from the dorsum of the hand throughout the arm, pitting and fibrosis with mild shape distortion, palliative care

Sizing



Circumferences cm		S	M	L
g	Axilla	25 - 30	30 - 40	40 - 50
d	Mid-point of forearm	23 - 26	26 - 30	30 - 38
c	Wrist	13 - 18	18 - 22	22 - 28
a	Pal at base of fingers	16 - 22	18 - 24	20 - 26
Lengths cm		SHORT	REGULAR	
c-g	Wrist to axilla	36 - 42	40 - 46	

Comfort

Elbow stretch zone

Pressure is evenly distributed to prevent tightness when the elbow bends.

“Excellent, very comfortable to sleep in. Having had to do massage and bandaging in the past, this is simple and works!”

Patient review, lymphshop.com



Codes

Also available on **Drug Tariff**

	S		M		L	
SHORT	CW1-AH2S 422-1214	■	CW1-AH3S 422-1206	■	CW1-AH4S 422-1222	■
REGULAR	CW1-AH2R 422-1230	■	CW1-AH3R 422-1248	■	CW1-AH4R 422-1255	■

Unit of Sale

Single

Comfiwave

Glove

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Positive Stemmer sign, mild shape distortion, pitting/fibrosis to the dorsum of the hand, palliative care

Sizing



Circumferences cm	S	M	L
b Base of thumb	17 - 22	20 - 25	23 - 28
c Wrist	14 - 20	16 - 24	22 - 28

Comfort

Reduced elastic at fingertips

Increased comfort and circulation

Knitted gussets

Gussets in between the fingers improves comfort and reduces chafing

Versatility

Can be used with Comfiwave Sleeve

Don and doff the glove as required without having to remove compression from your arm



Codes

	S		M		L	
REGULAR	CW1-ZC2		CW1-ZC3		CW1-ZC4	
Unit of Sale	Single					

Comfiwave

Glove to Elbow

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient's skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Oedema extending from the fingers (positive Stemmer sign) to just below the elbow, pitting and fibrosis to the forearm, mild shape distortion.

Sizing



Circumferences cm		S	M	L
b	Base of thumb webbing	17 - 22	20 - 25	23 - 28
c	Wrist	14 - 20	16 - 24	22 - 28
d	Mid-point of forearm	23 - 28	26 - 32	30 - 38
Lengths cm		REGULAR		
c-e	Wrist to elbow	20 - 24		

Comfort

Reduced elastic at fingertips

Increased comfort and circulation

3D finger webbing

Extra cotton in the finger gussets
relieves pressure and prevents chafing



Codes

	S		M		L	
REGULAR	CW1-ZE2		CW1-ZE3		CW1-ZE4	
Unit of Sale	Single					

Comfiwave

Glove to Axilla

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient's skin that help soften the skin while promoting lymphatic drainage.

Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Oedema extending from the fingers throughout the whole arm, pitting and fibrosis, palliative care, mild to moderate shape distortion

Sizing

	Circumferences cm	S	M	L
	g Axilla	25 - 30	30 - 40	40 - 50
	d Mid-point of forearm	23 - 26	26 - 30	30 - 38
	c Wrist	13 - 18	18 - 22	22 - 28
	a Palm at base of fingers	16 - 22	18 - 24	20 - 26
	Lengths cm	SHORT	REGULAR	
	c-g Wrist to Axilla	36 - 42	40 - 46	

Comfort

Reduced elastic at fingertips

Increased comfort and circulation

3D Finger webbing

Extra cotton in the finger gussets
relieves pressure and prevents chafing

“I’m very impressed with how this works I wear it mainly at night and my fingers look amazing.”

Patient review, lymphshop.com



Codes

	S		M		L
SHORT	CW1-ZH2S	■	CW1-ZH3S	■	CW1-ZH4S
REGULAR	CW1-ZH2R	■	CW1-ZH3R	■	CW1-ZH4R

Unit of Sale

Single

Comfiwave

Breast Band

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Breast discomfort, mild to severe lymphoedema with or without pitting tissues and shape distortion. Post breast surgery.

Sizing



Circumferences cm	S	M	L	XL
Above breast	70 - 88	80 - 102	90 - 116	100 - 130
Over breast	70 - 92	80 - 108	90 - 126	100 - 144
Below breast	70 - 80	80 - 92	90 - 104	100 - 116

Comfort

Gentle bust coverage

A softer knit around the chest for a comfortably reduced compression level.

Support

Underarm support

Increased elastic thread count under the arms for stronger compression where you need it most.

Use post-operatively

Following breast augmentation or for those who have had breast conserving surgery such as lumpectomy.



“WOW!! This has significantly increased the effective management of my breast lymphoedema. I only wear it overnight and, every morning, my ‘lymph boob’ feels very, very close to just ‘normal’.”

Lymphoedema patient

Codes

	S		M		L		XL	
REGULAR	CW1-TN2 436-1655		CW1-TN3 436-1648		CW1-TN4 436-1630		CW1-TN5 436-1770	
Unit of Sale	Single							

Comfiwave

Neck

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.


Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for
15-21	Submental and submandibular lymphoedema with or with pitting and shape distortion. Post surgery.

Sizing



Circumferences cm	S	M	L
a Base of neck	32 - 38	37 - 43	42 - 47
Lengths cm	SHORT	REGULAR	
a-b Base of neck to top of neck	5-7	7-9	

Comfort

Adaptive Knitting

Specially designed to accommodate drastic circumference changes from chin to neck.

Curved ear line

Features a subtle curve below the ears for a naturally contoured fit.

“The collar is keeping my lymphoedema completely under control without using the kinetic tape or self-massage. I can feel it working, move the skin, dispersing the fluid all of the time it’s on my neck.”

Lymphoedema patient

Codes

	S		M		L	
SHORT	CW1-HN2S	■	CW1-HN3S	■	CW1-HN4S	■
REGULAR	CW1-HN2R	■	CW1-HN3R	■	CW1-HN4R	■

Unit of Sale

—

Single

—

Comfiwave

Thigh + Hip

Material Composition: Cotton, Polyester, Polyurethane

A soft compressive fabric, with a distinctive wave pattern knit. The wave profile creates a pressure differential, that leaves lines/grooves on a patient’s skin that help soften the skin while promoting lymphatic drainage.

Comfiwave offers a unique additional option for patients to use in their self management regimen, which can be worn in a more relaxed setting at home and can even be worn overnight.



Compression

mmHg	Well suited for	
	LYMPHOEDEMA	LIPOEDEMA
15-21	Mild to severe oedema of the thigh extending into the hip with or without pitting and shape distortion.	Smooth or uneven skin with an increase of enlarged subcutaneous fat tissue or lipomas.

Sizing

	Circumferences cm	S	M	L	XL
	g Thigh top	40 - 48	46 - 56	54 - 65	62 - 75
	e Knee	27 - 35	33 - 39	41 - 50	48 - 60
	Lengths cm	SHORT		REGULAR	
e-g Knee to thigh top	20 - 26	25 - 31			

Comfort

Hip flap

allows for wearing under normal underwear, where an abdominal band may not be needed.

Versatility

Use in combination

Can be combined with Comfiwave abdominal band and/or below knee for whole leg treatment



Codes

	S		M		L		XL	
SHORT	CW1-ET2S	■	CW1-ET3S	■	CW1-ET4S	■	CW1-ET5S	■
REGULAR	CW1-ET2R	■	CW1-ET3R	■	CW1-ET4R	■	CW1-ET5R	■

Unit of Sale

Single

Comfiwave

Kids Below Knee

Material Composition: Cotton, Polyester, Polyurethane

A one-of-a-kind addition to treatment options for children with lymphoedema. Optimised for infants and children with a lighter, softer feel, prioritising comfort during wear.


Comfiwave Kids offers a unique option to assist parents in the management of children with lymphoedema, which can be worn in a more relaxed setting at home and even overnight.



Compression

mmHg	Well suited for
15-21	Primary lymphoedema all stages in toddlers and children

Sizing



Circumferences cm		PAEDIATRIC			
C	WIDEST CALF	20 - 33			
b	ANKLE	14 - 22			
Lengths cm		SHORT	REGULAR	LONG	EXTRA LONG
A-d	FLOOR TO BELOW KNEE	17 - 22	22 - 27	27 - 32	32 - 37

Comfort

Lighter

The cotton ply has been reduced by 20%, resulting in an even lighter and more comfortable fabric for children.

Friendly

New Kids on the box

Comfiwave’s packaging has had a makeover to appeal to children, with simple instructions to help them feel in control. Plus, shiny stickers in every box!



Codes

	P	
SHORT	CW1-AD0S	■
REGULAR	CW1-AD0R	■■
LONG	CW1-AD0L	■■■
EXTRA LONG	CW1-AD0X	■■■■
Unit of Sale	Single	



etogarments
Prairie Wear




etogarmments

Garments for head, neck, breast, genital, and torso oedema, as well as garments suitable for lipoedema and breast care.



Find out more

Turn the page to explore [Page 128](#) ►




eto Bra

Page 130 ►



eto Breast Band

Page 132 ►



eto Head & Neck

Page 133 ►




eto Vest

Page 134 ►



eto Scrotal Support

Page 135 ►



eto Shorts

Page 136 ►



eto Grace

Page 138 ►



eto Abdominal Pouch

Page 140 ►

Prairie Wear

Supportive, secure & comfortable compression bras for post-surgical recovery, lymphedema management or active recovery.



Find out more

Explore the range [Page 142](#) ►



HuggerPRIMA

Page 144 ►



HuggerVIDA

Page 146 ►



HuggerLUMA

Page 148 ►

etogarments

Eto garments encompass a broad range of supportive, breathable compression garments for treatment of oedema extending to the head, neck, breast, and torso, as well as garments suitable for lipoedema.



Comfort

Eto garments are made from a fabric that is silky smooth on the outside and lightly textured on the inside, for low friction and increased comfort.

External seams

Heat-sealed external seams reduce pressure marks on the skin.

“There isn't much out there that can even come near it in terms of what you can do with it - pockets, flies and seams.”

Mark Pearson, St George's Hospital London



Highly configurable

Eto garments are available custom made, meaning extensive sizing options for treating even the most severe cases of swelling and shape distortion.

Guaranteed
6 months of use.

Material composition
Perlon / Lycra

eto garments

Bras

Material Composition: Elastane

A high quality compression bra made from oeko-tex certified fabric, free from harmful chemicals and dyes.

Eto bras are configurable to the wearers specification - from altering the straps and closures, to specifying pockets for prosthesis to paddings for extra comfort.



“My patient absolutely loves this product and finds it incredibly comfortable and easy to wear.”


Feedback from lymphoedema therapist

Compression

mmHg	Well suited for
25-30	Post breast surgery. where there is pitting and some moderate skin changes to the breast.

Sizing

Custom made option available. Download form at www.hadhealth.com/f0004

	Circumferences cm	XS	S	M	L	XL	XXL	XXXL
	Bra size	A	A / B	B / C	C / D	D / E	E / F	F / G
	Above breast	78-82	83-87	88-92	93-97	98-102	103-107	108-112
	Below Breast	65-71	72-77	78-85	86-94	95-100	101-109	110-118
	Waist	60-65	66-72	73-79	80-90	91-98	99-106	107-113
	Hips H603L only	89-93	94-99	101-109	110-119	120-125	126-130	131-135

Configurable

Pockets

Add for a prosthesis or additional padding

Adjustable fastenings

Hook & eye fastenings to the front and straps allow support to be adjusted

Comfortable

Empire foam

Optional additional under-breast padding for extra comfort (pictured right)

Klima-comfort under arm gussets

Add for additional comfort under arm



ETO-H603L

Options

*Available for H601, H603 and H603 only.

Colour	 Beige	 Black	 White
Front fastening*	Hook & Eye Default	Zip -----	
Cups*	Empire foam under breast -----	Pocket	

Strap Width*	Narrow 5cm	Wide 8cm -----
Strap fastening*	Hook & Eye Default	Velcro Free
Openings*	Lower Neckline 4.5cm	Klima-comfort under arm gusset

Codes

				
25-30mmHg	ETO-H601	ETO-H603	ETO-H603L	ETO-H607
Neck Shape	V	Scoop	V	V
Strap width	Wide	Narrow	Narrow	Narrow
Strap fastening	Hook & Eye	Hook & Eye	Hook & Eye	None
Front fastening	Hook & Eye	Hook & Eye	Hook & Eye	None
Unit of Sale	Single	Single	Single	Single

eto garments

Breast Band

Material Composition: Elastane


High quality compression breast band made from oeko-tex certified fabric free from harmful chemicals and dyes.



Compression

Class	mmHg	Well suited for
-	25-30	Post surgery, where oedema is extensive and extra compression is required to reduce the swelling

Sizing

	Circumference cm	XS	S/M	L/XL	XXL
	Above breast	80-90	90-100	100-110	110-120

Options

Colour			
	Beige	Black	White

Codes



25-30mmHg	ETO-602
Unit of Sale	Single

Head & Neck Straps

Material Composition: Elastane

A high quality compression garment for head and neck oedema, made from oeko-tex certified fabric free from harmful chemicals and dyes.



Compression

Class	mmHg	Well suited for
-	25-30	Oedema of the neck, submental, submandibular region extending into the lower cheeks, moderate pitting with some shape distortion

Sizing

Custom made option available. Download form at www.hadhealth.com/f0055

	Circumferences cm	ETO 200, 201, 202, 203, 206		ETO209	
		STANDARD	XL	STANDARD	XL
1 ABOVE EYEBROWS		< 66	66 - 71	< 66	66 - 74
2 POINT OF CHIN TO CROWN OF HEAD		< 66	66 - 71	< 66	66 - 74
3 NECK		< 37	37 - 40	< 37	37 - 41

Options

Colour



Beige



White

Codes



ETO-200



ETO-201



ETO-202



ETO-203



ETO-206



ETO-209

25-30mmHg

Unit of Sale

Single

Single

Single

Single

Single

Single

Custom made by quotation

eto garments

Vest

Material Composition: Elastane

High quality compression vest made from oeko-tex certified fabric free from harmful chemicals and dyes.


Compression vests that provide support and compression in the chest and abdominal/ truncal area, with easy fastening and comfortable zip.



Compression

Class	mmHg	Well suited for
-	25-30	Bilateral mastectomy or unilateral mastectomy in men, pitting and mild to moderate skin changes.

Sizing

	Circumference cm	XS	S	M	L	XL	XXL
 Chest		79-83	84-89	90-99	100-110	111-120	121-125
Waist		75-79	80-85	86-95	96-106	107-115	116-122
Hips		85-89	90-95	96-105	106-116	117-123	124-130

Options

Colour			
	Beige	Black	White

Codes

	
25-30mmHg	
ETO-402	ETO-400
Unit of Sale	Single

eto garments

Scrotal Support

Material Composition: Elastane

High quality compression shorts made from oeko-tex certified fabric free from harmful chemicals and dyes.

Support pouch with an adjustable waistband and option for additional pad at perineum.



Compression


Class	mmHg	Well suited for
-	0	Oedema of the scrotum and penis

Sizing

Custom made option available. Pricing by quotation. Download form at www.hadhealth.com/f0005

The diagram shows a garment with three measurement points: A (waist), B (hip), and C (length). The size chart provides measurements for these points across different sizes.

	S/M	L/XL	XXL					
A	90-105	106-119	120+					
	1	2	3	4	5	6	7	8
B	22-24	24-26	26-28	28-30	30-32	32-34	34-36	36-38
C	23-26	26-29	29-32	32-35	35-38	38-41	41-44	44-47



Please Note: Point 'B' is measured from right to left of groin over widest part. Point 'C' is measured from pubic bone to perineum. We require both measurements in order to manufacture this garment.

Options

Colour			
	Beige	Black	White

Codes

	Scrotal Support
25-30mmHg	ETO-11P
Unit of Sale	Single
Price	£88.00
	Custom made by quotation

eto garments

Female shorts

Material Composition: Elastane

High quality compression shorts made from oeko-tex certified fabric free from harmful chemicals and dyes.

A female genital compression garment with a cotton lined crotch, with a pocket to add additional padding.




Compression

Class	mmHg	Well suited for
-	25-30	Lymphoedema of the suprapubic region and labia

Sizing

Custom Made available by quotation. Download form at www.hadhealth.com/f0106

	Circumferences cm	XS	S	M	L	XL	XXL	XXXL
	Waist	60-65	66-72	73-79	80-90	91-98	99-106	107-114
	Hip	89-93	94-99	100-109	110-119	120-125	126-130	131-135
	Left top	49-51	52-54	55-58	59-61	62-64	65-68	69-72

Options

Colour			
	Beige	Black	White

Codes



Sized

25-30mmHg	ETO-H106
Unit of Sale	Single

Custom made by quotation

eto garments

Male shorts

Material Composition: Elastane

High quality compression shorts made from oeko-tex certified fabric free from harmful chemicals and dyes.

A male genital compression garment with a cotton lined crotch with optional fly-opening.




Compression

mmHg	Well suited for
25-30	Scrotal swelling

Sizing

Custom Made available. Download form at www.hadhealth.com/f0099

	Circumferences cm	XS	S	M	L	XL	XXL	XXXL
	Waist	60-65	66-72	73-79	80-90	91-98	99-106	107-114
	Hip	89-93	94-99	100-109	110-119	120-125	126-130	131-135
	Left top	49-51	52-54	55-58	59-61	62-64	65-68	69-72

Options

Colour			
	Beige	Black	White

Codes



Sized Sized

25-30mmHg	ETO-H99	ETO-H96
-----------	---------	---------

Unit of Sale	Single	Single
--------------	--------	--------

Custom made by quotation

eto garments

Grace Leggings

Material Composition: Elastane

A high quality compression leggings made from oeko-tex certified fabric free from harmful chemicals and dyes.

Soft and stretchy to help donning and doffing. With a silky outer layer and lightly textured fabric against the skin to improve all-day comfort.

“I like the nice feeling and the security it gives me. I wear the Eto Grace more when I am sore and achy, as it is so comfy.”

Feedback from patient with lipoedema



ETO-G133

Compression

mmHg Well suited for

25-30 All types and stages of lipoedema, post liposuction

Sizing

Download measurement form at
www.hadhealth.com/f0060



Comfortable

Breathable

Warp knit with a light open weave to wick heat and moisture away

Advanced microfibres

Soft and stretchy to help donning and doffing.

Convenient

Lightweight

The fabric moves easily with you through your daily work/exercise

Machine washable

at 30 degrees



Codes

			
	Tights	Footless Tights	Capri
25-30mmHg	ETO-G134	ETO-G133A	ETO-G133
Unit of Sale	Single	Single	Single

Options

Colour	 Beige	 Black	 White
Front fastening	None Default	 Zip	 Hook & Eye
Crotch	Closed Default	Open Free	

eto garments

Abdominal Pouch

Material Composition: Elastane

High quality abdominal apron pouch made from oeko-tex certified fabric free from harmful chemicals and dyes.

Soft cotton, custom made abdominal support pouch with adjustable Velcro fastenings, with extra padding around the edge of the pouch.



Compression

Class	mmHg	Well suited for
-	25-30	Obesity related lymphoedema

Sizing

Custom made only. Download measurement form at www.hadhealth.com/f0008

Options

Colour
 Beige

Codes

Custom Made	
25-30mmHg	ETO-12PCM
Unit of Sale	Single
Price	by quotation

etogarments



Full range available in the etogarments brochure

Download www.hadhealth.com/etocat

Prices by quotation

Turn page for ►

Prairie Wear



Haddenham Healthcare proudly offers Prairie Wear support bras, featuring 360-degree, wire-free mapped compression for comfort and support.

Ideal for post-surgery recovery, managing lymphedema, and everyday active wear, these bras combine functionality with comfort.




Simple sizing

Just one circumference measurement is needed to find a perfect fit.

Full coverage

360 degree targeted and mapped compression that hugs the body.



Accommodating

Extenders allow for phased support post-surgery and to adapt to fluctuation in oedema.

Also includes pockets for prosthetics, compression padding, ice and more.

Attractive

Prairie Bras are available in a pleasing mixture of bright and natural tones.



Prairie Wear

HuggerPrima

Material Composition: Nylon, Spandex

Comfortable support bras with 360-degree, wire-free mapped compression. For after surgery, for Lymphoedema, or as everyday active wear.

The HuggerPrima is a medium-to-high compression bra with high coverage.



Comfort

Extenders

For phased support as swelling reduces


Adjustable

Padded straps, and pockets for foams or prosthesis

Compression

Class	mmHg	Well suited for
Med-High	-	Post breast surgery, where there is pitting and some moderate skin changes to the breast

Sizing




Circumferences cm	XS	S	M	L	XL	2XL	3XL
Above breast	66-71	76-81	86-91	96-102	107-112	117-122	127-132

2XL + L Extender



For **PRIMA**: Upsize if you're in-between sizes or you're a DD cup or higher.



Options

Colour



Codes



HuggerPrima



Extender

Bra sizes S-M Bra Sizes L-2XL

HUG

EXT-M

EXT-L

Unit of Sale

Single

Single

Single

Prairie Wear

HuggerVida

Material Composition: Nylon, Spandex

Comfortable support bras with 360-degree, wire-free mapped compression. For after surgery, for Lymphoedema, or as everyday active wear.

The HuggerVida is a low-to-medium compression bra, with low-to-medium coverage.



Comfort

Extenders

For phased support as swelling reduces


Adjustable

Padded straps, and pockets for foams or prosthesis

Compression

Class	mmHg	Well suited for
Low-med	-	As a step down from higher compression as maintenance post oedema reduction

Sizing




Circumferences cm	XS	S	M	L	XL	2XL	3XL
Above breast	66-71	76-81	86-91	96-102	107-112	117-122	127-132

2XL + L Extender



In-between sizes?
For **VIDA**, downsize.



Options

Colour



Codes



HuggerVida



Extender

Bra sizes S-M

Bra Sizes L-2XL

VIDA

VXT-M

VXT-L

Unit of Sale

Single

Single

Single

Prairie Wear

HuggerLuma

Material Composition: Nylon, Spandex

Comfortable support bras with 360-degree, wire-free mapped compression. For after surgery, for Lymphoedema, or as everyday active wear.

The HuggerLuma is a medium-to-high compression sports bra with low coverage.



Comfort

Extenders

For phased support as swelling reduces


Adjustable

Padded straps and back band

Compression

Class	mmHg	Well suited for
Med-High	-	Post-surgery daily wear where oedema fluctuates and higher compression is needed.

Sizing

	Circumferences cm	XS	S	M	L	XL	2XL
	Above breast	66-71	76-81	86-91	96-102	107-112	117-122



For **LUMA**: *Downsize* if you're in-between



Options

COLOUR



Big Sky Blue



Fawn



No Moon Black

Codes



HuggerLuma

LUMA

Unit of Sale

Single



LymphFlow
LymphScanner



LymphFlow LymphScanner

The latest technology in water content detection and pneumatic compression, to supplement a patients existing compression treatment.



**LymphFlow
Go**

Page 152 ►



**LymphFlow
Advance Pump**

Page 154 ►



**LymphFlow
Advance Garments**

Page 158 ►



Lymphscanner

Page 158 ►

LymphFlow Go

Leg Garment

Material Composition: Nylon, Polyester

The LymphFlow Go is a twelve chamber leg garment with an integrated pump module, which means cordless usage and IPC treatment on the go.

Haddenham LymphFlow products utilise pneumatic compression pumps paired with inflatable, multi-chambered garments to provide intermittent pressure, externally to an oedematous limb.



Compression

Class	mmHg	Well suited for
Low	40	Mild or watery oedema with slight to moderate pitting. Varicose veins or healed venous ulcer.
Medium	60	Moderate oedema with moderate shape distortion and woody tissues. Moderate to very deep pitting. Venous skin changes.
High	80	Severe oedema with extensive shape distortion, woody, fibrotic tissues.

Sizing

Circumferences cm

	Regular	Large
g Thigh Top	≤78	≤88

Straight Distance cm

A-g Leg Length	≤81	≤81
-----------------------	-----	-----

Leg

Codes

	REGULAR	LARGE
LFA Garment	LFG3	LFG4

Coming 2026



Go anywhere

Compact

The device packs neatly into its own carry bag, which fits within a car boot or a suitcase.

Cordless

Cordless operation means a simple setup for treatment anywhere, at anytime.

Two treatments

Lymphatic mode

works in a proximal to distal motion, intensively working the whole of the limb before then completing with a distal to proximal flush that lasts approximately 20-35 minutes on a medium intensity setting.

Vascular mode

offers distal to proximal treatment without a reverse flush, lasting for approximately 30-45 minutes on a medium intensity setting.



LymphFlow Advance

Pump Module

The LymphFlow Advance is a twelve chamber pump that connects to a variety of separate garments, which can be used for the treatment of upper and lower limb.

Haddenham LymphFlow products utilise a pneumatic compression pumps paired with inflatable, multi-chambered garments to provide intermittent pressure, externally to an oedematous limb.



Comprehensive

Seven treatment modes

For individualised care

Easy connections


For quick treatment setup

Compression

Class	mmHg	Well suited for
-	30-100	Lymphoedema, post-immobilisation oedema, venous insufficiency

Codes

Pump Module	Service
LFA5438	LFA SERVICE



How it works

The LymphFlow Advance has seven cycles to enable the clinician to focus treatment on a particular area. The LymphFlow Advance will work intensively on the whole limb, on the proximal or distal areas, or on the hand or foot. This allows the therapist to concentrate on the area that requires more intensive treatment.

1. Intensive Distal

This mode clears in a proximal to distal motion, focusing intensively on the distal problem area. The sequence completes with a sequential clearance from distal to proximal.

2. Intensive Proximal

This mode works in proximal to distal motion, focusing intensively on the proximal problem area. The sequence completes with a sequential clearance from distal to proximal.

3. Intensive Whole

This mode works in a proximal to distal motion, intensively working on the whole of the limb. The sequence completes with a sequential clearance from distal to proximal.

4. Intensive Hand/Foot

This mode clears the lymphatic channels prior to working and focusing on the affected area. The sequence focuses on the hand or foot and completes with a sequential clearance from distal to proximal.

5. Sequential

This mode focuses on inflating each chamber in turn, distal to proximal, and then deflates all chambers at the same time.

6. Peristaltic

This mode provides a sequential distal to proximal inflation mode but with a wave throughout, this is achieved by the deflation of distal chambers while the more proximal inflate. This setting is sometimes preferred as there is less static pressure to the distal limb which may offer greater comfort.

7. Short Intensive Whole

This is a shorter intensive whole limb mode which is useful when treating bilateral limb oedema using the trouser garment.

LymphFlow Advance

Garments

Material Composition: Nylon, Thermoplastic Polyurethane


The LymphFlow Advance is a twelve-chamber pump that connects to a variety of separate garments, which can be used for the treatment of upper and lower limb.

Haddenham LymphFlow products utilise a pneumatic compression pumps paired with inflatable, multi-chambered garments to provide intermittent pressure, externally to an oedematous limb.



Lower Limb Sizing

FOOT SECTION measurements shown in blue



Circumferences cm		Legging	Extra Wide Legging	Thigh High + Waist	Trousers + Expanders	Extra Wide Trousers
t	Waist	-	-	-	≤119	≤134
g	Thigh Top	≤78	≤88	≤78	≤65	≤75
b	Ankle	≤52	≤78	≤58 ≤48	- ≤48	- ≤48
Straight Distance cm						
A-g	Leg Length	≤82	≤82	≤81	≤81	≤81

Codes



Legging



Thigh High + Waist



Trousers

	Legging		Thigh High + Waist		Trousers	
	STANDARD	EXTRA WIDE	LEFT LEG	RIGHT LEG	WITH EXPANDERS	EXTRA WIDE
LFA Garment	LFA5812	LFA5912	LFA8012-L	LFA8012-R	LFA7612	LFA7812
Expanders	EXP5812GZ	EXP5912GZ	EXP8012GZ	EXP8012GZ	-	EXP7812GZ

Seamless

Adaptive garments


LymphFlow Advance garments adapt to a patient’s limb size

Upper and lower limb garments


For targeted treatment of oedema



Upper Limb Sizing

	Circumferences cm		
		Armsleeve	Armsleeve + Shoulder
	g Axilla	≤58	≤58
	c Wrist	≤48	≤48
Straight Distance cm			
	c-g Arm Length	≤76	≤76

Codes

		
	Armsleeve	Armsleeve + Shoulder
LFA Garment	LFA3612	LFA3812
Expanders	EXP3612GZ	EXP3812GZ

Lymphscanner

PWC Measurement Unit

The LymphScanner is an all in one Percentage Water Content (PWC) measurement unit that is composed of an integrated probe, a built-in contact pressure sensor and a display screen.



Codes

Lymphscanner Module	Service
MDL2001	MDL SERVICE



Precision

Tracks treatment effectiveness


Helps evaluate for response and intervention

Non-invasive measurements

Scanning is painless and can be repeated

The device is used for the assessment and detection of lymphoedema.

It measures localized tissue water content and displays a PWC value.



LymphScanner automatically calculates the interim ratio of PWC measurements

Measuring the percentage water content (PWC) of the tissue is useful in early detection of lymphoedema and in follow-up of lymphoedema treatments.

PWC values vary by measurement site and thus a measurement method capable of rapid skin scanning helps lymphologists and therapists record spatial information on the regional extent of the patient's lymphoedema.





**Accessories
+ Bandages**



Accessories & Bandaging

Compression bandages, kinesiology tape, foam paddings and skin adhesive for assisting existing compression treatment.



**Clinistretch
Bandaging**

Page 162 ►




**Lymphpad
Foam**

Page 164 ►



**Spaghetti
Foam**

Page 165 ►



**Malleolus
Pads**

Page 165 ►



**Kiniseology
Tape**

Page 166 ►



**Klebi Skin
Adhesive**

Page 167 ►



**Special Fabric
Wash**

Page 167 ►



**Abdominal
Binder**

Page 168 ►

Clinistretch

Cohesive Bandage

Material Composition: Cotton, Polyamide, Polyurethane, Latex

Clinistretch is a low profile, conforming cohesive short-stretch bandage that is highly effective in the treatment of lymphoedema, venous leg ulceration and chronic venous insufficiency.



Clinistretch Cohesive is a short-stretch compression bandage that contains a micro-coating of latex.

Strong

Cohesive

The bandage only sticks to itself, for easy application

Cut to length

Trim to any desired length, for individualised care

Codes

Available on Drug Tariff, in single rolls

	4cm x 6m	6cm x 6m	8cm x 6m	10cm x 6m	12cm x 6m
	CC04 403-9814	CC06 403-9806	CC08 403-9798	CC10 403-9772	CC12 403-9780
Box Price:	£92.34	£102.06	£59.94	£126.36	£77.76
10% Discount	(Box of 36)	(Box of 36)	(Box of 18)	(Box of 36)	(Box of 18)

100% Cotton Bandage

Material Composition: Cotton

Clinistretch is a low profile, conforming cohesive short stretch bandage that is highly effective in the treatment of lymphoedema, venous leg ulceration and chronic venous insufficiency.



Clinistretch Short is a 100% cotton short-stretch compression bandage.

Soft

100% Cotton

Gentler option for sensitive skin

Reusable and washable

Multiple applications from one bandage

Codes

Available on **Drug Tariff**, in single rolls

	6cm x 5m	8cm x 5m	10cm x 5m	12cm x 5m
	CS06 403-9855	CS08 403-9848	CS10 403-9830	CS12 403-9822
Box Price:	£204.12	£239.76	£129.60	£155.52
10% Discount	(Box of 72)	(Box of 72)	(Box of 36)	(Box of 36)

Lymphpad

Foam Padding

Material Composition: Styrene Butadiene, Polyamide, Polyester

Haddenham Lymphpad –when used in conjunction with compression therapy– can help to break down fibrotic tissue by creating changes in tissue pressure and enhancing the massaging effect of treatment.



Lymphpad comes in either a dot or a wave form pattern, backed by a woven fabric.

Supportive

Increase lymphatic drainage

Supplements existing compression treatments

Break down fibrotic tissue

Softening skin for improved mobility and comfort

Codes

Dot Pattern

20.5cm x 29.5cm

LPD

Wave Pattern

20.5cm x 29.5cm

LPW

Price: £34.00
(Pack of 3)

£34.00
(Pack of 3)

Spaghetti Foam

Large economical sheet that can be cut to size.

1.3m x 1m

87000



Malleolus Pads

Soft silicone pad for use at either side of the ankle.

Left Side

425

Right Side

426



K-Active

Kinesiology Tape

K-Active kinesiology tape can be used to assist with lymphatic drainage by gently lifting the skin, helping congested fluid to move deeper into the lymphatic system.



K-Active is available in three different strength categories: Classic, Elite and Gentle. Elite is suitable for moist or oily skin, as well as those active during the day. Gentle is recommended for children, elderly patients and those with sensitive skin.

“The tape is very good. After applying the tape it will last in place for 3-4 days easily, allowing for showers etc. It’s a great help with the lymphoedema I have in my foot.”

Patient review, lymphshop.com

Codes

	Classic	Elite	Gentle
	5cm x 5m	5cm x 5m	5cm x 5m
 Beige	KA-050	KAE-050	KAG-050
 Black	KA-050BL	-	-
 Blue	KA-050B	-	-
 Pink	KA-050R	-	-
 Green	KA-050GR	-	-

Klebi Skin Adhesive

Ingredients: Isopropyl Alcohol, Glycerin, Water

Water soluble liquid adhesive used for maintaining the position of hosiery, orthopaedic aids, bandages and other textile articles.

410



Fabric Wash

Formulated to maintain fabrics in good condition.

412



Abdominal Binder

Material Composition: Polyamide, Viscose, Elastane


The Haddenham abdominal binder supports the healing process of postoperative wounds, providing high lateral stability.



Supportive

High lateral stability
Supports the spine using long longitudinal elastic compression.

Sizing



Circumference cm

Waist

S

70-90

M

90-110

L

109-130

XL

125-150

XXL

140-170

Options

Colour



White

Codes



39270
Unisex

Terms and conditions

Settlement	Nett 30 days. Accounts payable monthly. Interest on amounts outstanding over 30 days at the current bank rate + 3% on a daily basis.
Delivery	Standard charge of £6.00 ex VAT per delivery. Additional charge for courier service. Delivery times are approximate. Not responsible for courier's failure due to adverse weather or exceptional circumstances.
Prices	Subject to change without notice. Exclude VAT.
Returns	There is a 90 day time limit in which to return all products. The following are required for return authorisation: Products code(s), Reason for return, applicable invoice number(s) or order reference(s). Returns will not be accepted without a Returns Authorisation Number. Customer liability returns for Credit Note may be subject to a £5.00 handling charge per item. Returned goods must be saleable and have a minimum of 12 months shelf life remaining. Custom made items are non-returnable unless faulty.
Non-delivery	Non-receipt of goods must be notified in writing within 21 days of the date on the despatch note/invoice.
Warranty	Products are guaranteed against defects in material and manufacture. Once a garment is put into use, a 6-month warranty applies, covering material defects only. The warranty does not cover normal wear and tear, misuse, improper care, or use beyond the expiration date. Haddenham Healthcare reserves the right to replace or repair at its discretion. Specific warranty terms may vary per product. Consult product's 'Instructions for Use' for details.
Shelf life	Products have a shelf life of 3 years under correct storage conditions (unless otherwise specified).
Liability	Responsibility for selection and application of products rests with professional staff treating the patient.
Privacy	Privacy policy can be viewed at www.hadhealth.com/privacy-policy
Colours	Product colours depicted as accurately as printing process allows. Sample swatches available on request.

Contraindications

Absolute contraindications for lower limb

1. Advanced peripheral arterial occlusive disease, 2. Decompensated Heart Disease, 3. Septic phlebitis & Phlegmasia cerulea dolens, 4. Diabetic patients with lack or reduced sensation in lower limb or foot, 5. Supportive dermatoses, intolerance of fabric, 6. Sensory disturbances of the limb, advanced peripheral neuropathy

Relative contraindications for lower limb

1. Primary chronic arthritis, 2. Advanced cancer to the root of the limb

Absolute contraindications for upper limb

1. Super vena cava obstruction, 2. Axillary vein thrombosis, 3. Supportive dermatoses, intolerance to fabric

Relative contraindications for upper limb

1. Brachial plexus neuropathy, 2. Rheumatoid arthritis, 3. Diabetes (peripheral neuropathy), 4. Advanced cancer to the root of the limb, 5. Cellulitis

Additional contraindications apply to particular garments. Always consult product instructions prior to use. Haddenham Healthcare Ltd assumes no liability for contraindicated use of garments.

Contact Us

Call: 01844 208842

Email: sales@hadhealth.com

Prices effective from

September 2025

Issue 10