

VIALS + TUBES

DROPPER BOTTLE & ASSEMBLIES

CULTURE TUBES

DIAGNOSTIC VIALS

CHROMATOGRAPHY, HPLC & TOC VIALS

CRIMP TOP & HEADSPACE VIALS

SCINTILLATING VIALS

TUBULAR GLASS VIALS

DENTAL CARTIRADGE



**PRIMARY GLASS
PACKAGING SOLUTIONS**

GLASSCO LABORATORY EQUIPMENTS PVT. LTD.

www.glasscolabs.com 

www.glasscolabs.in 

enquiry@glasscolabs.com 

Glassco Group is a producer of high-quality Lab Glassware and Benchtop Instruments for more than **53 years** and exporting to **85+ countries**. We manufacture our products from 3 different plants. **Glassco Laboratory Equipment Pvt Ltd (AMBALA PLANT)** From this plant we manufacture and export a wide range of Scientific glassware (VOLUMETRIC, INTERCHANGEABLE, FILTRATION & GENERAL GLASSWARE) according to DIN, USP, ISO, and ASTM standards from one production facility.

In **2019**, we started manufacturing Disposable Glassware (AMBALA PLANT) like Test Tubes, Pasteur Pipettes, Culture Tubes, a milestone achieved in Glassco Laboratories, having a capacity of producing one million pieces per day.

Our Latest acquisition **ABHISHEK Glass (VAPI PLANT)** Gujrat has opened new doors for not only our existing customers but also for the **Glassco Group** is now manufacturing more than **19000+ SKU's** from **3 different plants** in India.

We are the biggest exporters from India holding a **market share of 55%** and now also the largest Pharma Glass Product Manufacturers, which makes us ; "TRULY A GLASS APART"

- **DROPPER BOTTLE & ASSEMBLIES**
- **CULTURE TUBES**
- **DIAGNOSTIC VIALS**
- **CHROMATOGRAPHY, HPLC & TOC VIALS**
- **CRIMP TOP & HEADSPACE VIALS**
- **SCINTILLATING VIALS**
- **TUBULAR GLASS VIALS**
- **DENTAL CARTRIDGE**

This catalogue covers almost all the range of vials and closures that are used in Pharma, Diagnostic, Cosmetic and Food & Beverages Industries. We also manufacture custom made sizes for the vials and droppers in case you don't find the size which fits your needs.



FOR THE GLASS DROPPERS CODES

USE GUIDE BELOW:

Cat. No. Starts with 619 for Dropper Tip, Followed by (.)

**TO GENERATE PART NUMBER
FILL ITEM CODES FROM EACH STEPS IN BLANKS :**

619 . .

ST3



ST1

For Example :
619.82.13

STEP 01 (ST1)

SELECTION OF DROPPER TIP

Types	Item Code
LP (Long Press Bent Tip)	80
STB (Sharp Tip Bent)	82
LPPT (Long Press Pencil Tip)	84
BT (Ball Tip)	86
SST (Straight Sharp Tip)	88
BBT (Bent Ball Tip)	90
BT (Blunt Tip)	92




STEP 03 (ST3)


SELECTION OF DIA X LENGTH

Dia x Length	Item Code
5.25x35	01
5.25x37	02
5.25x45	03
5.25x65	04
5.75x40	05
5.75x45	06
5.75x31	07
6.25x45	08
6.25x59	09
6.25x60	10
6.25x54	11
7.25x42	12
7.25x44	13
7.25x48	14
7.25x54	15
7.25x58	16
7.25x60	17
7.25x62	18
7.25x64	19
7.25x65	20
7.25x66	21
7.25x67	22
7.25x68	23
7.25x70	24
7.25x71	25
7.25x75	26
7.25x77	27
7.25x80	28
7.25x82	29
7.25x88	30
7.25x89	31
7.25x92	32
7.25x94	33
7.25x104	34
7.25x107	35
7.25x108	36
7.25x109	37
7.25x115	38


GLASS DROPPERS




LPBT



SBT




LPPT




BT

Item Code Long Press Bent	Item Code Sharp Bent	Item Code Long Press Pencil	Item Code Ball Tip	Dia (mm)	Height (mm)	Inner Size	Case Size
619.80.01	619.82.01	619.84.01	619.86.01	5.25	35	6000	24000
619.80.02	619.82.02	619.84.02	619.86.02	5.75	37	6000	24000
619.80.03	619.82.03	619.84.03	619.86.03	5.75	45	-	-
619.80.04	619.82.04	619.84.04	619.86.04	5.75	65	-	-
619.80.05	619.82.05	619.84.05	619.86.05	6.25	40	4500	18000
619.80.06	619.82.06	619.84.06	619.86.06	6.25	45	-	-
619.80.07	619.82.07	619.84.07	619.86.07	6.25	31	-	-
619.80.08	619.82.08	619.84.08	619.86.08	6.25	45	3000	12000
619.80.09	619.82.09	619.84.09	619.86.09	7.25	59	2100	8400
619.80.10	619.82.10	619.84.10	619.86.10	7.25	60	2100	8400
619.80.11	619.82.11	619.84.11	619.86.11	7.25	54		-
619.80.12	619.82.12	619.84.12	619.86.12	7.25	42	3000	12000
619.80.13	619.82.13	619.84.13	619.86.13	7.25	44	3000	12000
619.80.14	619.82.14	619.84.14	619.86.14	7.25	48	2100	8400
619.80.15	619.82.15	619.84.15	619.86.15	7.25	54	1950	7800
619.80.16	619.82.16	619.84.16	619.86.16	7.25	58	1950	7800
619.80.17	619.82.17	619.84.17	619.86.17	7.25	60	1800	7200
619.80.18	619.82.18	619.84.18	619.86.18	7.25	62	1800	7200
619.80.19	619.82.19	619.84.19	619.86.19	7.25	64	1800	7200
619.80.20	619.82.20	619.84.20	619.86.20	7.25	65	1500	6000
619.80.21	619.82.21	619.84.21	619.86.21	7.25	66	1500	6000
619.80.22	619.82.22	619.84.22	619.86.22	7.25	67	1500	6000
619.80.23	619.82.23	619.84.23	619.86.23	7.25	68	1500	6000
619.80.24	619.82.24	619.84.24	619.86.24	7.25	70	1500	6000
619.80.25	619.82.25	619.84.25	619.86.25	7.25	71	1500	6000
619.80.26	619.82.26	619.84.26	619.86.26	7.25	75	1350	5400
619.80.27	619.82.27	619.84.27	619.86.27	7.25	77	1350	5400
619.80.28	619.82.28	619.84.28	619.86.28	7.25	80	1350	5400
619.80.29	619.82.29	619.84.29	619.86.29	7.25	82	1200	4800
619.80.30	619.82.30	619.84.30	619.86.30	7.25	88	1200	4800
619.80.31	619.82.31	619.84.31	619.86.31	7.25	89	1200	4800
619.80.32	619.82.32	619.84.32	619.86.32	7.25	92	1050	4200
619.80.33	619.82.33	619.84.33	619.86.33	7.25	94	1050	4200
619.80.34	619.82.34	619.84.34	619.86.34	7.25	104	1050	4200
619.80.35	619.82.35	619.84.35	619.86.35	7.25	107	1050	4200
619.80.36	619.82.36	619.84.36	619.86.36	7.25	108	1050	4200
619.80.37	619.82.37	619.84.37	619.86.37	7.25	109	1050	4200
619.80.38	619.82.38	619.84.38	619.86.38	7.25	115	1050	4200


GLASS DROPPERS



SST



BBT



BT

Item Code Straight Sharp	Item Code Bent Ball	Item Code Blunt	Dia (mm)	Height (mm)	Inner Size	Case Size
619.88.01	619.90.01	619.92.01	5.25	35	6000	24000
619.88.02	619.90.02	619.92.02	5.75	37	6000	24000
619.88.03	619.90.03	619.92.03	5.75	45	-	-
619.88.04	619.90.04	619.92.04	5.75	65	-	-
619.88.05	619.90.05	619.92.05	6.25	40	4500	18000
619.88.06	619.90.06	619.92.06	6.25	45	-	-
619.88.07	619.90.07	619.92.07	6.25	31	-	-
619.88.08	619.90.08	619.92.08	6.25	45	3000	12000
619.88.09	619.90.09	619.92.09	7.25	59	2100	8400
619.88.10	619.90.10	619.92.10	7.25	60	2100	8400
619.88.11	619.90.11	619.92.11	7.25	54		-
619.88.12	619.90.12	619.92.12	7.25	42	3000	12000
619.88.13	619.90.13	619.92.13	7.25	44	3000	12000
619.88.14	619.90.14	619.92.14	7.25	48	2100	8400
619.88.15	619.90.15	619.92.15	7.25	54	1950	7800
619.88.16	619.90.16	619.92.16	7.25	58	1950	7800
619.88.17	619.90.17	619.92.17	7.25	60	1800	7200
619.88.18	619.90.18	619.92.18	7.25	62	1800	7200
619.88.19	619.90.19	619.92.19	7.25	64	1800	7200
619.88.20	619.90.20	619.92.20	7.25	65	1500	6000
619.88.21	619.90.21	619.92.21	7.25	66	1500	6000
619.88.22	619.90.22	619.92.22	7.25	67	1500	6000
619.88.23	619.90.23	619.92.23	7.25	68	1500	6000
619.88.24	619.90.24	619.92.24	7.25	70	1500	6000
619.88.25	619.90.25	619.92.25	7.25	71	1500	6000
619.88.26	619.90.26	619.92.26	7.25	75	1350	5400
619.88.27	619.90.27	619.92.27	7.25	77	1350	5400
619.88.28	619.90.28	619.92.28	7.25	80	1350	5400
619.88.29	619.90.29	619.92.29	7.25	82	1200	4800
619.88.30	619.90.30	619.92.30	7.25	88	1200	4800
619.88.31	619.90.31	619.92.31	7.25	89	1200	4800
619.88.32	619.90.32	619.92.32	7.25	92	1050	4200
619.88.33	619.90.33	619.92.33	7.25	94	1050	4200
619.88.34	619.90.34	619.92.34	7.25	104	1050	4200
619.88.35	619.90.35	619.92.35	7.25	107	1050	4200
619.88.36	619.90.36	619.92.36	7.25	108	1050	4200
619.88.37	619.90.37	619.92.37	7.25	109	1050	4200
619.88.38	619.90.38	619.92.38	7.25	115	1050	4200

DROPPER WITH TEAT & CAP

FOR THE **STANDARD DROPPER, CAP & TEAT ASSEMBLY CODES**
USE GUIDE BELOW:

279 . [] [] [] [] [] [] [] [] [] []

ST1 ST2 ST3 ST4

For Example : 279.80BK.01W

STEP 02 (ST2)

SELECTION OF TEAT COLOUR

Colour Name	Item Code
Black	BK
White	W
Grey	GY

STEP 03 (ST3)


SELECTION OF DIA X LENGTH

Dia x Length	Thread Size	Item Code
5.25x35	13	01
5.25x37	13	02
5.25x45	13	03
5.25x65	13	04
5.75x40	13	05
5.75x45	13	06
5.75x31	13	07
6.25x45	16	08
6.25x59	16	09
6.25x60	16	10
6.25x54	16	11
7.25x42	18	12
7.25x44	18	13
7.25x48	18	14
7.25x54	18	15
7.25x58	18	16
7.25x60	18	17
7.25x62	18	18
7.25x64	18	19
7.25x65	18	20
7.25x66	18	21
7.25x67	18	22
7.25x68	18	23
7.25x70	18	24
7.25x71	18	25
7.25x75	18	26
7.25x77	18	27
7.25x80	18	28
7.25x82	18	29
7.25x88	18	30
7.25x89	18	31
7.25x92	18	32
7.25x94	18	33
7.25x104	18	34
7.25x107	18	35
7.25x108	18	36
7.25x109	18	37
7.25x115	18	38

STEP 04 (ST4)

SELECTION OF CAP COLOUR

Colour	Colour Code
White	W
Blue	B
Red	R
Yellow	Y
Green	G
Black	BK
Grey	GY
Sky Blue	SB
Pink	P
Violet	V
Orange	O
Silver	S



STEP 01 (ST1)

SELECTION OF DROPPER TIP


Types	Item Code
LP (Long Press Bent Tip)	80
STB (Sharp Tip Bent)	82
LPPT (Long Press Pencil Tip)	84
BT (Ball Tip)	86
SST (Straight Sharp Tip)	88
BBT (Bent Ball Tip)	90
BT (Blunt Tip)	92

CREATE YOUR OWN

DROPPER



SCAN TO SELECT



Cat. No. Starts with 279 for Dropper Assembly Followed by (.)

Standard Dropper Assembly with Long Press Bent Tip (LPBT)



Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Dia (mm)	Height (mm)	Thread Size (mm)	Inner Size	Case Size
279.80BK.02W	279.80W.02W	279.80GY.02W	○ White	5.25	37	13mm	400	4800
279.80BK.03W	279.80W.03W	279.80GY.03W	○ White	5.25	45	13mm	400	4800
279.80BK.04W	279.80W.04W	279.80GY.04W	○ White	5.25	65	13mm	-	-
279.80BK.08W	279.80W.08W	279.80GY.08W	○ White	6.25	45	16mm	250	3000
279.80BK.09W	279.80W.09W	279.80GY.09W	○ White	6.25	59	16mm	220	2640
279.80BK.10W	279.80W.10W	279.80GY.10W	○ White	6.25	60	16mm	220	2640
279.80BK.02B	279.80W.02B	279.80GY.02B	● Blue	5.25	37	13mm	400	4800
279.80BK.03B	279.80W.03B	279.80GY.03B	● Blue	5.25	45	13mm	400	4800
279.80BK.04B	279.80W.04B	279.80GY.04B	● Blue	5.25	65	13mm	-	-
279.80BK.08B	279.80W.08B	279.80GY.08B	● Blue	6.25	45	16mm	250	3000
279.80BK.09B	279.80W.09B	279.80GY.09B	● Blue	6.25	59	16mm	220	2640
279.80BK.10B	279.80W.10B	279.80GY.10B	● Blue	6.25	60	16mm	220	2640
279.80BK.02R	279.80W.02R	279.80GY.02R	● Red	5.25	37	13mm	400	4800
279.80BK.03R	279.80W.03R	279.80GY.03R	● Red	5.25	45	13mm	400	4800
279.80BK.04R	279.80W.04R	279.80GY.04R	● Red	5.25	65	13mm	-	-
279.80BK.08R	279.80W.08R	279.80GY.08R	● Red	6.25	45	16mm	250	3000
279.80BK.09R	279.80W.09R	279.80GY.09R	● Red	6.25	59	16mm	220	2640
279.80BK.10R	279.80W.10R	279.80GY.10R	● Red	6.25	60	16mm	220	2640
279.80BK.02Y	279.80W.02Y	279.80GY.02Y	● Yellow	5.25	37	13mm	400	4800
279.80BK.03Y	279.80W.03Y	279.80GY.03Y	● Yellow	5.25	45	13mm	400	4800
279.80BK.04Y	279.80W.04Y	279.80GY.04Y	● Yellow	5.25	65	13mm	-	-
279.80BK.08Y	279.80W.08Y	279.80GY.08Y	● Yellow	6.25	45	16mm	250	3000
279.80BK.09Y	279.80W.09Y	279.80GY.09Y	● Yellow	6.25	59	16mm	220	2640
279.80BK.10Y	279.80W.10Y	279.80GY.10Y	● Yellow	6.25	60	16mm	220	2640
279.80BK.02G	279.80W.02G	279.80GY.02G	● Green	5.25	37	13mm	400	4800
279.80BK.03G	279.80W.03G	279.80GY.03G	● Green	5.25	45	13mm	400	4800
279.80BK.04G	279.80W.04G	279.80GY.04G	● Green	5.25	65	13mm	-	-
279.80BK.08G	279.80W.08G	279.80GY.08G	● Green	6.25	45	16mm	250	3000
279.80BK.09G	279.80W.09G	279.80GY.09G	● Green	6.25	59	16mm	220	2640
279.80BK.10G	279.80W.10G	279.80GY.10G	● Green	6.25	60	16mm	220	2640
279.80BK.02BK	279.80W.02BK	279.80GY.02BK	● Black	5.25	37	13mm	400	4800
279.80BK.03BK	279.80W.03BK	279.80GY.03BK	● Black	5.25	45	13mm	400	4800
279.80BK.04BK	279.80W.04BK	279.80GY.04BK	● Black	5.25	65	13mm	-	-
279.80BK.08BK	279.80W.08BK	279.80GY.08BK	● Black	6.25	45	16mm	250	3000
279.80BK.09BK	279.80W.09BK	279.80GY.09BK	● Black	6.25	59	16mm	220	2640
279.80BK.10BK	279.80W.10BK	279.80GY.10BK	● Black	6.25	60	16mm	220	2640
279.80BK.08SB	279.80W.08SB	279.80GY.08SB	● Sky Blue	6.25	45	16mm	250	3000
279.80BK.09SB	279.80W.09SB	279.80GY.09SB	● Sky Blue	6.25	59	16mm	220	2640
279.80BK.02GY	279.80W.02GY	279.80GY.02GY	● Grey	5.25	37	13mm	400	4800
279.80BK.03GY	279.80W.03GY	279.80GY.03GY	● Grey	5.25	45	13mm	400	4800
279.80BK.04GY	279.80W.04GY	279.80GY.04GY	● Grey	5.25	65	13mm	-	-
279.80BK.08GY	279.80W.08GY	279.80GY.08GY	● Grey	6.25	45	16mm	250	3000
279.80BK.09GY	279.80W.09GY	279.80GY.09GY	● Grey	6.25	59	16mm	220	2640
279.80BK.10GY	279.80W.10GY	279.80GY.10GY	● Grey	6.25	60	16mm	220	2640
279.80BK.08O	279.80W.08O	279.80GY.08O	● Orange	6.25	45	16mm	250	3000
279.80BK.09O	279.80W.09O	279.80GY.09O	● Orange	6.25	59	16mm	220	2640
279.80BK.10O	279.80W.10O	279.80GY.10O	● Orange	6.25	60	16mm	220	2640
279.80BK.10SB	279.80W.10SB	279.80GY.10SB	● Sky Blue	6.25	60	16mm	220	2640
279.80BK.08P	279.80W.08P	279.80GY.08P	● Pink	6.25	45	16mm	250	3000
279.80BK.09P	279.80W.09P	279.80GY.09P	● Pink	6.25	59	16mm	220	2640
279.80BK.10P	279.80W.10P	279.80GY.10P	● Pink	6.25	60	16mm	220	2640
279.80BK.08V	279.80W.08V	279.80GY.08V	● Violet	6.25	45	16mm	220	3000
279.80BK.09V	279.80W.09V	279.80GY.09V	● Violet	6.25	59	16mm	220	2640
279.80BK.10V	279.80W.10V	279.80GY.10V	● Violet	6.25	60	16mm	200	2640

Standard Dropper Assembly with Sharp Bent Tip (SBT)



Item Code (Black Text)	Item Code (White Text)	Item Code (Grey Text)	Cap Colour	Dia (mm)	Height (mm)	Thread Size (mm)	Inner Size	Case Size
279.82BK.01W	279.82W.01W	279.82GY.01W	○ White	5.25	35	13mm	1500	6000
279.82BK.02W	279.82W.02W	279.82GY.02W	○ White	5.25	37	13mm	400	4800
279.82BK.03W	279.82W.03W	279.82GY.03W	○ White	5.25	45	13mm	400	4800
279.82BK.04W	279.82W.04W	279.82GY.04W	○ White	5.25	65	13mm	-	-
279.82BK.08W	279.82W.08W	279.82GY.08W	○ White	6.25	45	16mm	250	3000
279.82BK.09W	279.82W.09W	279.82GY.09W	○ White	6.25	59	16mm	220	2640
279.82BK.02B	279.82W.02B	279.82GY.02B	● Blue	5.25	37	13mm	400	4800
279.82BK.03B	279.82W.03B	279.82GY.03B	● Blue	5.25	45	13mm	400	4800
279.82BK.04B	279.82W.04B	279.82GY.04B	● Blue	5.25	65	13mm	-	-
279.82BK.08B	279.82W.08B	279.82GY.08B	● Blue	6.25	45	16mm	250	3000
279.82BK.09B	279.82W.09B	279.82GY.09B	● Blue	6.25	59	16mm	220	2640
279.82BK.02R	279.82W.02R	279.82GY.02R	● Red	5.25	37	13mm	400	4800
279.82BK.03R	279.82W.03R	279.82GY.03R	● Red	5.25	45	13mm	400	4800
279.82BK.04R	279.82W.04R	279.82GY.04R	● Red	5.25	65	13mm	-	-
279.82BK.08R	279.82W.08R	279.82GY.08R	● Red	6.25	45	16mm	250	3000
279.82BK.09R	279.82W.09R	279.82GY.09R	● Red	6.25	59	16mm	220	2640
279.82BK.02Y	279.82W.02Y	279.82GY.02Y	● Yellow	5.25	37	13mm	400	4800
279.82BK.03Y	279.82W.03Y	279.82GY.03Y	● Yellow	5.25	45	13mm	400	4800
279.82BK.04Y	279.82W.04Y	279.82GY.04Y	● Yellow	5.25	65	13mm	-	-
279.82BK.08Y	279.82W.08Y	279.82GY.08Y	● Yellow	6.25	45	16mm	250	3000
279.82BK.09Y	279.82W.09Y	279.82GY.09Y	● Yellow	6.25	59	16mm	220	2640
279.82BK.02G	279.82W.02G	279.82GY.02G	● Green	5.25	37	13mm	400	4800
279.82BK.03G	279.82W.03G	279.82GY.03G	● Green	5.25	45	13mm	400	4800
279.82BK.04G	279.82W.04G	279.82GY.04G	● Green	5.25	65	13mm	-	-
279.82BK.08G	279.82W.08G	279.82GY.08G	● Green	6.25	45	16mm	250	3000
279.82BK.09G	279.82W.09G	279.82GY.09G	● Green	6.25	59	16mm	220	2640
279.82BK.02BK	279.82W.02BK	279.82GY.02BK	● Black	5.25	37	13mm	400	4800
279.82BK.03BK	279.82W.03BK	279.82GY.03BK	● Black	5.25	45	13mm	400	4800
279.82BK.04BK	279.82W.04BK	279.82GY.04BK	● Black	5.25	65	13mm	-	-
279.82BK.08BK	279.82W.08BK	279.82GY.08BK	● Black	6.25	45	16mm	250	3000
279.82BK.09BK	279.82W.09BK	279.82GY.09BK	● Black	6.25	59	16mm	220	2640
279.82BK.08SB	279.82W.08SB	279.82GY.08SB	● Sky Blue	6.25	45	16mm	250	3000
279.82BK.09SB	279.82W.09SB	279.82GY.09SB	● Sky Blue	6.25	59	16mm	220	2640
279.82BK.10SB	279.82W.10SB	279.82GY.10SB	● Sky Blue	6.25	60	16mm	220	2640
279.82BK.10P	279.82W.10P	279.82GY.10P	● Pink	6.25	60	16mm	220	2640
279.82BK.10V	279.82W.10V	279.82GY.10V	● Violet	6.25	60	16mm	220	2640
279.82BK.02GY	279.82W.02GY	279.82GY.02GY	● Grey	5.25	37	13mm	400	4800
279.82BK.03GY	279.82W.03GY	279.82GY.03GY	● Grey	5.25	45	13mm	400	4800
279.82BK.04GY	279.82W.04GY	279.82GY.04GY	● Grey	5.25	65	13mm	-	-
279.82BK.08GY	279.82W.08GY	279.82GY.08GY	● Grey	6.25	45	16mm	250	3000
279.82BK.09GY	279.82W.09GY	279.82GY.09GY	● Grey	6.25	59	16mm	220	2640
279.82BK.08O	279.82W.08O	279.82GY.08O	● Orange	6.25	45	16mm	220	2640
279.82BK.09O	279.82W.09O	279.82GY.09O	● Orange	6.25	59	16mm	220	2640
279.82BK.08V	279.82W.08V	279.82GY.08V	● Violet	6.25	45	16mm	250	3000
279.82BK.09V	279.82W.09V	279.82GY.09V	● Violet	6.25	59	16mm	220	2640
279.82BK.08P	279.82W.08P	279.82GY.08P	● Pink	6.25	45	16mm	250	3000
279.82BK.09P	279.82W.09P	279.82GY.09P	● Pink	6.25	59	16mm	220	2640

Standard Dropper Assembly with Long Press Pencil Tip (LPPT)



Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Dia (mm)	Height (mm)	Thread Size (mm)	Inner Size	Case Size
279.84BK.10W	279.84W.10W	279.84GY.10W	○ White	6.25	60	16mm	220	2640
279.84BK.10B	279.84W.10B	279.84GY.10B	● Blue	6.25	60	16mm	220	2640
279.84BK.10R	279.84W.10R	279.84GY.10R	● Red	6.25	60	16mm	220	2640
279.84BK.10Y	279.84W.10Y	279.84GY.10Y	● Yellow	6.25	60	16mm	220	2640
279.84BK.10G	279.84W.10G	279.84GY.10G	● Green	6.25	60	16mm	220	2640
279.84BK.10BK	279.84W.10BK	279.84GY.10BK	● Black	6.25	60	16mm	220	2640
279.84BK.10SB	279.84W.10SB	279.84GY.10SB	● Sky Blue	6.25	60	16mm	220	2640
279.84BK.10GY	279.84W.10GY	279.84GY.10GY	● Grey	6.25	60	16mm	220	2640
279.84BK.10V	279.84W.10V	279.84GY.10V	● Violet	6.25	60	16mm	220	2640
279.84BK.10P	279.84W.10P	279.84GY.10P	● Pink	6.25	60	16mm	220	2640
279.84BK.10O	279.84W.10O	279.84GY.10O	● Orange	6.25	60	16mm	220	2640

Standard Dropper Assembly with Ball Tip (BT)



Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Dia (mm)	Height (mm)	Thread Size (mm)	Inner Size	Case Size
279.86BK.16S	279.86W.16S	279.86GY.16S	○ Silver	7.25	58	18mm	600	1800
279.86BK.19S	279.86W.19S	279.86GY.19S	○ Silver	7.25	64	18mm	600	1800
279.86BK.27S	279.86W.27S	279.86GY.27S	○ Silver	7.25	77	18mm	500	1500
279.86BK.32S	279.86W.32S	279.86GY.32S	○ Silver	7.25	92	18mm	450	1350
279.86BK.31LG	279.86W.31LG	279.86GY.31LG	● Light Green	7.25	89	18mm	450	1350
279.86BK.16G	279.86W.16G	279.86GY.16G	● Green	7.25	58	18mm	600	1800
279.86BK.19G	279.86W.19G	279.86GY.19G	● Green	7.25	64	18mm	600	1800
279.86BK.27G	279.86W.27G	279.86GY.27G	● Green	7.25	77	18mm	500	1500
279.86BK.32G	279.86W.32G	279.86GY.32G	● Green	7.25	92	18mm	450	1350

EU Dropper Assembly Bottle Set

TO GENERATE PART NUMBER
FILL ITEM CODES FROM EACH STEPS IN BLANKS :

V279. **ST1** **ST2** **ST3** **ST4**
8 0 B K . 0 9 G

For Example : **V279.80BK.01W**

CREATE YOUR OWN

DROPPER




SCAN TO SELECT

FILL ITEM CODES FROM
EACH STEPS IN BLANKS :

**STEP
01**

SELECTION OF DROPPER TIP



Use Code :

V279

**STEP
02**

SELECTION OF TEAT COLOUR



Use Code :

**STEP
03**

SELECTION OF DROPPER BOTTLE



Bottle Capacity

Use Code :



**STEP
04**

SELECTION OF CAP COLOUR

Colour	Colour Code	Colour	Colour Code
White	W	Sky Blue	SB
Blue	B	Pink	P
Red	R	Violet	V
Yellow	Y	Orange	O
Green	G	Silver	S
Black	BK		
Grey	GY		



DROPPER WITH VIALS

EU DROPPER BOTTLES & ASSEMBLIES

Glassco Dropper Bottle includes an extensive collection of fully customizable ranges that are made from low expansion clear glass and are the perfect primary packaging solution for numerous uses in various industries from Personal Care, Pharmaceutical, Diagnostics to Cosmetic industry.

With State-of-the-art production facility we ensure consistent and best quality products are delivered to the buyer.


10 ML DROPPER BOTTLES

Code	Product Description	Bottle Size (D x H mm)	Packaging / Inner Packs	Case Quantity
086.800.12	10ml Clear Glass Dropper Bottle	21.25x59	120	1080
086.800.13*	10ml Clear Glass Dropper Bottle	21.25x60	120	1080
086.800.14*	10ml Clear Glass Dropper Bottle	21.25x62	120	1080
086.800.15*	10ml Clear Glass Dropper Bottle	21.25x63	120	1080


*Compatible Vials



Dropper Assemblies with Long Press Bent Tip (LPBT) for 10ml Bottles




Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.80BK.09W	279.80W.09W	279.80GY.09W	○ White	220	2640
279.80BK.09B	279.80W.09B	279.80GY.09B	● Blue	220	2640
279.80BK.09R	279.80W.09R	279.80GY.09R	● Red	220	2640
279.80BK.09Y	279.80W.09Y	279.80GY.09Y	● Yellow	220	2640
279.80BK.09G	279.80W.09G	279.80GY.09G	● Green	220	2640
279.80BK.09BK	279.80W.09BK	279.80GY.09BK	● Black	220	2640
279.80BK.09GY	279.80W.09GY	279.80GY.09GY	● Grey	220	2640
279.80BK.09SB	279.80W.09SB	279.80GY.09SB	● Sky Blue	220	2640
279.80BK.09P	279.80W.09P	279.80GY.09P	● Pink	220	2640
279.80BK.09V	279.80W.09V	279.80GY.09V	● Violet	220	2640
279.80BK.09O	279.80W.09O	279.80GY.09O	● Orange	220	2640
279.80BK.09S	279.80W.09S	279.80GY.09S	○ Silver	220	2640




Dropper Assemblies with Sharp Bent Tip (SBT) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.82BK.09W	279.82W.09W	279.82GY.09W	○ White	220	2640
279.82BK.09B	279.82W.09B	279.82GY.09B	● Blue	220	2640
279.82BK.09R	279.82W.09R	279.82GY.09R	● Red	220	2640
279.82BK.09Y	279.82W.09Y	279.82GY.09Y	● Yellow	220	2640
279.82BK.09G	279.82W.09G	279.82GY.09G	● Green	220	2640
279.82BK.09BK	279.82W.09BK	279.82GY.09BK	● Black	220	2640
279.82BK.09GY	279.82W.09GY	279.82GY.09GY	● Grey	220	2640
279.82BK.09SB	279.82W.09SB	279.82GY.09SB	● Sky Blue	220	2640
279.82BK.09P	279.82W.09P	279.82GY.09P	● Pink	220	2640
279.82BK.09V	279.82W.09V	279.82GY.09V	● Violet	220	2640
279.82BK.09O	279.82W.09O	279.82GY.09O	● Orange	220	2640
279.82BK.09S	279.82W.09S	279.82GY.09S	○ Silver	220	2640




Dropper Assemblies with Long Press Pencil Tip (LPPT) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.84BK.09W	279.84W.09W	279.84GY.09W	○ White	220	2640
279.84BK.09B	279.84W.09B	279.84GY.09B	● Blue	220	2640
279.84BK.09R	279.84W.09R	279.84GY.09R	● Red	220	2640
279.84BK.09Y	279.84W.09Y	279.84GY.09Y	● Yellow	220	2640
279.84BK.09G	279.84W.09G	279.84GY.09G	● Green	220	2640
279.84BK.09BK	279.84W.09BK	279.84GY.09BK	● Black	220	2640
279.84BK.09GY	279.84W.09GY	279.84GY.09GY	● Grey	220	2640
279.84BK.09SB	279.84W.09SB	279.84GY.09SB	● Sky Blue	220	2640
279.84BK.09P	279.84W.09P	279.84GY.09P	● Pink	220	2640
279.84BK.09V	279.84W.09V	279.84GY.09V	● Violet	220	2640
279.84BK.09O	279.84W.09O	279.84GY.09O	● Orange	220	2640
279.84BK.09S	279.84W.09S	279.84GY.09S	○ Silver	220	2640




Dropper Assemblies with Ball Tip (BT) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.86BK.09W	279.86W.09W	279.86GY.09W	○ White	220	2640
279.86BK.09B	279.86W.09B	279.86GY.09B	● Blue	220	2640
279.86BK.09R	279.86W.09R	279.86GY.09R	● Red	220	2640
279.86BK.09Y	279.86W.09Y	279.86GY.09Y	● Yellow	220	2640
279.86BK.09G	279.86W.09G	279.86GY.09G	● Green	220	2640
279.86BK.09BK	279.86W.09BK	279.86GY.09BK	● Black	220	2640
279.86BK.09GY	279.86W.09GY	279.86GY.09GY	● Grey	220	2640
279.86BK.09SB	279.86W.09SB	279.86GY.09SB	● Sky Blue	220	2640
279.86BK.09P	279.86W.09P	279.86GY.09P	● Pink	220	2640
279.86BK.09V	279.86W.09V	279.86GY.09V	● Violet	220	2640
279.86BK.09O	279.86W.09O	279.86GY.09O	● Orange	220	2640
279.86BK.09S	279.86W.09S	279.86GY.09S	○ Silver	220	2640




Dropper Assemblies with Straight Sharp Tip (SST) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.88BK.09W	279.88W.09W	279.88GY.09W	○ White	220	2640
279.88BK.09B	279.88W.09B	279.88GY.09B	● Blue	220	2640
279.88BK.09R	279.88W.09R	279.88GY.09R	● Red	220	2640
279.88BK.09Y	279.88W.09Y	279.88GY.09Y	● Yellow	220	2640
279.88BK.09G	279.88W.09G	279.88GY.09G	● Green	220	2640
279.88BK.09BK	279.88W.09BK	279.88GY.09BK	● Black	220	2640
279.88BK.09GY	279.88W.09GY	279.88GY.09GY	● Grey	220	2640
279.88BK.09SB	279.88W.09SB	279.88GY.09SB	● Sky Blue	220	2640
279.88BK.09P	279.88W.09P	279.88GY.09P	● Pink	220	2640
279.88BK.09V	279.88W.09V	279.88GY.09V	● Violet	220	2640
279.88BK.09O	279.88W.09O	279.88GY.09O	● Orange	220	2640
279.88BK.09S	279.88W.09S	279.88GY.09S	○ Silver	220	2640



Dropper Assemblies with Bent Ball Tip (BBT) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.90BK.09W	279.90W.09W	279.90GY.09W	○ White	220	2640
279.90BK.09B	279.90W.09B	279.90GY.09B	● Blue	220	2640
279.90BK.09R	279.90W.09R	279.90GY.09R	● Red	220	2640
279.90BK.09Y	279.90W.09Y	279.90GY.09Y	● Yellow	220	2640
279.90BK.09G	279.90W.09G	279.90GY.09G	● Green	220	2640
279.90BK.09BK	279.90W.09BK	279.90GY.09BK	● Black	220	2640
279.90BK.09GY	279.90W.09GY	279.90GY.09GY	● Grey	220	2640
279.90BK.09SB	279.90W.09SB	279.90GY.09SB	● Sky Blue	220	2640
279.90BK.09P	279.90W.09P	279.90GY.09P	● Pink	220	2640
279.90BK.09V	279.90W.09V	279.90GY.09V	● Violet	220	2640
279.90BK.09O	279.90W.09O	279.90GY.09O	● Orange	220	2640
279.90BK.09S	279.90W.09S	279.90GY.09S	○ Silver	220	2640




Dropper Assemblies with Blunt Tip (BT) for 10ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.92BK.09W	279.92W.09W	279.92GY.09W	○ White	220	2640
279.92BK.09B	279.92W.09B	279.92GY.09B	● Blue	220	2640
279.92BK.09R	279.92W.09R	279.92GY.09R	● Red	220	2640
279.92BK.09Y	279.92W.09Y	279.92GY.09Y	● Yellow	220	2640
279.92BK.09G	279.92W.09G	279.92GY.09G	● Green	220	2640
279.92BK.09BK	279.92W.09BK	279.92GY.09BK	● Black	220	2640
279.92BK.09GY	279.92W.09GY	279.92GY.09GY	● Grey	220	2640
279.92BK.09SB	279.92W.09SB	279.92GY.09SB	● Sky Blue	220	2640
279.92BK.09P	279.92W.09P	279.92GY.09P	● Pink	220	2640
279.92BK.09V	279.92W.09V	279.92GY.09V	● Violet	220	2640
279.92BK.09O	279.92W.09O	279.92GY.09O	● Orange	220	2640
279.92BK.09S	279.92W.09S	279.92GY.09S	○ Silver	220	2640


7.5 ML DROPPER BOTTLES

Code	Product Description	Bottle Size (D x H mm)	Packaging / Inner Packs	Case Quantity
086.800.10	7.5ml Clear Glass Dropper Bottle	21.25x45	120	1440




Dropper Assemblies with Long Press Bent Tip (LPBT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.80BK.08W	279.80W.08W	279.80GY.08W	○ White	250	3000
279.80BK.08B	279.80W.08B	279.80GY.08B	● Blue	250	3000
279.80BK.08R	279.80W.08R	279.80GY.08R	● Red	250	3000
279.80BK.08Y	279.80W.08Y	279.80GY.08Y	● Yellow	250	3000
279.80BK.08G	279.80W.08G	279.80GY.08G	● Green	250	3000
279.80BK.08BK	279.80W.08BK	279.80GY.08BK	● Black	250	3000
279.80BK.08GY	279.80W.08GY	279.80GY.08GY	● Grey	250	3000
279.80BK.08SB	279.80W.08SB	279.80GY.08SB	● Sky Blue	250	3000
279.80BK.08P	279.80W.08P	279.80GY.08P	● Pink	250	3000
279.80BK.08V	279.80W.08V	279.80GY.08V	● Violet	250	3000
279.80BK.08O	279.80W.08O	279.80GY.08O	● Orange	250	3000
279.80BK.08S	279.80W.08S	279.80GY.08S	○ Silver	250	3000




Dropper Assemblies with Sharp Bent Tip (SBT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.82BK.08W	279.82W.08W	279.82GY.08GY	○ White	250	3000
279.82BK.08B	279.82W.08B	279.82GY.08B	● Blue	250	3000
279.82BK.08R	279.82W.08R	279.82GY.08R	● Red	250	3000
279.82BK.08Y	279.82W.08Y	279.82GY.08Y	● Yellow	250	3000
279.82BK.08G	279.82W.08G	279.82GY.08G	● Green	250	3000
279.82BK.08BK	279.82W.08W	279.82GY.08GY	● Black	250	3000
279.82BK.08GY	279.82W.08GY	279.82GY.08GY	● Grey	250	3000
279.82BK.08SB	279.82W.08SB	279.82GY.08SB	● Sky Blue	250	3000
279.82BK.08P	279.82W.08P	279.82GY.08P	● Pink	250	3000
279.82BK.08V	279.82W.08V	279.82GY.08V	● Violet	250	3000
279.82BK.08O	279.82W.08O	279.82GY.08O	● Orange	250	3000
279.82BK.08S	279.82W.08S	279.82GY.08S	○ Silver	250	3000




Dropper Assemblies with Long Press Pencil Tip (LPPT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.84BK.08W	279.84W.08W	279.84GY.08GY	○ White	250	3000
279.84BK.08B	279.84W.08B	279.84GY.08B	● Blue	250	3000
279.84BK.08R	279.84W.08R	279.84GY.08R	● Red	250	3000
279.84BK.08Y	279.84W.08Y	279.84GY.08Y	● Yellow	250	3000
279.84BK.08G	279.84W.08G	279.84GY.08G	● Green	250	3000
279.84BK.08BK	279.84W.08W	279.84GY.08GY	● Black	250	3000
279.84BK.08GY	279.84W.08GY	279.84GY.08GY	● Grey	250	3000
279.84BK.08SB	279.84W.08SB	279.84GY.08SB	● Sky Blue	250	3000
279.84BK.08P	279.84W.08P	279.84GY.08P	● Pink	250	3000
279.84BK.08V	279.84W.08V	279.84GY.08V	● Violet	250	3000
279.84BK.08O	279.84W.08O	279.84GY.08O	● Orange	250	3000
279.84BK.08S	279.84W.08S	279.84GY.08S	○ Silver	250	3000




Dropper Assemblies with Ball Tip (BT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.86BK.08W	279.86W.08W	279.86GY.08W	○ White	250	3000
279.86BK.08B	279.86W.08B	279.86GY.08B	● Blue	250	3000
279.86BK.08R	279.86W.08R	279.86GY.08R	● Red	250	3000
279.86BK.08Y	279.86W.08Y	279.86GY.08Y	● Yellow	250	3000
279.86BK.08G	279.86W.08G	279.86GY.08G	● Green	250	3000
279.86BK.08BK	279.86W.08BK	279.86GY.08BK	● Black	250	3000
279.86BK.08GY	279.86W.08GY	279.86GY.08GY	● Grey	250	3000
279.86BK.08SB	279.86W.08SB	279.86GY.08SB	● Sky Blue	250	3000
279.86BK.08P	279.86W.08P	279.86GY.08P	● Pink	250	3000
279.86BK.08V	279.86W.08V	279.86GY.08V	● Violet	250	3000
279.86BK.08O	279.86W.08O	279.86GY.08O	● Orange	250	3000
279.86BK.08S	279.86W.08S	279.86GY.08S	○ Silver	250	3000




Dropper Assemblies with Straight Sharp Tip (SST) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.88BK.08W	279.88W.08W	279.88GY.08W	○ White	250	3000
279.88BK.08B	279.88W.08B	279.88GY.08B	● Blue	250	3000
279.88BK.08R	279.88W.08R	279.88GY.08R	● Red	250	3000
279.88BK.08Y	279.88W.08Y	279.88GY.08Y	● Yellow	250	3000
279.88BK.08G	279.88W.08G	279.88GY.08G	● Green	250	3000
279.88BK.08BK	279.88W.08BK	279.88GY.08BK	● Black	250	3000
279.88BK.08GY	279.88W.08GY	279.88GY.08GY	● Grey	250	3000
279.88BK.08SB	279.88W.08SB	279.88GY.08SB	● Sky Blue	250	3000
279.88BK.08P	279.88W.08P	279.88GY.08P	● Pink	250	3000
279.88BK.08V	279.88W.08V	279.88GY.08V	● Violet	250	3000
279.88BK.08O	279.88W.08O	279.88GY.08O	● Orange	250	3000
279.88BK.08S	279.88W.08S	279.88GY.08S	○ Silver	250	3000



Dropper Assemblies with Bent Ball Tip (BBT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.90BK.08W	279.90W.08W	279.90GY.08W	○ White	250	3000
279.90BK.08B	279.90W.08B	279.90GY.08B	● Blue	250	3000
279.90BK.08R	279.90W.08R	279.90GY.08R	● Red	250	3000
279.90BK.08Y	279.90W.08Y	279.90GY.08Y	● Yellow	250	3000
279.90BK.08G	279.90W.08G	279.90GY.08G	● Green	250	3000
279.90BK.08BK	279.90W.08BK	279.90GY.08BK	● Black	250	3000
279.90BK.08GY	279.90W.08GY	279.90GY.08GY	● Grey	250	3000
279.90BK.08SB	279.90W.08SB	279.90GY.08SB	● Sky Blue	250	3000
279.90BK.08P	279.90W.08P	279.90GY.08P	● Pink	250	3000
279.90BK.08V	279.90W.08V	279.90GY.08V	● Violet	250	3000
279.90BK.08O	279.90W.08O	279.90GY.08O	● Orange	250	3000
279.90BK.08S	279.90W.08S	279.90GY.08S	○ Silver	250	3000




Dropper Assemblies with Blunt Tip (BT) for 7.5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.92BK.08W	279.92W.08W	279.92GY.08W	○ White	250	3000
279.92BK.08B	279.92W.08B	279.92GY.08B	● Blue	250	3000
279.92BK.08R	279.92W.08R	279.92GY.08R	● Red	250	3000
279.92BK.08Y	279.92W.08Y	279.92GY.08Y	● Yellow	250	3000
279.92BK.08G	279.92W.08G	279.92GY.08G	● Green	250	3000
279.92BK.08BK	279.92W.08BK	279.92GY.08BK	● Black	250	3000
279.92BK.08GY	279.92W.08GY	279.92GY.08GY	● Grey	250	3000
279.92BK.08SB	279.92W.08SB	279.92GY.08SB	● Sky Blue	250	3000
279.92BK.08P	279.92W.08P	279.92GY.08P	● Pink	250	3000
279.92BK.08V	279.92W.08V	279.92GY.08V	● Violet	250	3000
279.92BK.08O	279.92W.08O	279.92GY.08O	● Orange	250	3000
279.92BK.08S	279.92W.08S	279.92GY.08S	○ Silver	250	3000


5 ML DROPPER BOTTLES

Code	Product Description	Bottle Size (D x H mm)	Packaging / Inner Packs	Case Quantity
086.800.06	5ml Clear Glass Dropper Bottle	16.25x45	209	2508




Dropper Assemblies with Long Press Bent Tip (LPBT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.80BK.06W	279.80W.06W	279.80GY.06W	○ White	250	3000
279.80BK.06B	279.80W.06B	279.80GY.06B	● Blue	250	3000
279.80BK.06R	279.80W.06R	279.80GY.06R	● Red	250	3000
279.80BK.06Y	279.80W.06Y	279.80GY.06Y	● Yellow	250	3000
279.80BK.06G	279.80W.06G	279.80GY.06G	● Green	250	3000
279.80BK.06BK	279.80W.06BK	279.80GY.06BK	● Black	250	3000
279.80BK.06GY	279.80W.06GY	279.80GY.06GY	● Grey	250	3000
279.80BK.06SB	279.80W.06SB	279.80GY.06SB	● Sky Blue	250	3000
279.80BK.06P	279.80W.06P	279.80GY.06P	● Pink	250	3000
279.80BK.06V	279.80W.06V	279.80GY.06V	● Violet	250	3000
279.80BK.06O	279.80W.06O	279.80GY.06O	● Orange	250	3000
279.80BK.06S	279.80W.06S	279.80GY.06S	○ Silver	250	3000




Dropper Assemblies with Sharp Bent Tip (SBT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.82BK.06W	279.82W.06W	279.82GY.06GY	○ White	250	3000
279.82BK.06B	279.82W.06B	279.82GY.06B	● Blue	250	3000
279.82BK.06R	279.82W.06R	279.82GY.06R	● Red	250	3000
279.82BK.06Y	279.82W.06Y	279.82GY.06Y	● Yellow	250	3000
279.82BK.06G	279.82W.06G	279.82GY.06G	● Green	250	3000
279.82BK.06BK	279.82W.06W	279.82GY.06GY	● Black	250	3000
279.82BK.06GY	279.82W.06GY	279.82GY.06GY	● Grey	250	3000
279.82BK.06SB	279.82W.06SB	279.82GY.06SB	● Sky Blue	250	3000
279.82BK.06P	279.82W.06P	279.82GY.06P	● Pink	250	3000
279.82BK.06V	279.82W.06V	279.82GY.06V	● Violet	250	3000
279.82BK.06O	279.82W.06O	279.82GY.06O	● Orange	250	3000
279.82BK.06S	279.82W.06S	279.82GY.06S	○ Silver	250	3000




Dropper Assemblies with Long Press Pencil Tip (LPPT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.84BK.06W	279.84W.06W	279.84GY.06GY	○ White	250	3000
279.84BK.06B	279.84W.06B	279.84GY.06B	● Blue	250	3000
279.84BK.06R	279.84W.06R	279.84GY.06R	● Red	250	3000
279.84BK.06Y	279.84W.06Y	279.84GY.06Y	● Yellow	250	3000
279.84BK.06G	279.84W.06G	279.84GY.06G	● Green	250	3000
279.84BK.06BK	279.84W.06W	279.84GY.06GY	● Black	250	3000
279.84BK.06GY	279.84W.06GY	279.84GY.06GY	● Grey	250	3000
279.84BK.06SB	279.84W.06SB	279.84GY.06SB	● Sky Blue	250	3000
279.84BK.06P	279.84W.06P	279.84GY.06P	● Pink	250	3000
279.84BK.06V	279.84W.06V	279.84GY.06V	● Violet	250	3000
279.84BK.06O	279.84W.06O	279.84GY.06O	● Orange	250	3000
279.84BK.06S	279.84W.06S	279.84GY.06S	○ Silver	250	3000




Dropper Assemblies with Ball Tip (BT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.86BK.06W	279.86W.06W	279.86GY.06W	○ White	250	3000
279.86BK.06B	279.86W.06B	279.86GY.06B	● Blue	250	3000
279.86BK.06R	279.86W.06R	279.86GY.06R	● Red	250	3000
279.86BK.06Y	279.86W.06Y	279.86GY.06Y	● Yellow	250	3000
279.86BK.06G	279.86W.06G	279.86GY.06G	● Green	250	3000
279.86BK.06BK	279.86W.06BK	279.86GY.06BK	● Black	250	3000
279.86BK.06GY	279.86W.06GY	279.86GY.06GY	● Grey	250	3000
279.86BK.06SB	279.86W.06SB	279.86GY.06SB	● Sky Blue	250	3000
279.86BK.06P	279.86W.06P	279.86GY.06P	● Pink	250	3000
279.86BK.06V	279.86W.06V	279.86GY.06V	● Violet	250	3000
279.86BK.06O	279.86W.06O	279.86GY.06O	● Orange	250	3000
279.86BK.06S	279.86W.06S	279.86GY.06S	○ Silver	250	3000




Dropper Assemblies with Straight Sharp Tip (SST) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.88BK.06W	279.88W.06W	279.88GY.06W	○ White	250	3000
279.88BK.06B	279.88W.06B	279.88GY.06B	● Blue	250	3000
279.88BK.06R	279.88W.06R	279.88GY.06R	● Red	250	3000
279.88BK.06Y	279.88W.06Y	279.88GY.06Y	● Yellow	250	3000
279.88BK.06G	279.88W.06G	279.88GY.06G	● Green	250	3000
279.88BK.06BK	279.88W.06BK	279.88GY.06BK	● Black	250	3000
279.88BK.06GY	279.88W.06GY	279.88GY.06GY	● Grey	250	3000
279.88BK.06SB	279.88W.06SB	279.88GY.06SB	● Sky Blue	250	3000
279.88BK.06P	279.88W.06P	279.88GY.06P	● Pink	250	3000
279.88BK.06V	279.88W.06V	279.88GY.06V	● Violet	250	3000
279.88BK.06O	279.88W.06O	279.88GY.06O	● Orange	250	3000
279.88BK.06S	279.88W.06S	279.88GY.06S	○ Silver	250	3000



Dropper Assemblies with Bent Ball Tip (BBT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.90BK.06W	279.90W.06W	279.90GY.06W	○ White	250	3000
279.90BK.06B	279.90W.06B	279.90GY.06B	● Blue	250	3000
279.90BK.06R	279.90W.06R	279.90GY.06R	● Red	250	3000
279.90BK.06Y	279.90W.06Y	279.90GY.06Y	● Yellow	250	3000
279.90BK.06G	279.90W.06G	279.90GY.06G	● Green	250	3000
279.90BK.06BK	279.90W.06BK	279.90GY.06BK	● Black	250	3000
279.90BK.06GY	279.90W.06GY	279.90GY.06GY	● Grey	250	3000
279.90BK.06SB	279.90W.06SB	279.90GY.06SB	● Sky Blue	250	3000
279.90BK.06P	279.90W.06P	279.90GY.06P	● Pink	250	3000
279.90BK.06V	279.90W.06V	279.90GY.06V	● Violet	250	3000
279.90BK.06O	279.90W.06O	279.90GY.06O	● Orange	250	3000
279.90BK.06S	279.90W.06S	279.90GY.06S	○ Silver	250	3000




Dropper Assemblies with Blunt Tip (BT) for 5ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.92BK.06W	279.92W.06W	279.92GY.06W	○ White	250	3000
279.92BK.06B	279.92W.06B	279.92GY.06B	● Blue	250	3000
279.92BK.06R	279.92W.06R	279.92GY.06R	● Red	250	3000
279.92BK.06Y	279.92W.06Y	279.92GY.06Y	● Yellow	250	3000
279.92BK.06G	279.92W.06G	279.92GY.06G	● Green	250	3000
279.92BK.06BK	279.92W.06BK	279.92GY.06BK	● Black	250	3000
279.92BK.06GY	279.92W.06GY	279.92GY.06GY	● Grey	250	3000
279.92BK.06SB	279.92W.06SB	279.92GY.06SB	● Sky Blue	250	3000
279.92BK.06P	279.92W.06P	279.92GY.06P	● Pink	250	3000
279.92BK.06V	279.92W.06V	279.92GY.06V	● Violet	250	3000
279.92BK.06O	279.92W.06O	279.92GY.06O	● Orange	250	3000
279.92BK.06S	279.92W.06S	279.92GY.06S	○ Silver	250	3000


3 ML DROPPER BOTTLES

Code	Product Description	Bottle Size (D x H mm)	Packaging / Inner Packs	Case Quantity
086.800.03	3ml Clear Glass Dropper Bottle	16.25x35	209	3135
086.800.04	3ml Clear Glass Dropper Bottle	16.25x37	209	3135
086.800.05	3ml Clear Glass Dropper Bottle	16.25x36	209	3135




Dropper Assemblies with Long Press Bent Tip (LPBT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.80BK.01W	279.80W.01W	279.80GY.01W	○ White	1500	6000
279.80BK.01B	279.80W.01B	279.80GY.01B	● Blue	1500	6000
279.80BK.01R	279.80W.01R	279.80GY.01R	● Red	1500	6000
279.80BK.01Y	279.80W.01Y	279.80GY.01Y	● Yellow	1500	6000
279.80BK.01G	279.80W.01G	279.80GY.01G	● Green	1500	6000
279.80BK.01BK	279.80W.01BK	279.80GY.01BK	● Black	1500	6000
279.80BK.01GY	279.80W.01GY	279.80GY.01GY	● Grey	1500	6000
279.80BK.01SB	279.80W.01SB	279.80GY.01SB	● Sky Blue	1500	6000
279.80BK.01P	279.80W.01P	279.80GY.01P	● Pink	1500	6000
279.80BK.01V	279.80W.01V	279.80GY.01V	● Violet	1500	6000
279.80BK.01O	279.80W.01O	279.80GY.01O	● Orange	1500	6000
279.80BK.01S	279.80W.01S	279.80GY.01S	○ Silver	1500	6000




Dropper Assemblies with Sharp Bent Tip (SBT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.82BK.01W	279.82W.01W	279.82GY.01GY	○ White	1500	6000
279.82BK.01B	279.82W.01B	279.82GY.01B	● Blue	1500	6000
279.82BK.01R	279.82W.01R	279.82GY.01R	● Red	1500	6000
279.82BK.01Y	279.82W.01Y	279.82GY.01Y	● Yellow	1500	6000
279.82BK.01G	279.82W.01G	279.82GY.01G	● Green	1500	6000
279.82BK.01BK	279.82W.01W	279.82GY.01GY	● Black	1500	6000
279.82BK.01GY	279.82W.01GY	279.82GY.01GY	● Grey	1500	6000
279.82BK.01SB	279.82W.01SB	279.82GY.01SB	● Sky Blue	1500	6000
279.82BK.01P	279.82W.01P	279.82GY.01P	● Pink	1500	6000
279.82BK.01V	279.82W.01V	279.82GY.01V	● Violet	1500	6000
279.82BK.01O	279.82W.01O	279.82GY.01O	● Orange	1500	6000
279.82BK.01S	279.82W.01S	279.82GY.01S	○ Silver	1500	6000




Dropper Assemblies with Long Press Pencil Tip (LPPT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.84BK.01W	279.84W.01W	279.84GY.01GY	○ White	1500	6000
279.84BK.01B	279.84W.01B	279.84GY.01B	● Blue	1500	6000
279.84BK.01R	279.84W.01R	279.84GY.01R	● Red	1500	6000
279.84BK.01Y	279.84W.01Y	279.84GY.01Y	● Yellow	1500	6000
279.84BK.01G	279.84W.01G	279.84GY.01G	● Green	1500	6000
279.84BK.01BK	279.84W.01W	279.84GY.01GY	● Black	1500	6000
279.84BK.01GY	279.84W.01GY	279.84GY.01GY	● Grey	1500	6000
279.84BK.01SB	279.84W.01SB	279.84GY.01SB	● Sky Blue	1500	6000
279.84BK.01P	279.84W.01P	279.84GY.01P	● Pink	1500	6000
279.84BK.01V	279.84W.01V	279.84GY.01V	● Violet	1500	6000
279.84BK.01O	279.84W.01O	279.84GY.01O	● Orange	1500	6000
279.84BK.01S	279.84W.01S	279.84GY.01S	○ Silver	1500	6000




Dropper Assemblies with Ball Tip (BT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.86BK.01W	279.86W.01W	279.86GY.01W	○ White	1500	6000
279.86BK.01B	279.86W.01B	279.86GY.01B	● Blue	1500	6000
279.86BK.01R	279.86W.01R	279.86GY.01R	● Red	1500	6000
279.86BK.01Y	279.86W.01Y	279.86GY.01Y	● Yellow	1500	6000
279.86BK.01G	279.86W.01G	279.86GY.01G	● Green	1500	6000
279.86BK.01BK	279.86W.01BK	279.86GY.01BK	● Black	1500	6000
279.86BK.01GY	279.86W.01GY	279.86GY.01GY	● Grey	1500	6000
279.86BK.01SB	279.86W.01SB	279.86GY.01SB	● Sky Blue	1500	6000
279.86BK.01P	279.86W.01P	279.86GY.01P	● Pink	1500	6000
279.86BK.01V	279.86W.01V	279.86GY.01V	● Violet	1500	6000
279.86BK.01O	279.86W.01O	279.86GY.01O	● Orange	1500	6000
279.86BK.01S	279.86W.01S	279.86GY.01S	○ Silver	1500	6000




Dropper Assemblies with Straight Sharp Tip (SST) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.88BK.01W	279.88W.01W	279.88GY.01W	○ White	1500	6000
279.88BK.01B	279.88W.01B	279.88GY.01B	● Blue	1500	6000
279.88BK.01R	279.88W.01R	279.88GY.01R	● Red	1500	6000
279.88BK.01Y	279.88W.01Y	279.88GY.01Y	● Yellow	1500	6000
279.88BK.01G	279.88W.01G	279.88GY.01G	● Green	1500	6000
279.88BK.01BK	279.88W.01BK	279.88GY.01BK	● Black	1500	6000
279.88BK.01GY	279.88W.01GY	279.88GY.01GY	● Grey	1500	6000
279.88BK.01SB	279.88W.01SB	279.88GY.01SB	● Sky Blue	1500	6000
279.88BK.01P	279.88W.01P	279.88GY.01P	● Pink	1500	6000
279.88BK.01V	279.88W.01V	279.88GY.01V	● Violet	1500	6000
279.88BK.01O	279.88W.01O	279.88GY.01O	● Orange	1500	6000
279.88BK.01S	279.88W.01S	279.88GY.01S	○ Silver	1500	6000



Dropper Assemblies with Bent Ball Tip (BBT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.90BK.01W	279.90W.01W	279.90GY.01W	○ White	1500	6000
279.90BK.01B	279.90W.01B	279.90GY.01B	● Blue	1500	6000
279.90BK.01R	279.90W.01R	279.90GY.01R	● Red	1500	6000
279.90BK.01Y	279.90W.01Y	279.90GY.01Y	● Yellow	1500	6000
279.90BK.01G	279.90W.01G	279.90GY.01G	● Green	1500	6000
279.90BK.01BK	279.90W.01BK	279.90GY.01BK	● Black	1500	6000
279.90BK.01GY	279.90W.01GY	279.90GY.01GY	● Grey	1500	6000
279.90BK.01SB	279.90W.01SB	279.90GY.01SB	● Sky Blue	1500	6000
279.90BK.01P	279.90W.01P	279.90GY.01P	● Pink	1500	6000
279.90BK.01V	279.90W.01V	279.90GY.01V	● Violet	1500	6000
279.90BK.01O	279.90W.01O	279.90GY.01O	● Orange	1500	6000
279.90BK.01S	279.90W.01S	279.90GY.01S	○ Silver	1500	6000



Dropper Assemblies with Blunt Tip (BT) for 3ml Bottles

Item Code (Black Teat)	Item Code (White Teat)	Item Code (Grey Teat)	Cap Colour	Pack Size	Case Size
279.92BK.01W	279.92W.01W	279.92GY.01W	○ White	1500	6000
279.92BK.01B	279.92W.01B	279.92GY.01B	● Blue	1500	6000
279.92BK.01R	279.92W.01R	279.92GY.01R	● Red	1500	6000
279.92BK.01Y	279.92W.01Y	279.92GY.01Y	● Yellow	1500	6000
279.92BK.01G	279.92W.01G	279.92GY.01G	● Green	1500	6000
279.92BK.01BK	279.92W.01BK	279.92GY.01BK	● Black	1500	6000
279.92BK.01GY	279.92W.01GY	279.92GY.01GY	● Grey	1500	6000
279.92BK.01SB	279.92W.01SB	279.92GY.01SB	● Sky Blue	1500	6000
279.92BK.01P	279.92W.01P	279.92GY.01P	● Pink	1500	6000
279.92BK.01V	279.92W.01V	279.92GY.01V	● Violet	1500	6000
279.92BK.01O	279.92W.01O	279.92GY.01O	● Orange	1500	6000
279.92BK.01S	279.92W.01S	279.92GY.01S	○ Silver	1500	6000



DISPOSABLE CULTURE TUBE

Disposable Culture Tube - Boro Glass 5.1

Glassco Disposable Culture Tubes are resistant to chemical attack and weathering.

Salient Features:

- Borosilicate glass compliance to ASTM E-438, TYPE-1 Standard
- Offers excellent chemical resistance
- High heat tolerance
- Uniform dimensions
- Convenient dispensing packages

Manufactured from Type-I, Class B borosilicate glass, each tube is heat-treated at temperature exceeding 500°C (932°F) immediately before packing, which eliminates contamination, increases strength, and reduces breakage. Tubes feature well-rounded bottoms and fire polished tops free of rough edges. Dimensional specifications of these tubes are rigidly controlled for reliable, uniform results. Tubes are clean-sealed in a convenient dispensing package to keep them free from dust and other particulate contamination. Tube (O.D. 6 x 50 mm) is also known as a Durham Tube. Contents will not leach out and contaminate solutions. Borosilicate glass construction prevents pH change in laboratory tests.

US Sizes

Code	Size (mm)	Over Flow Volume (ml)	Inner Pack Size	Case Size
27729.566	6x50	0.75	100	2000
27729.568	10x75	3	250	1000
27729.570	12x75	5	250	1000
27729.572	13x100	9	250	1000
27729.574	15x85	10	250	1000
27729.576	16x100	14	250	1000
27729.578	16x125	18	250	1000
27729.580	16x150	22	250	1000
27729.583	18x150	28	125	500
27729.584	20x150	34	250	500
27729.586	25x150	60	125	500



Test Tube Boro Glass 5.1 Expansion

- With Rim
- EU Sizes

Code	Size (mm)	Over Flow Volume (ml)	Inner Pack Size	Case Size
27729.667	8x70	0.4-0.5	100	2000
27729.669	10x100	0.4-0.5	100	1000
27729.671	12x100	0.4-0.5	100	1000
27729.682	16x160	0.5-0.6	100	1000
27729.683	18x180	0.5-0.6	100	400
27729.685	20x180	0.5-0.6	100	400
27729.668	10x75	0.4-0.5	100	1000
27729.670	12x75	0.4-0.5	100	1000
27729.673	14x130	0.4-0.5	100	1000
27729.679	16x130	0.5-0.6	100	1000
27729.686	25x150	0.6-0.7	50	500
27729.688	25x200	0.6-0.7	50	500
27729.692	30x200	0.7-0.8	50	300



Disposable Culture Tube - Boro Glass 7.0

These tubes are chemically clean and well-suited for general use in chemistry, blood banks, hematology and other laboratory procedures.

Salient Features:

- Fire-polished top
- Shrink-film sealed to prevent contamination
- Well-rounded bottoms

Each tube is annealed at temperatures over 426.7°C (800°F) to increase strength and reduce breakage. Tubes feature smooth, well-rounded bottoms with no irregularities. Tops are fire-polished to eliminate rough edges. Tubes are packed in lint-free packaging and are shrink-film sealed to prevent any contamination from outside sources.

US Sizes

Code	Size (mm)	Over Flow Volume (ml)	Inner Pack Size	Case Size
20825.402	10 x 75	3	250	1000
20825.406	12 x 75	5	250	1000
20825.414	13 x 100	9	250	1000
20825.425	16 x 100	11	250	1000
20825.430	16 x 125	15	250	1000
20825.435	16 x 150	20	250	1000
20825.443	18 x 150	28	125	500
20825.447	20 x 150	34	125	500



Disposable Culture Tube without rim

- Boro 7.0 Expansion
- EU Sizes

Code	Size (mm)	Wall Thickness Volume (mm)	Inner Pack Size	Case Size
20825.390	8 x40	0.8-1	250	2000
20825.400	10x70	0.8-1	250	1000
20825.406A	12x75	0.8-1	250	1000
20825.403	10x100	0.8-1	250	1000
20825.408	12x100	0.8-1	250	1000
20825.425A	16x100	0.8-1	250	1000
20825.418	14x100	0.8-1	250	1000
20825.420	14x130	0.8-1	250	1000
20825.438	16x160	0.8-1	250	1000
20825.444	18x180	0.8-1	250	500
20825.402	10x75	0.5-0.6	250	1000
20825.406	12x75	0.5-0.6	250	1000
20825.408A	12x100	0.5-0.6	250	1000
20825.414	13x100	0.5-0.6	250	1000
20825.418A	14x100	0.5-0.6	250	1000
20825.425	16x100	0.6-0.7	250	1000
20825.430	16x125	0.6-0.7	250	1000
20825.435	16x150	0.6-0.7	250	1000
20825.438A	16x160	0.6-0.7	250	1000



INNER SIZE



CASE SIZE



PACKAGING

INNER SIZE



CASE SIZE



Disposable Culture Tube with Screw Cap Finish Boro Glass 5.1

Tubes have a marking spot and reference mark for orientation during use and a threaded neck to accept screw cap.

Salient Features :

- Borosilicate glass 5.1 expansion construction
- Marking spot and orientation reference mark on all tubes

Code	O.D. x Height (mm)	Thread Finish	Over Flow Volume (ml)	Case Size
23283.800	13 x 100	13-415	8	1000
23283.802	16 x 100	15-415	12	1000
23283.804	16 x 125	15-415	16	1000
23283.806	16 x 150	15-415	20	1000
23283.685	20 x 125	18-415	30	500
23283.686	20 x 150	18-415	33	500



Ps: Any other sizes can be made on request

Phenolic Caps, With Bonded Rubber Liner

• For Disposable Culture Tube

Code	O.D. x Height (mm)	Thread Finish	Case Size
23283.800C	13 x 100	13-415	1000
23283.804C	16 x 125	15-415	1000
23283.685C	20 x 125	18-415	1000



Screw Neck Tubes and Culture Media Tubes

The screw neck tubes are made out of clear Borosilicate glass can be printed with ceramic paint. In this page you will find our standard dimensions. It is however also possible to produce tubes according to your specifications.

Code	Dia (mm)	Length (mm)	Wall Thickness	Nominal volume (ml)	Bottom (mm)	Screw Thread	Case Size
23284.500	16	50	0.95	5	Flat	Capalu 16	480
23284.501	16	75	0.95	8	Flat	Capalu 16	480
23284.502	16	112.5	0.95	15	Flat	Capalu 16	299
23284.503	16	125	0.95	18	Flat	Capalu 16	299
23284.504	16	160	0.95	20	Round	Capalu 16	299
23284.505	18	125	0.95	22	Flat	Capalu 16	252
23284.506	18	145	1.05	24	Flat	Capalu 16	252
23284.507	18	180	1.05	31	Flat	Capalu 16	252
23284.508	20	200	1.05	42	Round	Capalu 16	215



Screw Caps

Our standard screw neck tubes are produced with the screw thread Capalu 16.

For this thread we have aluminium screw caps available. These caps come with Silicone Liner

Code	Screw Thread	Material	Colour	Inlay	Qty per box	Qty per Pallet
23284.CA	Capalu 16	Aluminium	Gold	Rubber	5.000	150.000



DISPOSABLE PASTEUR PIPETTE

Pasteur Pipette

Glassco Pasteur pipettes are manufactured to conform strict specifications laid down by research and industrial laboratories and are suitable for various laboratory applications.

Salient features

- Made up of Borosilicate Glass
- Facilitate Controlled Drop
- Smooth, Straight-cut, Fire-polished tips (Short & Long Tip)
- Ideal for Clinical & Industrial use
- Disposable
- Supplied in Convenient Packing

Item Code	Product Description	Length (cm)	Case Size
24672.200	Boro 7.0 Expansion	14.6	1000
24672.608	Boro 7.0 Expansion with Controlled Drop**	14.6	1000
24672.380	Boro 7.0 Expansion	22.9	1000
24673.010	Boro (3.3) Pipette - Short Tip (2 ml)	14.6	1000
23283.916	Boro 3.3 Pipette - Short Tip, Controlled Drop (2 ml)	14.6	1000
24673.043	Boro (3.3) Pipette - Long Tip (2 ml)	22.9	1000

****Pipette body O.D. : 7 mm. | Capillary stem O.D. : Approx. 1.5mm, Delivers 22 to 28 Drops per ml**



Glassco Pasteur Pipettes with Cotton Plugs

Glassco Pasteur pipettes are manufactured to conform strict specifications laid down by research and industrial laboratories and are suitable for various laboratory applications.

Salient Features:

- Fire-polished top
- Made up of Borosilicate Glass
- Smooth, Straight-cut, Fire-polished tips (Short & Long Tip)
- Ideal for Clinical & Industrial use
- Disposable
- Supplied in Convenient Packing

Item Code	Product Description	Length (cm)	Case Size
24672.400	Boro 7.0 Expansion	14.6	1000
24672.405	Boro7.0 Expansion	22.9	1000
24672.410	Boro 3.3	14.6	1000
24672.412	Boro 3.3	22.9	1000



DIAGNOSTIC VIALS

Manufactured from Low Expansion Borosilicate TYPE-1 Glass and controlled bore ensures excellent fit with rubber stoppers. Screw neck vials are used in the pharmaceutical, diagnostic and food industries.

- Perfect for in-vitro diagnostic procedures
- Vials with screw thread design eliminates need for crimping tools
- Vial available in clear and amber color
- Available from stock in convenient pack sizes
- Suitable for autoclave testing machines
- Caps available in Polypropylene or Urea material

DIAGNOSTIC VIALS

Code	Capacity (ml)	Vial DxH (mm)	Neck Size (mm)	Pack Size
GW.DVCG.01	3	18x36	15-425	320
GW.DVAG.01	3	18x36	15-425	320
GW.DVCG.02	5	18x41	15-425	320
GW.DVAG.02	5	18x41	15-425	320
GW.DVCG.03	10	25x50	22-400	154
GW.DVAG.03	10	25x50	22-400	154
GW.DVCG.04	20	27.50x63	22-400	130
GW.DVAG.04	20	27.50x63	22-400	130
GW.DVCG.05	25	30x63	22-400	99
GW.DVAG.05	25	30x63	22-400	99

RUBBER STOPPERS

Code	Glass Type	Stopper Material	Type	Pack Qty
GW.DVRS.01	13mm Stopper	Rubber	Igloo Bromobutyl Slotted	1000
GW.DVRS.02	20mm Stopper	Rubber	Igloo Bromobutyl Slotted	1000



RUBBER STOPPERS

Code	Glass Type	Stopper Material	Type	Pack Qty
GW.DVRS.01	13mm Stopper	Rubber	Igloo Bromobutyl Slotted (RFU)	1000
GW.DVRS.02	20mm Stopper	Rubber	Igloo Bromobutyl Slotted (RFU)	1000



CAPS

Code	Cap Colour	Capacity Material	Liner	Neck Size (ml)	Pack Qty
GW.DVWCPP.01	White	PP	PE	15.425	1000
GW.DVBCPP.02	Blue	PP	PE	15.425	1000
GW.DVGCPP.03	Green	PP	PE	15.425	1000
GW.DVRCPP.04	Red	PP	PE	15.425	1000
GW.DVYCPP.05	Yellow	PP	PE	15.425	1000
GW.DVWCU.01	White	Urea	PE	15.425	1000
GW.DVWCPP.06	White	PP	PE	22.400	1000
GW.DVBCPP.07	Blue	PP	PE	22.400	1000
GW.DVGCPP.08	Green	PP	PE	22.400	1000
GW.DVRCPP.09	Red	PP	PE	22.400	1000
GW.DVYCPP.10	Yellow	PP	PE	22.400	1000
GW.DVWCU.06	White	Urea	PE	22.400	1000



CHROMATOGRAPHY VIALS

SCREW NEW VIALS

Vial are made from ASTM 438 TYPE-1 Class B Borosilicate Glass that are crafted to fit on almost all autosamplers that ensures long-term sample storage when used with top closure caps.

Available in Two Sizes 12x32 mm (ND8/ND9) & 15x45 mm (ND13)



Screw New Vials 12 X 32 mm (ND 8)

Screw neck vials ND8 are available in clear and amber Type 1 Glass. vials comes with 8-425 thread for used as standard GC & HPLC applications.

Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Pack Unit
GW.SVC8.01	Clear	-	ND8	11.6 x 32	1.5	100
GW.SVC8.02	Clear	Yes	ND8	11.6 x 32	1.5	100
GW.SVA8.01	Amber	-	ND8	11.6 x 32	1.5	100
GW.SVA8.02	Amber	Yes	ND8	11.6 x 32	1.5	100

PP Caps for Screw Neck Vials 12 x 32 mm (ND 8)

Screw Caps ND8 are made of Polypropylene - PP with fitted septa. Caps have 8-425 thread, with open top centre hole

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.SV8.BC1	Black	Open Top	ND8	Silicone White / PTFE red	100
GW.SV8.BC2	White	Open Top	ND8	Silicone Red / PTFE White	100



Screw New Vials 12 X 32 mm (ND 9)

Screw neck vials ND9 are available in clear and amber Type 1 Glass. vials come with wide opening 9mm thread and can be used on almost any autosampler



Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Pack Unit
GW.SVC9.01	Clear	-	ND9	11.6 x 32	1.5	100
GW.SVC9.02	Clear	Yes	ND9	11.6 x 32	1.5	100
GW.SVA9.01	Amber	-	ND9	11.6 x 32	1.5	100
GW.SVA9.02	Amber	Yes	Nd9	11.6 x 32	1.5	100

PP Caps for Screw Neck Vials 12 x 32 mm (ND 9)

Screw Caps ND9 are made of Polypropylene - PP with fitted septa. Caps have 9mm thread, with open top centre hole

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.SV9.BLC2	Blue	Open Top	ND9	PTFE Red / White Silicone	100
GW.SV9.BLC4	Blue	Open Top	ND9	PTFE Blue/ White Silicone/ pre-slit (single)	100
GW.SV9.BLC5	Blue	Open Top	ND9	PTFE Red / Silicone White pre-slit	100



Screw New Vials 15 x 45 mm (ND 13)

Screw neck vials ND13 are available in clear and amber Type 1 Glass. vials comes with 13-425 thread wide opening for easy filling of viscous liquids.



Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Pack Unit
GW.SVC13.01	Clear	-	ND13	14.7 x 45	4	100
GW.SVC13.02	Clear	Yes	ND13	14.7 x 45	4	100
GW.SVA13.01	Amber	-	ND13	14.7 x 45	4	100
GW.SVA13.02	Amber	Yes	ND13	14.7 x 45	4	100

PP Caps for Screw Neck Vials 15 x 45 mm (ND 13)

Screw Caps ND13 are made of Polypropylene - PP with fitted septa. Caps have 13-425 thread, with open top centre hole.

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.SV13.BC1	Black	Open Top	ND13	Silicone White / PTFE Red pre-slit	100
GW.SV13.BC2	Black	Open Top	ND13	Silicone Red / PTFE White	100
GW.SV13.BC4	Black	Open Top	Nd13	Silicone White / PTFE Red	100




CRIMP TOP Vials

- Designed to use for High Pressure / Performance Liquid Chromatography (HPLC) and Gas Chromatography (GC) Applications.
- Crimp top vials consist of an Aluminum cap and PTFE / Silicon / Rubber Septum thus forms a superior seal and reduces evaporation.
- Optional with White Ceramic Writing Patch and fill marks for easy labelling.

Crimp Top Vials 12 x 32 mm (ND 11)

Crimp Neck vials ND11 are available in clear and amber Type 1 Glass for used as standard GC & HPLC applications.

Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Pack Unit
GW.CTC11.01	Clear	-	ND11	11.6 x 32	1.5	100
GW.CTC11.02	Clear	Yes	ND11	11.6 x 32	1.5	100
GW.CTA11.01	Amber	-	ND11	11.6 x 32	1.5	100
GW.CTA11.02	Amber	Yes	ND11	11.6 x 32	1.5	100



Aluminium Caps for Crimp Top Vials 12 x 32 mm

Crimp Caps ND11 are made of Aluminium with septa. Caps have open top centre hole.


Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.V11.SC1	Silver	Open Top	ND11	Silicone White / PTFE Red	100
GW.V11.SC2	Silver	Open Top	ND11	PTFE Red / Silicone White / PTFE Red	100
GW.V11.SC3	Silver	Open Top	ND11	Transparent TEF/Natural Rubber Red-orange Septa	100
GW.V11.SC4	Silver	Open Top	Nd11	Silicone Red / PTFE Clear	100
GW.V11.SC5	Silver	Open Top	Nd11	Silicone Red / PTFE White	

SNAP RING/Cap Vials

- A snap cap provides a moderate seal and doesn't require special tools.
- Snap Caps can be snapped onto a vial and removed by hand thus eliminate the need for crimping and decapping tools.
- With an opening in the top, snap caps provide easy access for autosamplers.

Snap Ring Vials 12 x 32 mm (ND 11)

- Snap ring vials ND11 are available in clear and amber type-1 glass can be used on almost any autosamplers
- Also suitable for robotic handling.
- These vials can also be used as crimp neck vial.
- Snap ring vials are suitable for HPLC applications



Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Pack Unit
GW.SRVC11.01	Clear	-	ND11	11.6 x 32	1.5	100
GW.SRVC11.02	Clear	Yes	ND11	11.6 x 32	1.5	100
GW.SRVA11.01	Amber	-	ND11	11.6 x 32	1.5	100
GW.SRVA11.02	Amber	Yes	ND11	11.6 x 32	1.5	100

PP Caps for Snap Ring/Cap Vials 12 x 32 mm

Snap ring caps ND11 are made of PP & Aluminium and are supplied with various septa. Snap ring caps are very easy to use and time-saving.

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.V11.CC2	Clear	Open Top	ND11	Silicone White / PTFE Red	100
GW.V11.CC3	Clear	Open Top	ND11	Silicone White / PTFE Red/pre-slit	100
GW.V11.CC4	Clear	Open Top	ND11	PTFE Red / Silicone White / PTFE Red	100

HEAD SPACE Vials

- Designed to use for High Pressure / Performance Liquid Chromatography (HPLC) and Gas Chromatography (GC) Applications.
- Crimp top vials consist of an Aluminum cap and PTFE / Silicon / Rubber Septum thus forms a superior seal and reduces evaporation.
- Optional with White Ceramic Writing Patch and fill marks for easy labelling.

Head Space Vials ND 20

Headspace vials with Crimp neck ND20 available with round or flat bottom are made of clear and amber Type 1 glass. The vials are available with a flat DIN crimp neck or with a bevelled neck (HS neck).

On the flat DIN crimp neck, liners have a greater contact surface, which ensures a better seal. A bevelled HS neck is required if the overpressure safety seal.



Code	Glass Color	Patch	Neck Type	D x H (mm)	Working Volume (ml)	Bottom Design	Finish	Special	Pack Unit
GW.HPVC20.01	Clear	-	ND20	22 X 38.5	5	Round	HS	-	100
GW.HPVC20.02	Clear	-	ND20	22 X 38	5	Flat	HS	-	100
GW.HPVC20.03	Clear	-	ND20	20 X 38	5	Flat	DIN	-	100
GW.HPVC20.04	Clear	-	ND20	22.5 X 46	10	Round	DIN	-	100
GW.HPVC20.05	Clear	-	ND20	22.5 X 46	10	Flat	DIN	Long Neck	100
GW.HPVC20.06	Clear	-	ND20	20 X 54.5	10	Flat	DIN	-	100
GW.HPVC20.07	Clear	-	ND20	22.5 X 75.5	20	Flat	DIN	Long Neck	100
GW.HPVC20.08	Clear	-	ND20	22.5 X 75.5	20	Round	DIN	Long Neck	100
GW.HPVC20.09	Clear	-	ND20	22.75 X 75.5	20	Round	HS	-	100
GW.HPVC20P.10	Clear	Yes	ND20	22.75 X 75.5	20	Round	HS	-	100

Standard Caps for Head Space Vials ND 20

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.V20.SC1	Aluminum	Open Top	ND20	Moulded Grey PTFE / Butyl	100
GW.V20.SC3	Aluminum	Open Top	ND20	Moulded Grey PTFE / Butyl	100
GW.V20.SC4	Aluminum	Open Top	ND20	Moulded Grey PTFE (Pharma Fix) / Butyl	100

Magnetic Crimp Cap ND 20

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.V20.GC1	Gold	Open Top	ND20	White PTFE/Blue Silicone Septa	100
GW.V20.GC2	Gold & Silver	Open Top	ND20	Natural PTFE/Natural Silicone Septa	100



Pressure Release Cap for Headspace Vials ND 20

Code	Cap Color	Cap Design	Neck Type	Liner	Pack Unit
GW.V20.SC5	Silver	Open Top	ND20	Moulded Grey PTFE/Butyl Septa	100



SCINTILLATION VIALS

Sample Vials (With Unattached Caps)

• Boro 5.0 Expansion

These are general-purpose vials are available various sizes & capacities. these are made from clear type -1 glass. cap is flushed with vial shoulder. Supplied with black phenolic caps with polyvinyl-faced pulp liner. The cap on the vial is the same diameter as the body, fitting flush to the shoulder.



Item Code	Size	Cap Colour	Cap Material	Linear Material	Thread Size	Case Size
086.211.01	12x35	Black	Phenolic	Polyvinyle faced Pulp	8-425	5184
086.211.02	15x45	Black	Phenolic	Polyvinyle faced Pulp	13-425	2304
086.211.03	16x50	Black	Phenolic	Polyvinyle faced Pulp	13-425	2592
086.211.04	17x60	Black	Phenolic	Polyvinyle faced Pulp	15-425	1728
086.211.05	19x65	Black	Phenolic	Polyvinyle faced Pulp	15-425	1728
086.211.06	21x70	Black	Phenolic	Polyvinyle faced Pulp	18-400	1152
086.211.07	23x85	Black	Phenolic	Polyvinyle faced Pulp	20-400	864
086.211.08	25x95	Black	Phenolic	Polyvinyle faced Pulp	22-400	576
086.211.09	28x57	Black	Phenolic	Polyvinyle faced Pulp	24-400	432
086.211.10	28x70	Black	Phenolic	Polyvinyle faced Pulp	24-400	432
086.211.11	28x95	Black	Phenolic	Polyvinyle faced Pulp	24-400	432
086.211.12	28x108	Black	Phenolic	Polyvinyle faced Pulp	24-400	432

Scintillation Vial (With Unattached Caps)

• Boro 3.3 Expansion

Item Code	Size	Cap Colour	Cap Material	Linear Material	Thread Size	Case Size
086.202.06*	17x57	No Cap	-	-	15-425	1000
086.209.01	17x57	White	Urea	Foil Cork	15-425	1000
086.202.20*	28x61	No Cap	-	-	22-400	500
086.208.01	28x61	White	Urea	Cone liner	22-400	500
086.210.01	28x61	White	PP	Foil Pulp	22-400	500
086.210.02	28x61	White	PP	PE	22-400	500
086.208.02	28x61	White	PP	Foil Pulp	22-400	500
086.208.03	28x61	White	Urea	Foil Cork	22-400	500
086.208.04	28x61	White	Urea	Foil Cork	24-400	500

Ps : * Is offered without cap

TOC Vials, Amber & Clear Glass

- Borosilicate Glass Exp. 5.1
- With Bonded Cap & PTFE Lined Silicone Septa having flip cover

Item Code	Capacity (ml)	O.D. Height	Glass	Septa	Pack
083.205.02	20	28x57	Clear	Bonded	1x72
083.205.04	40	28x95	Clear	Bonded cap with flip cover	1x72
083.206.02	20	28x57	Amber	Bonded	1x72
083.206.04	40	28x95	Amber	Bonded cap with flip cover	1x72



Plain Neck Flat Bottom

Code	Tube Dia x Height (mm)	Neck Dia (mm)	Cap. (ml)	Pcs/Inner	Inner/Case	Case Qty (Pcs.)
086.802.03	16.25x35	13	3ml	209	15	3135
086.804.14	22.50x75	20	20ml	100	9	900
086.804.17	24.00x33	20	5ml	91	15	1365
086.804.18	24.00x45	20	13ml	91	12	1092
086.804.19	24.00x60	20	18ml	91	9	819
086.804.20	24.00x74	20	20ml	145	5	725
086.804.21	30.00x62	20	30ml	60	9	540

For Ordering Amber Vials, Please Add 'A' After Cat. Number & For Light Amber Vails, Please Add 'LA' After Cat. Number



Plain Neck Round Bottom

Item Code	Tube Dia (mm)	Plain Neck	1 Inner Pcs	Inner	Case
086.806.16	22.50	20ml	100	9	900

Dental Cartridge

Item Code	Tube Dia (mm)	1 Inner Pcs	Inner	Case
619.DC.01	8.65x62	1230	4920	7200

Dental Cartridge Closure Cap

Item Code	Cap Colour	1 Inner Pcs	Inner	Case
619.DCC.01	Silver	1230	4920	7200

Dental Cartridge Plunger

Item Code	Cap Colour	1 Inner Pcs	Inner	Case
619.DCP.01	Grey	1230	4920	7200





GLASSCO LABORATORY EQUIPMENTS PVT. LTD.

Head Off.: Manglai, P. O. Khuda Kalan, Ambala Cantt. 133 004 (INDIA)



 www.glasscolabs.com

 www.glasscolabs.in

 enquiry@glasscolabs.com